


neria™ soft 90

Intended use:

The product is intended for the subcutaneous infusion of medication such as immunoglobulins for the treatment of Primary Immune Deficiency (for solutions up to 16 % only), apomorphine for Parkinson's Disease and morphine for pain management therapy.

The neria™ soft 90 infusion set has been tested to be compatible with the above-mentioned drugs approved for subcutaneous infusion.


Product range:

Product code	Description	Cannula length	Tubing length	EAN	Pcs per inner box	Pcs per outer box
neria™ soft 90						
722030-5226	neria™ soft 90 infusion set	6 mm	30 cm	05705244003637	10	100
722030-5229	neria™ soft 90 infusion set	9 mm	30 cm	05705244003651	10	100
722060-5226	neria™ soft 90 infusion set	6 mm	60 cm	05705244003675	10	100
722060-5229	neria™ soft 90 infusion set	9 mm	60 cm	05705244003699	10	100
722110-5226	neria™ soft 90 infusion set	6 mm	110 cm	05705244003712	10	100
722110-5229	neria™ soft 90 infusion set	9 mm	110 cm	05705244003736	10	100
neria™ soft 90 cannula						
722000-5226	neria™ soft 90 cannula	6 mm	N/A	05705244003835	10	100
722000-5229	neria™ soft 90 cannula	9 mm	N/A	05705244003859	10	100
neria™ soft 90 tubing						
722030-5220	neria™ soft 90 tubing	N/A	30 cm	05705244003774	10	200
722060-5220	neria™ soft 90 tubing	N/A	60 cm	05705244003798	10	200
722110-5220	neria™ soft 90 tubing	N/A	110 cm	05705244003811	10	200

Product specifications:

Insertion angle 90 degrees
Soft cannula length 6 mm and 9 mm
Soft cannula outer diameter ~ 0,68 mm
Introducer needle 27 gauge
Tubing length in cm 30, 60 and 110
Tubing inner diameter 0,385 mm
Tubing outer diameter 1,50 mm

Material specifications:

PVC free Yes
Made with natural rubber latex No
Raw materials or derivation of animal origin No
The product contains phthalates No
Introducer & connector needles may contain . . ≤ 10% nickel

Certifications:

CE-marking Yes
Product classification Class IIB (EU)
Quality system certified to ISO 13485
Tariff code HS9018390000
GMDN 35833

Sterilization:

Sterile Yes
Sterilized by Ethylene oxide
Shelf life 3 years
Usage time (drug depended) Maximum 72 hours

Drug compatibility/stability tested with:

Apomorphine
Morphine (hydromorphone, morphine sulphate & morphine chloride)

Storage and handling:

	AVOID DIRECT SUNLIGHT - Keep the product away from direct sunlight and any humid environment
	KEEP DRY - Store in a dry place at room temperature
	DO NOT RE-USE
	DO NOT USE - If package is damaged
	CAUTION - See Instruction for Use



Unomedical a/s A ConvaTec Company

Aaholmvej 1-3
 Østved, 4320 Lejre, Denmark
www.neria.com