

## Behandlerstrømpe sæt

Juzo Ulcer System • Juzo Ulcer System Add

### Materiale

#### Kompressionsstrømpe

- 75% Polyamid
- 25% Elastan

#### Liner

- 73% Polyamid
- 27% Elastan



Vaskegrader



Anvend ikke skyllemiddel



Må ikke tørretumbles




Lufttørres



Stryges ved meget lav temperatur



### Omkreds på ben - mål i cm



	XS	S	M	L	XL	XXL
<b>cD</b>	27,0 - 29,5	29,5 - 32,5	32,5 - 35,5	35,5 - 39,0	38,5 - 45,0	41,5 - 45,5
<b>cC</b>	29,0 - 31,5	31,5 - 34,5	34,5 - 37,5	37,5 - 41,0	40,5 - 44,5	43,0 - 47,0
<b>cB<sup>1</sup></b>	23,5 - 25,5	25,5 - 28,0	27,5 - 30,0	30,0 - 33,0	32,5 - 35,5	35,0 - 38,5
<b>cB</b>	18,5 - 20,0	20,0 - 22,0	22,0 - 24,0	24,0 - 26,0	26,0 - 28,5	28,0 - 30,5
<b>cY</b>	28,0 - 30,5	29,0 - 31,5	31,0 - 34,0	33,0 - 36,0	35,0 - 38,5	37,0 - 40,5
<b>cA</b>	18,5 - 20,0	20,0 - 22,0	22,0 - 24,0	24,0 - 26,0	26,0 - 28,5	28,0 - 30,5

### Længde på fod

**lZ** ————— Liner: 22,5 - 27,5 cm

**lA** ————— Kompressionsstrømpe: 15,0 - 21,0 cm

### Længde på ben **lD** - mål i cm

Kort model	Middel model
34,0 - 39,0	40,0 - 46,0