

evercare[®] Medical face mask with earloops, type IIR

REF: 2881

evercare[®] Medical face mask with earloops, type II R provides a barrier to minimize the direct transmission of infective agents between medical staff and patient. It is a splash resistant mask offering protection against potentially contaminated fluid splashes. Recommended to be used during wet surgery. The mask is slightly shorter in width than many other face masks but at the same time it is also wider when in open size. The design of the side sealing and the anchoring of the earloops help to reduce the gaping at the cheeks and makes the mask more conforming to the face. This overall design results in a deeper cup-shape, which is better conforming and provides a fitter coverage of the nose and mouth.

Key Features

- ▶ EN 14683:2019
- ▶ Bacterial Filtration Efficiency (BFE) ≥ 98%
- ▶ Single use
- ▶ Flexible nose wire for proper fit
- ▶ Splash resistant
- ▶ Designed to give an excellent /improved fit of the face mask for the user



Ear-loop



Splash resistant

REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
2881	Medical face mask with earloops, type IIR, non-sterile, White	165 x 90 mm	50 / 2500

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-5 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR Class I
Compliance with product standards:	EN 14683:2019
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417

Material

Product Part	Raw Material & Characteristics
Outer facing	Spunbond PP nonwoven, white
Filter media	Meltblown PP nonwoven, white
Inner facing	Spunbond PP nonwoven, white
Elastic earloops	Spandex, Polyester and Polymide, latex-free, white
Nose wire	PP coated Galvanized iron wire, 105 mm

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	PE bag, Carton, TCF or ECF quality	188 x 98 x 92
Transport Carton	Carton, TCF or ECF quality	515 x 395 x 485
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
2881	N/A	06438129324563	06438129324567

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.