

PRODUCT DATA SHEET

LUOFUCON® Foam Dressing, which is EO sterilized, is composed of soft and highly absorbent polyurethane. It can bulkily absorb wound exudates and maintain moist environment, which accelerates wound healing. It is especially suitable for excessive exuding wounds.



Benefits:

- Various in forms, shapes and sizes;
- Highly absorbent (more than 10 times of fluid);
- Transmitting moisture vapour and oxygen and providing moist environment for wound healing;
- Occlusive for bacteria proof;
- Comfortable and conformable.



Indications:

LUOFUCON® Foam Dressing is used for heavily exuding wounds, especially during the inflammatory phase following debridement and desloughing, when drainage is at its peak, including:

- Deep cavity wounds;
- Pressure ulcers;
- Leg ulcers;
- Donor sites;
- Lacerations and abrasions.

ITEM CODE	PRODUCT		SIZE	THICKNESS	GLUE	Unit	qty/ ctn	Pcs/ Dispenser
FD050050	PU FOAM NORMAL	NON ADHESIVE	5x5CM	5MM	N/A	pc	300	10
FD100100	PU FOAM NORMAL	NON ADHESIVE	10X10CM	5MM	N/A	pc	400	10
FD100200	PU FOAM NORMAL	NON ADHESIVE	10X20CM	5MM	N/A	pc	100	5
FD150150	PU FOAM NORMAL	NON ADHESIVE	15X15CM	5MM	N/A	pc	300	10
FD200200	PU FOAM NORMAL	NON ADHESIVE	20X20CM	5MM	N/A	pc	100	5
FD050050U	PU FOAM THICK/PU	NON ADHESIVE	10X10CM	5MM	N/A	pc	300	10
FD100100U	PU FOAM THICK/PU	NON ADHESIVE	10X10CM	5MM	N/A	pc	400	10
FD100200U	PU FOAM THICK/PU	NON ADHESIVE	10X20CM	5MM	N/A	pc	100	5
FD150150U	PU FOAM THICK/PU	NON ADHESIVE	15X15CM	5MM	N/A	pc	300	10
FD200200U	PU FOAM THICK/PU	NON ADHESIVE	20X20CM	5MM	N/A	pc	100	5
FD170100UH	PU FOAM THICK/PU	NON ADHESIVE	10,5X17CM	5MM	N/A	pc	50	5
FD70110UH	PU FOAM THICK/PU	NON ADHESIVE	7X11CM	5MM	N/A	pc	300	10
FD100090UX	PU FOAM THICK/PU	NON ADHESIVE	10X9CM	5MM	N/A	pc	300	10
FD065065UX	PU FOAM THICK/PU	NON ADHESIVE	6,5X6,5	5MM	N/A	pc	600	10
FD100100B	PU FOAM THICK/PU	SELF ADHESIVE	10X10CM	5MM	Acrylic	pc	300	10
FD150150B	PU FOAM THICK/PU	SELF ADHESIVE	15X15CM	5MM	Acrylic	pc	300	10
FD200200B	PU FOAM THICK/PU	SELF ADHESIVE	20X20CM	5MM	Acrylic	pc	100	5

Storage

Store the product in a cool and dry place away from sunlight exposure. Recommended storage temperature 5°C - 40°C.

Shelf-life

3 years.

Batch numbering method: YYY MM AAA

YYY- Year of production

MM- Month of production

AAA-Serial number of batch, based on annual production

Example: 201605087 - Manufactured 2016-05 as 87th batch in year 2016.

Compliance with Standards & Regulations

Products are in full compliance with following standards	
ISO 10993-5	Biological evaluation of medical devices – Part 5: Tests for in vitro cytotoxicity
ISO 10993-7	Biological evaluation of medical devices – Part 7: Ethylene oxide sterilization residuals
ISO 10993-10	Biological evaluation of medical devices – Part 10: Test for irritation and
ISO 11607-1	Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems
ISO 11135	Sterilization of health-care products -- Ethylene oxide -- Requirements for the development, validation and routine control of a sterilization process for medical devices

Packaging System

ITEM CODE	SIZE	Unit	qty/ ctn	Pcs/ Dispenser	pouch size(mm)		Dispenser box size(mm)			Carton size(cm)		
FD050050	5x5CM	pc	300	10	115	95	95	60	120	48,8	37,4	13
FD100100	10X10CM	pc	400	10	160	145	146	60	163	60,8	32,5	35,6
FD100200	10X20CM	pc	100	5	265	150	151	30	268	32,3	33	29,6
FD150150	15X15CM	pc	300	10	213	185	186	60	216	58,1	32,5	46,5
FD200200	20X20CM	pc	100	5	265	250	251	30	268	52,3	33	29,6
FD050050U	10X10CM	pc	300	10	115	95	95	60	120	48,8	37,4	13
FD100100U	10X10CM	pc	400	10	160	145	146	60	163	60,8	32,5	35,6
FD100200U	10X20CM	pc	100	5	260	145	151	30	268	32,3	33	29,6
FD150150U	15X15CM	pc	300	10	210	195	186	60	216	58,1	32,5	46,5
FD200200U	20X20CM	pc	100	5	265	250	251	30	268	52,3	33	29,6
FD170100UH	10,5X17CM	pc	50	5	240	175	177	60	245	51,2	32,5	20,5
FD70110UH	7X11CM	pc	300	10	110	160	177	60	245	44,3	37,5	35,4
FD100090UX	10X9CM	pc	300	10	157	135	140	70	162	44,3	37,5	35,4
FD065065UX	6,5X6,5	pc	600	10	112	110	115	70	127	48,4	37,5	41,3
FD100100B	10X10CM	pc	300	10	160	145	146	80	163	46,2	42,5	35,6
FD150150B	15X15CM	pc	300	10	213	185	186	80	216	58,1	42,5	46,5
FD200200B	20X20CM	pc	100	5	265	250	251	40	268	52,3	43	29,6

HUIZHOU FORYOU MEDICAL DEVICES CO., LTD.

ADD: No.1, Shangxia North Road, Dongjiang Hi-tech Industrial Park, Huizhou, Guangdong, China

TEL: +86-752-5302188 FAX: +86-752-5302088 E-MAIL: medical2@foryougroup.com

www.foryoumedical.com

DISTRIBUTOR WITHIN EU:

OneMed Sverige AB

adress: Tagenevägen 29; 425 37 Hisings Kärra

website: www.onemed.se