


evercare® Hand/foot set

REF: 1232-50

evercare Hand/foot set, sterile is a draping set for hand/foot operations.

Key Features

- ▶ Elastic fenestration
- Absorbent patch
- ▶ Tube holders
- Adhesive OP Tape
- ▶ Main drape complies to EN 13795-1 High performance critical product area requirements




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
1232-50	Hand/foot set, sterile	See Material table	6/12

Sales Unit: Transport Carton

Page 1 2023-12-05

Regulations


The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.























REACH regulation:

The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation:

EN ISO 10993-01 EN ISO 10993-5 EN ISO 10993-7 EN ISO 10993-10

Compliance with legal requirements:

Medical Device Regulation (EU) 2017/745, MDR class Is

Compliance with product standards:

EN13795-1

Quality standards:

EN ISO 13485

Label information and packaging compliance:

EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

Material

Product Part	Raw Material & Characteristics	
1 Mayo cover	78x145 cm, reinforced, Laminate: PP nonwoven/PE-film/PE-film, blue, absorbent patch $65x85$ cm	
1 Adhesive Op-tape	9 x 50 cm, Nonwoven polyester, blue, adhesive coated	
4 Hand towel	30 x 34 cm, Pulp, natural white, thread reinforced	
1 Hand/foot drape	230 x 300 cm, Laminate: blue PP nonwoven /blue PE-film, blue absorbent patch75 x 80 cm, PP nonwoven, blue tube holders PP, aperture 3.5 cm, elastic material SEBS based thermoplastic elastomer	
1 Instrument table cover	150 x 190 cm, Laminate: PP nonwoven/PE-film, blue	

Page 2 2023-12-05

Warnings & Recommendations

Warnings conta	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur. Do not use if package is damaged.		
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.		
Disposal Recommendation	Dispose in accordance with local regulations.		
Shelf Life	5 years		

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, PE/PET film	420 x 603
Dispenser Pack	TCF or ECF quality	485 x 285 x 287
Transport Carton	TCF or ECF quality	500 x 310 x 590
Packing Tape	Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
1232-50	06438129118261	06438129318265	06438129518269

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.