

# evercare<sup>®</sup> Adhesive Plasters nonwoven

REF: CMT352177, CMT352184, CMT352191, CMT352207

evercare Adhesive plaster nonwoven is a dressing intended to protect smaller wounds from impurities in outer surrounding. The absorbent patch attached to the plaster absorbs excudate from the covered wound. The plaster material is soft and comfortable for the skin. The plaster is typically used for post-injection site protection.

## **Key Features**

- Soft, high breathable nonwoven plaster
- Solvent-free adhesive
- Non-sticking absorbent pad
- ▶ Two different shapes: stripe and rounded
- Individual pouch for each plaster



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
CMT352177	Adhesive Plaster, nonwoven, strip, Beige	2,4x7cm	200 / 10200
CMT352184	Adhesive Plaster, nonwoven, stripe, Beige	1,2x4cm	200 / 10200
CMT352191	Adhesive Plaster, nonwoven, stripe, Beige	3,8x3,8cm	200 / 10200
CMT352207	Adhesive Plaster, nonwoven, stripe, Beige	2x4cm	200 / 10200

Sales Unit: Transport Carton

Page 1 2024-12-12

### Regulations

The product complies with legal requirements for medical devices.







**REACH** regulation:

The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

EN ISO 10993-01 EN ISO 10993-5 **Biological evaluation:** 

EN ISO 10993-10

Compliance with legal requirements:

Medical Device Regulation (EU) 2017/745, MDR Class I

**Compliance with product** standards:

EN 14079, ASTM D3330

EN ISO 13485 **Quality standards:** 

Label information and packaging compliance:

EN ISO 15223-1, EN ISO 20417

#### Material

Product Part	Raw Material & Characteristics	
Main material	Nonwoven, 100 Polyester, beige	
Adhesive	Hot-melt acrylic glue, UV cross linked	
Absorbent pad	Viscose/PP/PE, layered with non-stick LDPE	
Peel-paper	Silicone paper, one sided	

Page 2 2024-12-12

## Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.				
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.				
Disposal Recommendation	Dispose in accordance with local regulations.				
Shelf Life	5 years				

### **Packaging Information**

Package	Material	Size ( L x W x H ) mm
Single Pack	N/A	N/A
Dispenser Pack	Corrugated cardboard, TCF or ECF quality	93 x 165 x 35
Transport Carton	Corrugated cardboard, TCF or ECF quality	580 x 185 x 300
Packing Tape	Transparent, PVC-free	N/A

#### Barcodes / UDI-DI

REF	Single Pack	Dispenser Pack	Transport Carton
CMT352177	N/A	8017990117337	18017990117334
CMT352184	NA	8017990117344	18017990117341
CMT352191	N/A	8017990117368	18017990117365
CMT352207	N/A	8017990117382	18017990117389

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.

Page 3