

SELEFA® Optimia Gauze swabs

REF: 222001, 222002, 222003, 222004, 222005

SELEFA® Optimia Gauze swabs are produced from environmental friendly bleached cotton fibres, woven in the amount of 17 threads per square cm. They are produced in four measurements with folded edges to avoid threads loose. Quality cotton gauze swabs are produced in four to sixteen layers that gives its softness, increasing its absorbent properties.

Key Features

- ▶ Wide area of use: pre-operative skin cleansing, wound cleansing, drying, absorption, protection etc
- ▶ Effective in moisture and liquid absorption, can be sterilized
- ▶ Peroxide bleached cotton gauze without chlorine
- ▶ Woven edges prevent loose threads drawing
- ▶ 100 % cotton, low linting



REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
222001	SELEFA® Optimia Gauze swabs, White	5x5cm,8ply	120 / 10800
222002	SELEFA® Optimia Gauze swabs, White	7,5x7,5cm,8ply	120 / 7200
222003	SELEFA® Optimia Gauze swabs, White	10x10cm,8ply	120 / 7200
222004	SELEFA® Optimia Gauze swabs, White	10x20cm,12ply	50 / 1000
222005	SELEFA® Optimia Gauze swabs, White	10x10cm 16ply	60 / 3600

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-5 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class I
Compliance with product standards:	EN 14079
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417

Material

Product Part	Raw Material & Characteristics
100 % cotton gauze, white	100 % cotton gauze, white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	TCF quality, ECF quality	222001: 162 x 60 x 54 222002: 235 x 60 x 77 222003: 204 x 85 x 104 222004: 204 x 90 x 104 222005: 204 x 85 x 104
Transport Carton	TCF quality, ECF quality	635 x 545 x 365
Packing Tape	Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
222001	N/A	06438129353921	06438129553925
222002	N/A	06438129353938	06438129553932
222003	N/A	06438129353945	06438129553949
222004	N/A	06438129353952	06438129553956
222005	N/A	06438129353969	06438129553963

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.