

Blood collection devices and accessories

Safety blood collection sets

BD Vacutainer® Safety-Lok™ Blood Collection Set

The BD Vacutainer® Safety-Lok™ Blood Collection Set for venous blood collection has a fully integrated safety shield, which once activated, protects against needlestick injuries.

- Indication of successful venepuncture:
When the vein has been successfully penetrated, blood can be seen in the device.
- Versatile:
For taking blood samples and for short-term infusions of up to 2 hours.
- Single-handed activation possible:
Activating the safety mechanism with one hand allows you to focus more attention on the patient and the venepuncture site.
- Protection against needlestick injuries:
Following successful venepuncture, the integrated safety shield is pushed over the needle, covering it completely, indicated by an audible click.








Further information

BD recommends for all needles without a pre-attached holder that a BD Vacutainer® One-Use Holder is used.

Clinical and technical information is available on request.

BD Vacutainer® Safety-Lok™ Blood Collection Sets

Cat. no.	Size	Needle length	Length of tubing	With Luer adapter	Colour code	Box/Case
367282	21 G (0.8 mm)	19 mm	178 mm	Yes		50/200
367286	21 G (0.8 mm)	19 mm	305 mm	Yes		50/200
367246	21 G (0.8 mm)	19 mm	305 mm	No		50/200
367284	23 G (0.6 mm)	19 mm	178 mm	Yes		50/200
367288	23 G (0.6 mm)	19 mm	305 mm	Yes		50/200
367247	23 G (0.6 mm)	19 mm	305 mm	No		50/200
367295	25 G (0.5 mm)	19 mm	178 mm	Yes		50/200
368383	25 G (0.5 mm)	19 mm	305 mm	No		50/200

Blood collection devices and accessories

Safety blood collection sets

BD Vacutainer® Safety-Lok™ Blood Collection Set with Pre-Attached Holder

With this Safety Blood Collection Set, the holder is already pre-attached, so it is not necessary to manually assemble the needle and holder. This ready-for-use blood collection set is individually packaged in a sterile blister.





The sterile closed system is ideally suited for the taking of samples using BD BACTEC™ Blood Culture Bottles.

Further information

Technical information is available on request.



BD Vacutainer® Safety-Lok™ Blood Collection Sets with Pre-Attached Holder

Cat. no.	Size	Needle length	Length of tubing	Colour code	Box/Case
368654	21 G (0.8 mm)	19 mm	178 mm		25/200
368652	21 G (0.8 mm)	19 mm	305 mm		25/200
368655	23 G (0.6 mm)	19 mm	178 mm		25/200
368653	23 G (0.6 mm)	19 mm	305 mm		25/200

**Sterile
closed
system**