



Technical datasheet

Item code: 730220

MEDI-INN Disposable Underpad SAP-Plus Premium

Description

Absorbent body made of soft cellulose-flake-superabsorbent mixture (carrier tissue), full surface covering with textile dry fleece to the edges; closed all around; impermeable, non-slip outer film; pre-folded ready for cabinet; different qualities for economical usage.

Applicable Standards

Class I medical product in accordance with Regulation (EU) 2017/745, suitable for work in the medical/nursing sector.

Product Information

Item code	730220
Weight incl. SAP	~ 15 gr.
Absorb volume	200 ml
Color	blue
Dimensions in cm	
Width	20
Length	50
Material	
Laundry protection	Polyethylene film, blue, non-slip
Absorbent	SAP-mixture (tissue and fluff bleached without chlorine)
Fleece	Polypropylene, hydrophilic (water absorbing or water permeable)
Glue	Hot glue

Packaging / Storage

Item code	730220
Unit	Piece
Quantity of single packaging	15
Quantity in carton	40
Type of encase packaging	Vacuum bags
Carton size	50 x 22 x 47 cm
Shelf life	5 years under the conditions below
Labeling	LOT Date of manufacture Best before date Single use
Disposal	Disposal in accordance with applicable regulations. Products contaminated by chemical or biological substances must be disposed of in accordance with the regulations for the relevant chemicals and/or infectious material.
Storage	Dry between 10-30° C. Avoid direct sun exposure.



Quality Assurance

Product specific standards	EN 1041 EN ISO 10993-10 EN 15233-1 EN ISO 14971
Quality certificates	EN ISO 9001:2015 EN ISO 13485:2016

Notes

Before use, check the disposable underpads critically for its suitability for the intended activity, areas of application and associated risks. Check the disposable underpads for damage before use and after putting it on under no circumstances use unsuitable or damaged disposable underpads! The disposable underpads are for single use certainly.

Contact us for further information.

Created by: BODY PRODUCTS RELAX Pharma und Kosmetik GmbH
Alfred-Nobel-Straße 1 – 3
D-50226 Frechen
Germany

Date, Place Frechen, 09.07.2021