


# 613100 Apron

Protective Apron

## Product details

**Size :** 145cm  
**Colour :** White  
**Sterile :** Non-sterile  
**Brand :** BARRIER

## Images



## Delivered items

613100-03

**Country of origin:** Czechia

**Shelf life:** 5 years

**Sterilization method:** Non-sterile

**Production Responsibility:** Intergaze, s.r.o., Ul. Svobody, 73301 Karviná, Czech Republic

**Packing information:** First packaging layer is a plastic liner/bag. Second layer is a corrugated board transport box.

**Is suitable for Tray:** No

Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Piece	1	7332430090683			

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.  
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2022 Mölnlycke Health Care AB. All rights reserved.



Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Transport box	150	7323190211760	389x399x278 mm	43.1 dm3	9.2 / 8.6 kg
Pallet	5400	7323190211753	800x1200x1812 mm		

613100-04

Country of origin: Czechia

Shelf life: 5 years

Sterilization method: Non-sterile

Production Responsibility: Intergaze, s.r.o., Ul. Svobody, 73301 Karviná, Czech Republic

Packing information: First packaging layer is a plastic liner/bag. Second layer is a corrugated board transport box.

Is suitable for Tray: No

Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Piece	1	7332430090683			
Transport box	150	7323190248759	389x399x278 mm	43.1 DECIMETER3	9.2 / 8.6 KILOGRAM
Pallet	5400	7323190248742	800x1200x1812 mm		

## Material

Animal tissues :	No
Human blood derivatives :	No
Natural rubber latex :	No
Medicinal substances :	No
Phthalates :	No
Polyvinyl chloride :	No

## Product Composition Wearing Apparel, Apron

Product Component	Composition
Main material	Tissue nonwoven 28 g/m <sup>2</sup> , Polyethylene film 25 µm

## Technical

## Dimension

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.  
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2022 Mölnlycke Health Care AB. All rights reserved.



Dimension text	Dimension value
Length	145 cm
Width	75 cm

### Disposal instructions

Non-hazardous waste used BARRIER products and sterility barriers should, in the majority of cases, be classified as non-hazardous waste. They contain high amounts of energy and are well suited for incineration. BARRIER products do not contain any hazardous substances that can leach out if the products are land filled. Transport boxes are designed to fit existing recovery systems. The new BARRIER packaging system complies with the Packaging Waste Directive of the European Union.

### Storage instructions

Mölnlycke Health Care recommends that BARRIER products are stored under normal storage conditions. All layers of packaging should be kept intact until access to the underlying layers is needed. Storage facilities for products only protected by the sterility barrier should be kept under conditions where low level of particulate air contamination prevail, so that it would not constitute a risk to the patient when the package is opened and the product is used.

### Classifications

Regulation type(s)	MDR Class I ns	Locally Regulated
Intended Purpose :	The product should be worn in situations where there is a risk that infectious agents can be transferred throughout the hospital. Secondly the product also protects the clothes of the user from different kind of fluids.	The product should be worn in situations where there is a risk that infectious agents can be transferred throughout the hospital. Secondly the product also protects the clothes of the user from different kind of fluids.
MDR Classification Rule :	1	
Conformity Annexes :	IV	
Measuring Function :	No	
Notified body medical devices/PPE :	MPA	Local Authority Mexico
510(k) clearance number :		213300401J0021

**Applied standards :** The standards presented below is a selection of the most essential standards that are adhered to.

EN 1041, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 10993-5, EN ISO 15223-1, ISO 15223-2, EN ISO 10993-10, ISO 14001

### GMDN Code (Global Medical Device Nomenclature)

40511 Hygiene apron, single-use

### UNSPSC

42131601 Medical staff aprons or bibs

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

## Commodity Code

4818500000 Wearing apparel mainly made of tissue/pulp paper

## CE Responsibility

Mölnlycke Health Care AB, Gamlestadsvägen 3c, Box 13080, SE-402 52 Göteborg Sweden

## Basic UDIDI

733243000000000097KH

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.  
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.  
©2022 Mölnlycke Health Care AB. All rights reserved.

