


708850 Adhesive Aperture Drape

Product details

| | |
|------------------------------|-------------------------------------|
| Size : | 150cm x 180cm, Aperture 5x15cm |
| Descriptive feature : | Adhesive, Centered ellipse aperture |
| Sterile : | Sterile |
| Brand : | BARRIER |
| Equivalences : | 968850 |

Images



Delivered items

708850-21

Sales released in: Australia, Austria, Bahrain, Belgium, Czechia, Denmark, Estonia, Finland, France, Germany, Hungary, Iceland, India, Ireland, Italy, Latvia, Lithuania, Luxembourg, Macedonia (the former Yugoslav Republic of), Malaysia, Martinique, Morocco, Netherlands, New Zealand, Norway, Poland, Portugal, Qatar, Romania, Saudi Arabia, Singapore, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey, United Kingdom of Great Britain and Northern Ireland

Shelf life: 5 years

Sterilization method: Irradiation

Production Responsibility: Mölnlycke Health Care (Thailand) Ltd, 160 Moo 17 Bangplee Industrial Estate, Bangna-Trad Road, Km 23 Tambol Bangsaothong, Amphur Bangsaothong, Samutprakarn 10570, Thailand

Packing information: First packaging layer is a peel-open sterile barrier, plastic/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a plastic liner/bag. Third layer is a corrugated board dispenser box. Fourth layer is a corrugated board transport box.

Is suitable for Tray: No

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2022 Mölnlycke Health Care AB. All rights reserved.



| Packing level | Quantity | GS1 Code / UDI-DI | Width x Length x Height | Vol | Weight gross / net |
|---------------|----------|-------------------|-------------------------|------------------------|--------------------|
| Piece | 1 | 7332551041922 | | | |
| Consumer pack | 1 | 7323190294152 | | | |
| Dispenser box | 28 | 7323190293186 | | | |
| Transport box | 56 | 7323190293179 | 397x580x329 mm | 75.8 dm ³ | 12.6 / 10.9 kg |
| Pallet | 1120 | 7323190293162 | 800x1200x1789 mm | 1717.4 dm ³ | 277.8 / 200.0 kg |

Technical

Dimension

| Dimension text | Dimension value |
|-----------------|-----------------|
| Outer dimension | 150 cm x 180 cm |
| Aperture | 5 cm x 15 cm |

Disposal instructions

Used Mölnlycke products and sterile barriers should, in the majority of cases, be classified as non-hazardous waste. They contain high amounts of energy and are well suited for incineration. Mölnlycke products do not contain any hazardous substances that can leach out if the products are land filled. Transport- and dispenser boxes are designed to fit existing recovery systems. The Mölnlycke packaging system complies with the EU Directive on packaging and packaging waste.

Storage instructions

Mölnlycke Health Care recommends that Mölnlycke products are stored under normal storage conditions. All layers of packaging should be kept intact until access to the underlying layers is needed. Storage facilities for products only protected by the sterile barrier should be kept under conditions where low level of particulate air contamination prevail, so that it would not constitute a risk to the patient when the sterile barrier is opened and the product is used.

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2022 Mölnlycke Health Care AB. All rights reserved.



Classifications

| Regulation type(s) | MDR Class IS | Locally Regulated | Locally Regulated | AU Class IS |
|-------------------------------------|--|--|--|--|
| CE Certificate Number : | MDR 722028 | | | MDR 722028 |
| Intended Purpose : | The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure. | The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure. | The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure. | The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure. |
| MDR Classification Rule : | 1 | | | |
| Conformity Annexes : | IX | | | |
| Measuring Function : | Yes | | | |
| Notified body medical devices/PPE : | BSi 2797 | Local Authority India | Local Authority New Zealand | TGA (Australia) |

Applied standards : The standards presented below is a selection of the most essential standards that are adhered to.

EN 1041, EN 556-1, EN 13795, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 11607-1, EN ISO 11607-2, EN ISO 15223-1, ISO 15223-2, ISO 14001

Removable label

Yes

GMDN Code (Global Medical Device Nomenclature)

47783 Patient surgical drape, single-use

UNSPSC

42131701 Surgical drapes

Commodity Code

3005100000 Wadding, Gauze, Dressings, drapes singlepacked - adhesive articles, Sets mainly consisting thereof

CE Responsibility

Molnlycke Health Care AB, Gamlestadsvägen 3c, Box 13080, SE-402 52 Goteborg Sweden

Basic UDIDI

73324300000000160JX

Find out more at www.molnlycke.com

