

84541251 Baby Blanket


Product details

Size : 77cm x 77cm

Sterile : Sterile

Brand : BARRIER

Images



Delivered items

84541251-02

Country of origin: Thailand

Shelf life: 5 years

Sterilization method: Irradiation

Production Responsibility: Mölnlycke Health Care (Thailand) Ltd, 160 Moo 17 Bangplee Industrial Estate, Bangna-Trad Road, Km 23 Tambol Bangsaothong, Amphur Bangsaothong, Samutprakarn 10570, Thailand

Packing information: First packaging layer is a peel-open sterile barrier, plastic/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a corrugated board dispenser box. Third layer is a corrugated board transport box.

Is suitable for Tray: No

Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Piece	1	7332430081919			
Consumer pack	1	7323190039913			
Dispenser box	18	7323190039173			

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2022 Mölnlycke Health Care AB. All rights reserved.

Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Transport box	108	7323190039166	366x753x295 mm	81.3 dm3	8.4 / 6.4 kg
Pallet	1944	7323190039159	800x1200x1914 mm		

84541251-21

Country of origin: Thailand**Shelf life:** 5 years**Sterilization method:** Irradiation**Production Responsibility:** Mölnlycke Health Care (Thailand) Ltd, 160 Moo 17 Bangplee Industrial Estate, Bangna-Trad Road, Km 23 Tambol Bangsaothong, Amphur Bangsaothong, Samutprakarn 10570, Thailand**Packing information:** First packaging layer is a peel-open sterile barrier, plastic/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a corrugated board dispenser box. Third layer is a corrugated board transport box.**Is suitable for Tray:** No

Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Piece	1	7332430081919			
Consumer pack	1	7323190292424			
Dispenser box	20	7323190289561			
Transport box	80	7323190289554	365x560x268 mm	54.8 dm3	5.7 / 4.3 kg
Pallet	1920	7323190289547	800x1200x1751 mm	1681.0 dm3	161.1 / 103.7 kg

Material

Animal tissues :	No
Human blood derivatives :	No
Natural rubber latex :	No
Medicinal substances :	No
Phthalates :	No
Polyvinyl chloride :	No

Composition Drape Other

Product Component	Composition
Main Material	Viscose/Polyester nonwoven 71 g/m ²

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2022 Mölnlycke Health Care AB. All rights reserved.



Technical

Dimension

Dimension text	Dimension value
Outer dimension	77 cm x 77 cm

Disposal instructions

Used Mölnlycke products and sterile barriers should, in the majority of cases, be classified as non-hazardous waste. They contain high amounts of energy and are well suited for incineration. Mölnlycke products do not contain any hazardous substances that can leach out if the products are land filled. Transport- and dispenser boxes are designed to fit existing recovery systems. The Mölnlycke packaging system complies with the EU Directive on packaging and packaging waste.

Storage instructions

Mölnlycke Health Care recommends that Mölnlycke products are stored under normal storage conditions. All layers of packaging should be kept intact until access to the underlying layers is needed. Storage facilities for products only protected by the sterile barrier should be kept under conditions where low level of particulate air contamination prevail, so that it would not constitute a risk to the patient when the sterile barrier is opened and the product is used.

Classifications

Regulation type(s)	MDD Class IS	Locally Regulated	Locally Regulated	Locally Regulated
CE Certificate Number :	CE 01966			
Intended use MDD :	Intended use is present in the classification assessment connected to this classification object			
Intended Purpose :		Intended use is present in the classification assessment connected to this classification object	Intended use is present in the classification assessment connected to this classification object	Intended use is present in the classification assessment connected to this classification object
Conformity Annexes :	V;VII			
Measuring Function :	No			
Notified body medical devices/PPE :	BSi 2797	Local Authority United Arab Emirates	Local Authority APAC	Local Authority Americas

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2022 Mölnlycke Health Care AB. All rights reserved.



Applied standards : The standards presented below is a selection of the most essential standards that are adhered to.

EN 1041, EN 556-1, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 11607-1, EN ISO 11607-2, EN ISO 15223-1, ISO 15223-2, ISO 14001

Removable label

Yes

GMDN Code (Global Medical Device Nomenclature)

14077 Surgical towel

UNSPSC

42131701 Surgical drapes

Commodity Code

3005905000 Articles of Wadding, Gauze, Dressings, drapes singlepacked and Sets - textile

CE Responsibility

Mölnlycke Health Care AB, Gamlestadsvägen 3c, Box 13080, SE-402 52 Goteborg Sweden

Basic UDIDI

733243000000000065K4

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2022 Mölnlycke Health Care AB. All rights reserved.

