


906542 Adhesive Aperture Drape

Product details

Size :	75cm x 90cm, Aperture 6x8cm
Descriptive feature :	Adhesive, Centered ellipse aperture
Sterile :	Sterile
Brand :	BARRIER
Equivalences :	904758

Images



Delivered items

906542-60

Sales released in: Austria, Bahrain, Belgium, Bosnia and Herzegovina, Chile, Colombia, Croatia, Cyprus, Czechia, Denmark, Estonia, Finland, France, Germany, Greece, Hong Kong, Hungary, Iceland, Ireland, Italy, Japan, Latvia, Lebanon, Lithuania, Luxembourg, Mexico, Moldova (the Republic of), Montenegro, Morocco, Netherlands, Norway, Oman, Panama, Poland, Portugal, Qatar, Russian Federation, Saudi Arabia, Serbia, Singapore, Slovakia, Slovenia, Spain, Sweden, Switzerland, Syrian Arab Republic, Taiwan (Province of China), Thailand, United Arab Emirates, United Kingdom of Great Britain and Northern Ireland

Shelf life: 5 years

Sterilization method: Irradiation

Production Responsibility: Mölnlycke Health Care (Thailand) Ltd, Amata Nakorn (Bang Pakong), Industrial Estate, 700/461 Moo 7 Bangna-Trad Road, Km 57 Tambol Donhuaroh, Amphur Muang Chonburi 20000, Thailand

Packing information: First packaging layer is a peel-open sterile barrier, paper/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a corrugated board dispenser box. Third layer is a corrugated board transport box.

Is suitable for Tray: No

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2022 Mölnlycke Health Care AB. All rights reserved.



Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Piece	1	7323190258352			
Consumer pack	1	7323190272556			
Dispenser box	40	7323190262205			
Transport box	80	7323190262199	332x428x288 mm	40.9 dm3	4.5 / 3.6 kg
Pallet	2880	7323190262182	800x1200x1872 mm	1797.1 dm3	188.1 / 129.6 kg

906542-20

Sales released in: Algeria, Austria, Bahrain, Belgium, Bosnia and Herzegovina, Chile, Colombia, Croatia, Cyprus, Czechia, Denmark, Estonia, Finland, France, Germany, Greece, Hong Kong, Hungary, Iceland, Ireland, Italy, Japan, Latvia, Lebanon, Lithuania, Luxembourg, Mexico, Moldova (the Republic of), Montenegro, Morocco, Netherlands, Norway, Oman, Panama, Poland, Portugal, Qatar, Russian Federation, Saudi Arabia, Serbia, Singapore, Slovakia, Slovenia, Spain, Sweden, Switzerland, Syrian Arab Republic, Taiwan (Province of China), Thailand, United Arab Emirates, United Kingdom of Great Britain and Northern Ireland

Shelf life: 5 years

Sterilization method: Irradiation

Production Responsibility: Mölnlycke Health Care (Thailand) Ltd, Amata Nakorn (Bang Pakong), Industrial Estate, 700/461 Moo 7 Bangna-Trad Road, Km 57 Tambol Donhuaroh, Amphur Muang Chonburi 20000, Thailand

Packing information: First packaging layer is a peel-open sterile barrier, paper/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a corrugated board dispenser box. Third layer is a corrugated board transport box.

Is suitable for Tray: No

Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Piece	1	7323190258352			
Consumer pack	1	7323190179718			
Dispenser box	50	7323190178575			
Transport box	200	7323190178568	398x594x444 mm	105.0 dm3	11.0 / 7.2 kg
Pallet	3200	7323190178551	800x1200x1920 mm		

Technical

Dimension

Dimension text	Dimension value
Outer dimension	75 cm x 90 cm

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2022 Mölnlycke Health Care AB. All rights reserved.



Dimension text	Dimension value
Aperture	6 cm x 8 cm

Disposal instructions

Used Mölnlycke products and sterile barriers should, in the majority of cases, be classified as non-hazardous waste. They contain high amounts of energy and are well suited for incineration. Mölnlycke products do not contain any hazardous substances that can leach out if the products are land filled. Transport- and dispenser boxes are designed to fit existing recovery systems. The Mölnlycke packaging system complies with the EU Directive on packaging and packaging waste.

Storage instructions

Mölnlycke Health Care recommends that Mölnlycke products are stored under normal storage conditions. All layers of packaging should be kept intact until access to the underlying layers is needed. Storage facilities for products only protected by the sterile barrier should be kept under conditions where low level of particulate air contamination prevail, so that it would not constitute a risk to the patient when the sterile barrier is opened and the product is used.

Classifications

Regulation type(s)	MDR Class IS	Locally Regulated	Locally Regulated	Locally Regulated
CE Certificate Number :	MDR 722028			
Intended Purpose :	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.
MDR Classification Rule :	1			
Conformity Annexes :	IX			
Measuring Function :	Yes			
Notified body medical devices/PPE :	BSi 2797	Local Authority Americas	Local Authority Mexico	Local Authority Thailand
510(k) clearance number :			203300CT070363	

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
 The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
 ©2022 Mölnlycke Health Care AB. All rights reserved.



Classifications

Regulation type(s)	Locally Regulated	Class I
CE Certificate Number :		
Intended Purpose :	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.
MDR Classification Rule :		
Conformity Annexes :		
Measuring Function :		
Notified body medical devices/PPE :	Local Authority Hong Kong	PMDA (Japan)
510(k) clearance number :		

Applied standards : The standards presented below is a selection of the most essential standards that are adhered to.

EN 1041, EN 556-1, EN 13795, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 11607-1, EN ISO 11607-2, EN ISO 15223-1, ISO 15223-2, ISO 14001

Removable label

Yes

GMDN Code (Global Medical Device Nomenclature)

47783 Patient surgical drape, single-use

UNSPSC

42131701 Surgical drapes

Commodity Code

3005100000 Wadding, Gauze, Dressings, drapes singlepacked - adhesive articles, Sets mainly consisting thereof

CE Responsibility

Mölnlycke Health Care AB, Gamlestadsvägen 3c, Box 13080, SE-402 52 Goteborg Sweden

Basic UDIDI

733243000000000160JX

Find out more at www.molnlycke.com

