


920496 Adhesive OP-Sheet

Product details

Size :	175cm x 175cm
Descriptive feature :	3Ply, Adhesive
Sterile :	Sterile
Brand :	BARRIER
Equivalences :	967775

Images



Delivered items

920496-21

Sales released in: Australia, Austria, Belgium, Bosnia and Herzegovina, Croatia, Cyprus, Czechia, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, India, Ireland, Israel, Italy, Japan, Kuwait, Latvia, Lithuania, Luxembourg, Malta, Netherlands, New Zealand, Norway, Poland, Qatar, Saudi Arabia, Serbia, Slovakia, Slovenia, Sweden, Switzerland, Thailand, United Arab Emirates, United Kingdom of Great Britain and Northern Ireland

Shelf life: 5 years

Sterilization method: Irradiation

Production Responsibility: Mölnlycke Health Care (Thailand) Ltd, Amata Nakorn (Bang Pakong), Industrial Estate, 700/461 Moo 7 Bangna-Trad Road, Km 57 Tambol Donhuaroh, Amphur Muang Chonburi 20000, Thailand

Packing information: First packaging layer is a peel-open sterile barrier, plastic/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a plastic liner/bag. Third layer is a corrugated board dispenser box. Fourth layer is a corrugated board transport box.

Is suitable for Tray: No

Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Piece	1	7323190258468			

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2022 Mölnlycke Health Care AB. All rights reserved.



Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Consumer pack	1	7323190294275			
Dispenser box	18	7323190294084			
Transport box	36	7323190294077	397x580x329 mm	75.8 dm3	10.2 / 9.2 kg
Pallet	720	7323190294060	800x1200x1789 mm	1717.4 dm3	229.4 / 180.0 kg

Technical

Dimension

Dimension text	Dimension value
Outer dimension	175 cm x 175 cm

Disposal instructions

Used Mölnlycke products and sterile barriers should, in the majority of cases, be classified as non-hazardous waste. They contain high amounts of energy and are well suited for incineration. Mölnlycke products do not contain any hazardous substances that can leach out if the products are land filled. Transport- and dispenser boxes are designed to fit existing recovery systems. The Mölnlycke packaging system complies with the EU Directive on packaging and packaging waste.

Storage instructions

Mölnlycke Health Care recommends that Mölnlycke products are stored under normal storage conditions. All layers of packaging should be kept intact until access to the underlying layers is needed. Storage facilities for products only protected by the sterile barrier should be kept under conditions where low level of particulate air contamination prevail, so that it would not constitute a risk to the patient when the sterile barrier is opened and the product is used.

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2022 Mölnlycke Health Care AB. All rights reserved.



Classifications

Regulation type(s)	MDR Class IS	MDR Class IS	MDR Class IS	MDR Class IS
CE Certificate Number :	MDR 722028			
Intended Purpose :	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.
MDR Classification Rule :	1			
Conformity Annexes :	IX			
Measuring Function :	No			
Notified body medical devices/PPE :	BSi 2797	BSi 2797	BSi 2797	BSi 2797

Classifications

Regulation type(s)	Locally Regulated	Locally Regulated	Locally Regulated	AU Class IS
CE Certificate Number :				CE 01966
Intended Purpose :	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.
MDR Classification Rule :				
Conformity Annexes :				
Measuring Function :				
Notified body medical devices/PPE :	Local Authority India	Local Authority Thailand	Local Authority APAC	TGA (Australia)

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
 The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
 ©2022 Mölnlycke Health Care AB. All rights reserved.



Classifications

Regulation type(s)	Class I
CE Certificate Number :	
Intended Purpose :	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.
MDR Classification Rule :	
Conformity Annexes :	
Measuring Function :	
Notified body medical devices/PPE :	PMDA (Japan)

Applied standards : The standards presented below is a selection of the most essential standards that are adhered to.

EN 1041, EN 556-1, EN 13795, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 11607-1, EN ISO 11607-2, EN ISO 15223-1, ISO 15223-2, ISO 14001

Removable label

Yes

GMDN Code (Global Medical Device Nomenclature)

47783 Patient surgical drape, single-use

UNSPSC

42132103 Patient barrier drapes

Commodity Code

3005905000 Articles of Wadding, Gauze, Dressings, drapes singlepacked and Sets - textile

CE Responsibility

Mölnlycke Health Care AB, Gamlestadsvägen 3c, Box 13080, SE-402 52 Goteborg Sweden

Basic UDIDI

73324300000000177KG

Find out more at www.molnlycke.com

