

evercare® Medical face mask Duckbill with tie-bands, type II and IIR

REF: 2852

evercare[®] Medical face mask Duckbill with tie-bands is a type II and IIR face mask that provides a barrier to minimize the direct transmission of infective agents between medical staff and patient. The upper part of the face mask is type IIR and the lower part is type II. Pleasantly off-the face and splash resistant. Recommended for all types of surgery.

Key Features

- ► EN 14683:2019
- ▶ Bacterial Filtration Efficiency (BFE) ≥ 98%
- ▶ Single use
- ▶ Flexible nose wire for proper fit







REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2852	Medical face mask Duckbill with tie-bands, type IIR, non-sterile, Blue	187 x 107 mm	50 / 300

Sales Unit: Transport Carton

Page 1 2022-11-30

Technical Summary

Regulations

The product complies with legal requirements for medical devices.















REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).		
Biological evaluation:	EN ISO 10993-05 EN ISO 10993-10		
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class I		
Compliance with product standards:	EN 14683:2019, Type II and Type IIR, Bacterial Filtration Efficiency (BFE) ≥98%		
Quality standards:	EN ISO 13485		
Label information and packaging compliance:	EN ISO 15223-1, EN 1041		

Material

Product Part	Raw Material & Characteristics	
Outer facing, upper part	Spunbond PP, blue	
Outer facing, lower part	Meltblown PP, white	
Filter media, upper part	Meltblown PP, white	
Filter media, lower part	Meltblown PP, white	
Inner facing	Meltblown PP, white	
Ties	Spunbond PP, white, upper/lower length 430/355 mm	
Nose wire	PE coated galvanize wire, 150 mm	

2022-11-30 Page 2

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.	
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.	
Disposal Recommendation	Dispose in accordance with local regulations.	
Shelf Life	3 Years	

Packing information

Package	Material	Size (L x W x H) mm	
Single Pack	N/A	N/A	
Dispenser Pack	Coated paper	185 x 110 x 130	
Transport Carton	Carton, TCF or ECF quality	345 x 197 x 280	
Packing Tape	PVC-free, Transparent	N/A	

Barcodes

REF	Single Pack	Dispenser	Transport Carton
2852	N/A	6438129315615	6438129515619

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.

Page 3