

evercare® Medical face mask with tie-bands, type IIR

REF: 2853

evercare® Medical face mask with tie-bands, type IIR provides a barrier to minimize the direct transmission of infective agents between medical staff and patient. It is a splash resistant mask offering protection against potentially contaminated fluid splashes. Recommended to be used during wet surgery.

Key Features

- ▶ EN 14683:2019
- ▶ Bacterial Filtration Efficiency (BFE) ≥ 98%
- ▶ Single use
- ▶ Flexible nose wire for proper fit
- ▶ Splash resistant
- ▶ Anti-fog
- ▶ Soft inner layer for comfort



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2853	Medical face mask with tie-bands, anti-fog, extra wide, type IIR, non-sterile, White	190 x 85 mm	50 / 300

Sales Unit: Transport Carton

Technical Summary

Regulations

The product complies with legal requirements for medical devices.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-05 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class I
Compliance with product standards:	EN 14683:2019, Type IIR, Bacterial Filtration Efficiency (BFE) ≥98%
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN 1041

Material

Product Part	Raw Material & Characteristics
1st outer facing	Spunbond PP, white
2nd outer facing	Spunbond PP, white
Filter	Meltblown PP, white
Inner facing	Hydrophobic tissue, white
Ties and side bindings	Spunbond PP, white, left/right length 405 mm
Nose wire	Aluminium wire, 133 mm
Anti-fog outside and inside	Spunbond PP, white + Coated LDPE, opaque white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	3 Years

Packing information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	Coated paper	203 x 127 x 184
Transport Carton	Carton, TCF or ECF quality	395 x 215 x 385
Packing Tape	PVC-free, Transparent	N/A

Barcodes

REF	Single Pack	Dispenser	Transport Carton
2853	N/A	6438129315622	6438129515626

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.