

evercare® Medical face mask Multifix, standard, with earloops, type II

REF: 2856

evercare[®] Medical face mask Multifix, standard with earloops, type II provides a barrier to minimize the direct transmission of infective agents between medical staff and patient. It is an allround face mask with high filtration capacity and low breathing resistance. Recommended to be used during dry surgery, with no risk of splashes.

Key Features

- PEN 14683:2019
- ▶ Bacterial Filtration Efficiency (BFE) ≥ 98%
- ▶ Single use
- ▶ Flexible nose wire for proper fit
- ▶ Color: Outer facing green, inner facing white





REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2856	Medical face mask Multifix Standard with earloops, type II, non-sterile, Green	175 x 90 mm	50 / 300

Sales Unit: Transport Carton

Page 1 2022-11-30

Technical Summary

Regulations

The product complies with legal requirements for medical devices.













EN ISO 15223-1, EN 1041



as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of **REACH** regulation: the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA). EN ISO 10993-05 **Biological evaluation:** EN ISO 10993-10 Compliance with legal Medical Device Regulation (EU) 2017/745, MDR class I requirements: Compliance with product EN 14683:2019, Type II, Bacterial Filtration Efficiency (BFE) ≥98% standards: **Quality standards:** EN ISO 13485

The products covered by this data sheet do not contain substances classified

Material

Label information and

packaging compliance:

Product Part	Raw Material & Characteristics	
Outer facing	Spunbond PP, green	
Filter media	Meltblown PP, white	
Inner facing	Thermalbonded PP, white	
Bindings	Spunbond PP, white	
Elastic earloops	Latex-free, white	
Nose wire	PE coated galvanize wire, 125 mm	

Page 2 2022-11-30

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.	
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.	
Disposal Recommendation	Dispose in accordance with local regulations.	
Shelf Life	5 Years	

Packing information

Package	Material	Size (L x W x H) mm	
Single Pack	N/A	N/A	
Dispenser Pack	Coated paper	185 x 110 x 80	
Transport Carton	Carton, TCF or ECF quality	213 x 197 x 260	
Packing Tape	PVC-free, Transparent	N/A	

Barcodes

REF	Single Pack	Dispenser	Transport Carton
2856	N/A	6438129312935	6438129512939

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.

Page 3