



evercare® Medical face mask Multifix Splash, extra, with earloops, type IIR

REF: 2875

evercare® Medical face mask Multifix Splash, extra with earloops type IIR provides a barrier to minimize the direct transmission of infective agents between medical staff and patient. It is a splash resistant mask offering protection against potentially contaminated fluid splashes. Recommended to be used during wet surgery.

Key Features

- ▶ EN 14683:2019
- ▶ Bacterial Filtration Efficiency (BFE) $\geq 98\%$
- ▶ Single use
- ▶ Flexible nose wire for proper fit
- ▶ Splash resistant
- ▶ Color: Outer facing green, inner facing white










REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2875	Medical face mask Multifix Splash Extra with earloops, type IIR, non-sterile, Green	175 x 90 mm	50 / 300

Sales Unit: Transport Carton

Technical Summary

Regulations

The product complies with legal requirements for medical devices.

REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-05 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class I
Compliance with product standards:	EN 14683:2019, Type IIR, Bacterial Filtration Efficiency (BFE) ≥98%
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN 1041

Material

Product Part	Raw Material & Characteristics
Outer facing	Spunbond PP, green
Filter media	Meltblown PP, white
Inner facing	Spunbond PP, white
Bindings	Spunbond PP, white
Elastic earloops	Latex-free, white
Nose wire	PE coated galvanize wire, 125 mm

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 Years

Packing information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	Coated paper	185 x 110 x 130
Transport Carton	Carton, TCF or ECF quality	345 x 197 x 280
Packing Tape	PVC-free, Transparent	N/A

Barcodes

REF	Single Pack	Dispenser	Transport Carton
2875	N/A	6438129312980	6438129512984

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.