

## evercare<sup>®</sup> Medical face mask with earloops, type IIR

REF: 2876

evercare<sup>®</sup> Medical face mask, white with earloops type IIR provides a barrier to minimize the direct transmission of infective agents between medical staff and patient. It is a splash resistant mask offering protection against potentially contaminated fluid splashes. Recommended to be used during wet surgery.

### Key Features

- ▶ EN 14683:2019
- ▶ Bacterial Filtration Efficiency (BFE)  $\geq 98\%$
- ▶ Single use
- ▶ Flexible nose wire for proper fit
- ▶ Splash resistant
- ▶ Anti-fog
- ▶ Soft inner layer for comfort



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2876	Medical face mask with earloops, anti-fog, type IIR, non-sterile, White	175 x 95 mm	50 / 300

Sales Unit: Transport Carton

# Technical Summary

## Regulations

The product complies with legal requirements for medical devices.

**REACH regulation:**

The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

**Biological evaluation:**

EN ISO 10993-05  
EN ISO 10993-10

**Compliance with legal requirements:**

Medical Device Regulation (EU) 2017/745, MDR class I

**Compliance with product standards:**

EN 14683:2019, Type IIR, Bacterial Filtration Efficiency (BFE) ≥98%

**Quality standards:**

EN ISO 13485

**Label information and packaging compliance:**

EN ISO 15223-1, EN 1041

## Material

Product Part	Raw Material & Characteristics
1st outer facing	Spunbond PP, white
2nd outer facing	Spunbond PP, white
Filter	Meltblown PP, white
Inner facing	Hydrophobic tissue, white
Bindings	Spunbond PP, white
Nose wire	PE coated galvanize wire, 125 mm
Elastic earloops	Latex-free, white
Anti-fog outside and inside	Spunbond PP, white + Coated PE, opaque white

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	3 Years

## Packing information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	N/A	N/A
<b>Dispenser Pack</b>	Coated paper	185 x 100 x 80
<b>Transport Carton</b>	Carton, TCF or ECF quality	213 x 197 x 260
<b>Packing Tape</b>	PVC-free, Transparent	N/A

## Barcodes

REF	Single Pack	Dispenser	Transport Carton
<b>2876</b>	N/A	6438129315646	6438129515640

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.