



evercare[®] Essentia medical face mask with tie-band, type II

REF: 2877

evercare[®] Essentia Medical face mask with tie-bands type II provides a barrier to minimize the direct transmission of infective agents between medical staff and patient. It is an allround face mask with high filtration capacity and low breathing resistance. Recommended to be used during dry surgery, with no risk of splashes.

Key Features

- ▶ EN 14683:2019
- ▶ Bacterial Filtration Efficiency (BFE) ≥ 98%
- ▶ Single use
- ▶ Flexible nose wire for proper fit
- ▶ Color: Outer facing blue, inner facing white
- ▶ Soft inner layer for comfort
- ▶ Each mask separately folded; easy to take out one by one



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2877	Essentia medical face mask with tie-band, non-sterile, Blue	175 x 90 mm	50 / 900

Sales Unit: Transport Carton

Technical Summary

Regulations

The product complies with legal requirements for medical devices.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-05 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class I
Compliance with product standards:	EN 14683:2019, Type II, Bacterial Filtration Efficiency (BFE) ≥98%
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN 1041

Material

Product Part	Raw Material & Characteristics
Outer facing	Spunbond PP nonwoven, blue
Filter media	Meltblown PP nonwoven, white
Inner facing	Spunbond PP nonwoven, white
Tie-bands	Spunbond PP nonwoven, white, upper/lower length 400 mm
Nose wire	PE coated ferrum/titanium wire, 110 mm

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 Years

Packing information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	Coated paper	185 x 100 x 80
Transport Carton	Carton, TCF or ECF quality	213 x 197 x 260
Packing Tape	PVC-free, Transparent	N/A

Barcodes

REF	Single Pack	Dispenser	Transport Carton
2877	N/A	6438129322316	6438129522310

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.