

evercare[®] Medikit Catheterization set, sterile REF: 400934

evercare[®] MediKit Single-use procedure small kits are assembled to provide healthcare professionals with the necessary supplies for specific medical interventions. The primary purpose of catheterization sets is to facilitate urological catheterization. Available in a single, sterile package, ensuring optimal hygiene conditions and patient safety.

Key Features

- 1 Utility drape 70x90cm
- 1 Utility drape 75x90cm slitted
- 1 Kidney-shaped basin 700ml
- > 7 Nonwoven swab 10x10cm 4ply 30g
- ▶ 1 Plastic tweezer forceps with double grooves 12,5cm
- 1 Bowl 120ml graduated






10 × 10 cm, 4 ply 30 g



1 Kidney-shaped basin 700 ml

1 Plastic tweezer forceps with double grooves 12.5 cm

Πn



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
400934	Cath. Set		14 / 28

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-7
Compliance with legal requirements:	Art.22/23 of Medical Device Regulation (EU) 2017/745, MDR Class IIa
Compliance with product standards:	EN ISO 11135
Label information and packaging compliance:	EN ISO 20417, EN ISO 11607-1, EN ISO 15223-1

Material

Product Part	Raw Material & Characteristics
1 Utility drape 70x90cm	SPP/PE
1 Utility drape 75x90cm slitted	SPP/PE
1 Kidney-shaped basin 700ml	PE blue
7 Nonwoven swab 10x10cm 4ply 30g	VIS+PES white
1 Plastic tweezer forceps with double grooves 12,5cm	PP blue
1 Bowl 120ml graduated	PE translucent

Warnings & Recommendations

Warnings	VarningsThe product is single use. If reused, the performance of product may deteriorate and cross contamination may occur. Do not use if package is damaged.		
Storage Recommendation		Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.	
Disposal Recommendation		Dispose in accordance with local regulations.	
Shelf Life		5 years	

Packaging Information

Package	Material	Size(L x W x H)mm
Single Pack	Medical grade paper	420 x 250 x 10
Dispenser Pack	Cardboard	370 x 282 x 305
Transport Carton	TCF or ECF quality	597 x 405 x 375
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser Pack	Transport Carton
400934	03700433504683	03700433504676	03700433504669

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.