

# Terumo Syringe

3-part Syringe



## Product Information

NEW PRODUCT CODES!



### Disposable syringes for single use

Every single detail is the result of Terumo's long time experience in this field which is appreciated by users throughout the world.

- Siliconisation for optimal plunger mobility
- The integral gasket ensures a firm and distinct end-stop
- Clear graduated scale allows accurate dosage
- Transparent barrel for visualisation
- Finger grips to support easy handling

Description	Raw Materials
1 Barrel	Polypropylene
2 Plunger	Polypropylene
3 Gasket	Thermoplastic elastomer (for ,SS'-codes) Polyisoprene rubber (for ,MDSS'-codes)
Lubricant	Silicone oil

Does not contain components made of PVC or natural rubber latex.

#### Sterilization

With radiation (e-beam) (for ,SS'-codes) or ethylene oxide (for ,MDSS'-codes) to a Sterility Assurance Level of at least  $10^{-6}$

#### Storage

- Shelf life of 5 years
- The batch number and expiry date are printed on the paper (for ,MDSS'-codes) or embossed in the film of the blister pack (for ,SS'-codes)
- Do not store at extreme temperature and humidity

#### Precautions

- Do not use if package is damaged
- Do not re-use
- Dispose of sharps after single use according to facility protocol to avoid risk of infection

Always consult the product label for a complete overview of warnings, cautions and/or precautions prior to actual use.



## Your syringes of choice NEW PRODUCT CODES! in different sizes and tip types

### Product Range

Product Code	Nominal Capacity	Tip	Replacement for
<b>MDSS01SE</b>	1 ml	Concentric Luer Slip	SS+01T1
SS*02LE1	2,5 ml	Concentric Luer Lock	-
SS*02SE1	2,5 ml	Concentric Luer Slip	-
<b>MDSS03SE</b>	3 ml	Concentric Luer Slip	SS*02SE1, SS+03S
<b>MDSS03LE</b>	3 ml	Concentric Luer Lock	SS*03L, SS+03L1
<b>MDSS05SE</b>	5 ml	Concentric Luer Slip	SS*05SE1, SS+05S
SS*05LE1	5 ml	Concentric Luer Lock	-
SS*05SE1	5 ml	Concentric Luer Slip	-
<b>MDSS10SE</b>	10 ml	Concentric Luer Slip	SS+10ES1
<b>MDSS10ESE</b>	10 ml	Eccentric Luer Slip	SS+10ES1
SS*10LE1	10 ml	Concentric Luer Lock	-
<b>MDSS20ESE</b>	20 ml	Eccentric Luer Slip	SS+20ES1
SS+20L1	20 ml	Concentric Luer Lock	-
SS+30L1MP	30 ml	Concentric Luer Lock	-
SS+30ES1	30 ml	Eccentric Luer Slip	-
SS+50L1	50 ml	Concentric Luer Lock	-
SS+50ES1	50 ml	Eccentric Luer Slip	-
SS+50C1	50 ml	Catheter Tip	-

### Packaging

Unit Box	Quantity	Shipping Carton	
		Gross Weight	Dimensions (L x W x H)
100	3000	16.60 kg	635 x 295 x 505 mm
100	2400	6.00 kg	397 x 581 x 444 mm
100	2400	6.00 kg	394 x 581 x 444 mm
100	1800	10.90 kg	445 x 447 x 449 mm
100	1800	10.90 kg	445 x 447 x 449 mm
100	1800	14.90 kg	656 x 419 x 426 mm
100	1600	6.25 kg	412 x 573 x 444 mm
100	1600	6.25 kg	411 x 576 x 444 mm
100	1200	14.60 kg	569 x 473 x 459 mm
100	1200	14.60 kg	565 x 425 x 485 mm
100	800	15.00 kg	359 x 581 x 383 mm
50	400	14.60 kg	468 x 390 x 455 mm
50	600	10.15 kg	573 x 410 x 377 mm
25	200	6.32 kg	560 x 398 x 255 mm
25	200	6.20 kg	558 x 395 x 254 mm
25	100	5.10 kg	464 x 320 x 275 mm
25	100	5.10 kg	472 x 324 x 282 mm
25	100	5.40 kg	472 x 324 x 282 mm



Single unit: Blister pack  
(film/film for ,SS'-codes and film/paper for ,MDSS'-codes)

### Regulatory Status

CE Certification – Class IIa medical device

**Notified body** TÜV Rheinland 0197

**Manufacturer**   
 for ,SS+'-codes TERUMO (Philippines) CORPORATION,   
 124 East Main Avenue, Laguna Technopark, Binan, Philippines   
 for ,SS\*'-codes TERUMO CORPORATION, 2-44-1 Hatagaya,   
 Shibuya-Ku, Toyko 151-0072, Japan   
 for ,MDSS'-codes Yangzhou Medline Industry Co., Ltd. (China)

**EC representative**   
 for ,SS'-codes TERUMO EUROPE N.V. – Leuven, Belgium   
 for ,MDSS'-codes Shanghai International Corp. GmbH (Europe) – Hamburg, Germany

**Distributor** TERUMO EUROPE N.V., Interleuvenlaan 40, 3001 Leuven, Belgium

Choose from a range of volumes from 1 ml to 50 ml

Syringes are available in all common tip types.



**Terumo Corporation**  
+81 3 6742 8500

**Terumo Europe N.V.**  
+32 16 38 12 11

Your Terumo partner

All brand names are trademarks or registered trademarks of TERUMO CORPORATION and their respective owners. Published by Terumo Europe N.V.

Please note that some products may not be available in all regions, check with your local Terumo contact.