


## HOSPITAL CARE

# Versatus - W

## Winged I.V. Catheter Product Information



- The flashback chamber visually confirms that the catheter is correctly inserted in the vein
- The sharp needle tip minimizes trauma to the vessel during use
- Wide range from 14G up to 26G










Description		Raw Materials
1	Catheter tube	Polyurethane
2	Winged Catheter hub Caulking pin	Polypropylene + master batch Polyoxymethylene
3	Cannula (internal needle)	Stainless steel
4	Needle hub	Styrene butadiene copolymer
5	Luer-lock cap	Polypropylene
6	Vent fitting filter Vent fitting body	Polyethylene Styrene butadiene copolymer
	Protector	Polyethylene
	Unit pack	Adhesive coated lid paper. Thermoformed blister film: laminated polyethylene/ polypropylene complex

Does not contain components made of PVC or natural rubber latex.

Designed and  
manufactured by skilled  
technicians in Europe

## Your catheters of choice for a wide range of I.V. applications

### Product Range

Product Code	Colour	External Catheter Ø		Internal Catheter Ø	Catheter Length	Flow Rate ml/min.
SR+DU2619WD	Violet 	26G	0.64 mm	0.45 mm	19 mm	19
SR+DU2419WD	Yellow 	24G	0.74 mm	0.55 mm	19 mm	29
SR+DU2225WD	Blue 	22G	0.90 mm	0.65 mm	25 mm	42
SR+DU2032WD	Pink 	20G	1.00 mm	0.75 mm	32 mm	59
SR+DU1832WD	Green 	18G	1.30 mm	0.95 mm	32 mm	103
SR+DU1845WD	Green 	18G	1.30 mm	0.95 mm	45 mm	96
SR+DU1745WD	White 	17G	1.50 mm	1.15 mm	45 mm	155
SR+DU1645WD	Grey 	16G	1.75 mm	1.35 mm	45 mm	225
SR+DU1445WD	Orange 	14G	2.00 mm	1.55 mm	45 mm	290



### Packaging

**Single Unit** Blister pack

**Unit Box** 50 pieces

**Shipping Carton** 500 pieces, Gross Weight ± 4 kg  
Dimensions ± 385 x 779 x 170 mm

### Regulatory Status

CE Certification – Class IIa medical device

#### Notified body

TÜV SÜD Product Service GmbH (0123)

#### Manufacturer

Delta Med S.p.A. (Italy)

#### Distributor

TERUMO EUROPE N.V., Interleuvenlaan 40, 3001 Leuven, Belgium

### Sterilization

With ethylene oxide to a Sterility Assurance Level of at least 10<sup>-6</sup>

### Storage

- Shelf life of 5 years
- The batch number and expiry date are printed on the paper of the unit pack
- Keep away from direct sunlight and heat.
- Keep dry
- Fragile, handle with care

### Precautions

- Do not use if the package is damaged and consult instructions for use
- Do not re-use / do not re-sterilize
- Do not attempt to re-insert a partially or completely withdrawn needle
- Dispose of safely after single use to avoid risk of infection
- Application time should be consistent with accepted standards of practice

Always consult the product label and the instructions for use for a complete overview of warnings, cautions and/or precautions prior to actual use.

This range complements the wider Terumo product offering including a comprehensive portfolio of both safety and conventional I.V. products which can be found at our web page on [www.terumo-europe.com](http://www.terumo-europe.com)

Your Terumo partner



#### TERUMO EUROPE N.V.

Terumo Corporation Terumo Europe N.V.  
+81 3 6742 8500 +32 16 38 12 11

[www.terumo-europe.com](http://www.terumo-europe.com)

All brand names are trademarks or registered trademarks owned by TERUMO CORPORATION, its affiliates, or unrelated third parties. Published by Terumo Europe N.V.

Image for visualisation purpose only