


# evercare® Angiography set

REF: 1239-01

evercare R Angiography set, sterile is a set for draping patients and operating units in angiographic procedures.

## **Key Features**

- Absorbent patch
- Adhesive OP Tape
- Main drape complies to EN 13795-1 High performance critical product area requirements



REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
1239-01	Angiography set, sterile	See Material table	10 / 10

Sales Unit: Transport Carton

Page 1 2023-12-05

### Regulations


The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.























**REACH regulation:** 

The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

EN ISO 10993-01 EN ISO 10993-5 **Biological evaluation:** EN ISO 10993-7

EN ISO 10993-10

Compliance with legal requirements:

Medical Device Regulation (EU) 2017/745, MDR class Is

Compliance with product standards:

EN13795-1

**Quality standards:** 

EN ISO 13485

Label information and packaging compliance:

EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

#### **Material**

Product Part	Raw Material & Characteristics		
1 Instrument table cover	100 x 150 cm, Laminate: PP nonwoven/PE film, blue		
1 Absorbent sheet,	56 x 75 cm, Spunlaced nonwoven of polyester/pulp viscose, white		
2 Hand towel	30 x 34 cm, Pulp, natural white, thread reinforced		
1 Adhesive drape sheet	50 x 60 cm, Laminate: PP nonwoven/PE film, blue		
1 Angiography drape	230 x 300 cm, Laminate: blue PP nonwoven/transparent PE film, absotbent patch 75 x 70 cm: PP nonwoven		
1 Instrument table cover, unfolded	150 x 190 cm, Laminate: PP nonwoven/PE film, blue		

Page 2 2023-12-05

## Warnings & Recommendations

Warnings contain	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.  Do not use if package is damaged.		
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.		
Disposal Recommendation	Dispose in accordance with local regulations.		
Shelf Life	5 years		

### **Packing Information**

Package	Material	Size ( L x W x H ) mm
Single Pack	Medical grade paper, PE/PET film	420 x 603
Dispenser Pack	TCF or ECF quality	485 x 285 x 287
Transport Carton	TCF or ECF quality	500 x 310 x 590
Packing Tape	Transparent	N/A

#### Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
1239-01	06438129110210	06438129323634	06438129523638

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.