

NasoGastric Tube Holders



The Dale NasoGastric Tube Holders are ideal for holding nasogastric feeding and aspiration tubes such as: Salem Sump™ tubes, Levin stomach tubes, nasoenteric feeding tubes and naso-gastric intestinal tubes.

The holder is an improvement over tape and other types of adhesive backed holders, because it features a skin-friendly adhesive that adheres to the skin of the nose quickly and comfortably. The holder can remain in place for up to 72 hours. The holder centers the tube exiting either from the left or right nostril, allowing the tube to “hang from” the nose in a more natural position.



Key Features

- ◆ Easy to use
- ◆ Securely holds all nasogastric tubes in place
- ◆ Prevents nasal irritation, erosion or necrosis
- ◆ Helps prevent inadvertent dislocation or removal of tube
- ◆ Provides maximum tube security
- ◆ Nose pad stretches and conforms to various nose shapes
- ◆ Non-adherent tips allow holder to be removed easily without sticking to exam gloves
- ◆ Not made with natural rubber latex

Product #	Description	Quantity
160	NasoGastric Tube Holder Large/Adult	50/box
161	NasoGastric Tube Holder Small/Young Adult/Pediatric/Geriatric	50/box

NasoGastric Tube Holders

Features	Clinical Advantages
Skin Friendly Adhesive	<ul style="list-style-type: none">▪ Provides maximum tube security▪ Securely holds feeding and/or aspiration tubes in place▪ Gentle on delicate skin▪ Resists humidified air of oxygen mask▪ Leaves no sticky residue on tube or nose after removal
Woven Fabric	<ul style="list-style-type: none">▪ Breathable material can be left in place for up to 3 days▪ Material stretches to conform to any shape nose▪ More comfortable than traditional tape
Dual Tabs	<ul style="list-style-type: none">▪ Maximum tube security▪ Prevents inadvertent dislocation or removal of the tube▪ Reduction of tube movement prevents nasal irritation▪ Centralizes tube in nostril with counterpulling action of over-wrap and underwrap tabs▪ Can accommodate one tube▪ More secure than adhesive tape
Non-Adherent Tips	<ul style="list-style-type: none">▪ “Starter” tab allows for easy application and removal from tube▪ No cutting required means no fear of accidentally cutting tube▪ Ease of use with gloved hands
Comfortable & Convenient	<ul style="list-style-type: none">▪ Allows nurse to independently change holder without waiting for additional staff
Not Made with Natural Rubber Latex	<ul style="list-style-type: none">▪ Eliminates patient and caregiver latex reactions

Learn More

