



CADD™ and Portex™ Portfolio with the NRFit® Connector

Product Catalogue

CADD and Portex Portfolio with the NRFit Connector

In 2014, The Joint Commission released Sentinel Event Alert 53: Managing Risk During Transition to ISO Tubing Connector Standards. The ISO 80369-6 neuraxial connector standard was published on March 2016. The goal is to reduce the incidence of medical device misconnections by creating unique non-Luer connectors for neuraxial procedures. The neuraxial non-Luer connectors will be branded by all major manufacturers as NRFit connectors.

ICU Medical would like to help you to switch to the NRFit neuraxial non-Luer connector. We'll help make the transition easier by providing support, guidance and product expertise.

A Familiar Way to Connect

The hubs and tapers of the NRFit connectors are smaller than those of standard Luer connectors, used commonly on many medical devices, but they operate much the same way. Like standard Luer connectors, NRFit connectors come together with a simple twist.




Security From Pump To Patient

The CADD and Portex Portfolio with the NRFit connector is an ISO 80369-6 compliant connection system intended to enhance patient safety through a reduction in the risk of misconnection, while minimising any change to clinical technique.



Portex Portfolio with the NRFIT Connector

Complete Solutions for Epidural Procedures



Epidural Needles					
Needle Gauge	Needle Type	Length	Connector	Product Code	Unit
● 16G	Tuohy	80 mm (3.1")	NRFit	24-1010-64	20/box
● 16G	Tuohy	110 mm (4.3")	NRFit	24-1011-64	10/box
○ 17G	Tuohy	90 mm (3.5")	NRFit	24-1014-64	10/box
● 18G	Tuohy	50 mm (2.0")	NRFit	24-1354-64	20/box
● 18G	Tuohy	80 mm (3.1")	NRFit	24-1012-64	20/box
● 18G	Tuohy	90 mm (3.5")	NRFit	24-1015-64	10/box
● 18G	Hustead	90 mm (3.5")	NRFit	24-1016-64	10/box
● 18G	Tuohy	110 mm (4.3")	NRFit	24-1013-64	10/box
● 19G	Tuohy	50 mm (2.0")	NRFit	24-1355-64	20/box
● 20G	Tuohy	90 mm (3.5")	NRFit	24-1017-64	10/box

Portex Epidural Catheter Kits with NRFit Connector

Clear Catheter Type	Length	Catheter ID Label	Catheter Guide	Catheter Connector	Outer Diameter	Product Code	Unit
19G Closed-end 3 Eyes	920 mm	•	•	NRFit	1.05 mm	24-1033-22	10/box
21G Closed-end 3 Eyes	920 mm	•	•	NRFit	0.85 mm	24-1034-22	10/box
19G Closed-end 3 Closer Eyes	920 mm	•	•	NRFit	1.05 mm	24-1035-22	10/box
21G Closed-end 3 Closer Eyes	920 mm	•	•	NRFit	0.85 mm	24-1036-22	10/box

Portex Epidural Single-Shot Packs with NRFit Connector

Tuohy Needle Gauge	Length	Connector	L.O.R. Syringe with NRFit Connector	Needle Filter with NRFit Connector	Product Code	Unit
● 16G	80 mm (3.1")	NRFit	Plastic, 10 mL, Slip	•	24-1206-22	10/box
● 18G	80 mm (3.1")	NRFit	Plastic, 10 mL, Slip	•	24-1207-22	10/box

Portex Epidural Minipack Systems 1 with NRFit Connector

Tuohy Needle Gauge	Length	With Attachable Wings	Connector	Clear Catheter with NRFit Connector & Guide	Plastic L.O.R. Slip with NRFit Connector	Catheter ID Label	With LockIt Plus™ Device	0.2 µm Epidural Filter with NRFit Connector	Product Code	Unit
● 16G	80 mm (3.1")	•	NRFit	Closed-End 3 Eyes	10 mL	•		•	24-1300-22	10/box
● 18G	80 mm (3.1")	•	NRFit	Closed-End 3 Eyes	10 mL	•		•	24-1301-22	10/box
● 16G	80 mm (3.1")	•	NRFit	Closed-End 3 Eyes	10 mL	•	•	•	24-1302-22	10/box
● 18G	80 mm (3.1")	•	NRFit	Closed-End 3 Eyes	10 mL	•	•	•	24-1303-22	10/box
● 16G	110 mm (4.3")	•	NRFit	Closed-End 3 Eyes	10 mL	•		•	24-1323-22	10/box
● 18G	110 mm (4.3")	•	NRFit	Closed-End 3 Eyes	10 mL	•		•	24-1324-22	10/box
● 16G	80 mm (3.1")	•	NRFit	Closed-End 3 Closer Eyes	10 mL	•		•	24-1328-22	10/box
● 18G	80 mm (3.1")	•	NRFit	Closed-End 3 Closer Eyes	10 mL	•		•	24-1325-22	10/box
● 18G	50 mm (2.0")	•	NRFit	Closed-End 3 Closer Eyes	10 mL	•		•	24-1335-22	10/box
● 19G	50 mm (2.0")	•	NRFit	Open-End	10 mL	•		•	24-1334-22	10/box

Portex Epidural Minipack Systems 2 with NRFit Connector

Tuohy Needle Gauge	Length	With Attachable Wings	Connector	Clear Catheter with NRFit Connector & Guide	Catheter ID Label	0.2 µm Epidural Filter with NRFit Connector	Product Code	Unit
● 16G	80 mm (3.1")	•	NRFit	Closed-End 3 Eyes	•	•	24-1326-22	10/box
● 18G	80 mm (3.1")	•	NRFit	Closed-End 3 Eyes	•	•	24-1327-22	10/box

Portex Epidural Maxipacks with NRFit Connector

	Tuohy Needle (80 mm/3") with NRFit Connector	With Attachable Wings	L.O.R. Slip with NRFit Connector	Clear Catheter with NRFit Connector	Catheter Identification Label	0.2 µm Epidural Filter with NRFit Connector	Product Code	Unit
●	16G	•	10 mL Plastic	Closed-end 3 Eyes	•	•	24-1304-22	10/box
●	18G	•	10 mL Plastic	Closed-end 3 Eyes	•	•	24-1305-22	10/box

Prep components: 2 hand towels, 5 gauze swabs, 2 applicator sponges, and solution prep wells.

Other procedural components: fenestrated drape, 20 mL lock syringe with NRFit connector, 19G x 38 mm (1.5") filter needle with NRFit connector, catheter fixation sponge, filter device (Luer), 3 mL plastic syringe, and hypodermic needles (22G x 40 mm – 1.5" and 25G x 16 mm – 5/8").



Portex Epidural Flat Filter with NRFit Connector

Filter Type	Connector	Membrane Pore Size	Priming Volume	Product Code	Unit
Flat	NRFit	0.2 µm hydrophilic 2-way	0.8 mL	24-1000-64	10/box



Portex Filter Needle with NRFit Connector

Type	Connector	Product Code	Unit
Filter Needle	NRFit	24-1356-64	10/box



Caps with NRFit Connector

Cap Type	Connector	Product Code	Unit
Yellow female	NRFit	24-1009-64	50/box
Yellow male non-vented	NRFit	24-1004-64	50/box



Portex Catheter Connector Key for use with NRFit Connectors

Description	Product Code	Unit
Key for catheter connector	24-1055-64	50/box



EpiFuse Catheter Connectors with NRFit Connector

Size	Connector	Product Code	Unit
19G-21G	NRFit	24-1005-64	25/box
24G	NRFit	24-1047-64	25/box



Loss of Resistance (L.O.R.) Syringes with NRFit Connector

Syringe Type	Connection	Connector	Capacity	Product Code	Unit
Plastic	Slip	NRFit	10 mL	24-1018-64	50/box
Glass	Slip	NRFit	10 mL	24-1019-64	10/box
Glass	Lock	NRFit	10 mL	24-1020-64	10/box
Plastic	Slip	NRFit	7 mL	24-1021-64	50/box
Plastic	Lock	NRFit	7 mL	24-1022-64	50/box
Glass	Slip	NRFit	5 mL	24-1023-64	10/box
Glass	Lock	NRFit	5 mL	24-1024-64	10/box



LockIt Plus Epidural Catheter Securement Device		
For catheter used with	Product Code	Unit
16G/17G needle	100/399/216	50/box
18G needle	100/399/218	50/box

NRFit Syringes		
Description	Product Code	Unit
NRFit Syringe, Lock 1 mL	15011010	100/box
NRFit Syringe, Lock 3 mL	15011030	100/box
NRFit Syringe, Lock 5 mL	15011050	100/box
NRFit Syringe, Lock 10 mL	15011100	100/box
NRFit Syringe, Lock 20 mL	15011200	50/box
NRFit Syringe, Lock 60 mL	15011600	25/box
NRFit Syringe, Slip 1 mL	15021010	100/box
NRFit Syringe, Slip 3 mL	15021030	100/box
NRFit Syringe, Slip 5 mL	15021050	100/box
NRFit Syringe, Slip 10 mL	15021100	100/box
NRFit Syringe, Slip 20 mL	15021200	50/box

CADD Portfolio with the NRFit Connector

Complete NRFit Solutions from the Pump to Patient

CADD-Solis™ Infusion Pump



CADD extension set with NRFit connector



CADD Administration Set with NRFit Connector



CADD Medication Cassette Reservoirs with NRFit Connector

Description	Length	Product Code	Unit
100 mL yellow CADD medication cassette reservoir with NRFit female connector, flow stop, clamp and non-vented stopper	20 cm (8 in)	21-7600-24	12/box
250 mL yellow CADD medication cassette reservoir with NRFit female connector, flow stop, clamp and non-vented stopper	20 cm (8 in)	21-7609-24	12/box



CADD Extension Sets with NRFit Connector

Description	Length	Priming Volume	Product Code	Unit
CADD yellow striped extension set with NRFit male/male connector, integral anti-siphon valve and clamp	229 cm (90 in)	2.3 mL	21-7605-24	12/box
CADD yellow-striped extension set with NRFit male/male connector, 0.2-micron air-eliminating filter and integral anti-siphon valve with clamp	229 cm (90 in)	2.4 mL	21-7608-24	50/box



CADD Administration Sets with NRFit Connector

Description	Length	Priming Volume	Product Code	Unit
CADD yellow striped administration set, bag spike, flow stop, clamp and one-way checkvalve with NRFit male connector	312 cm (123 in)	3.2 mL	21-7624-24	12/box
CADD yellow-striped administration set, bag spike, flow stop, 0.2-micron air-eliminating filter, clamp and one-way checkvalve with NRFit male connector	330 cm (130 in)	6.0 mL	21-7649-24	12/box
CADD yellow-striped high volume administration set, bag spike, flow stop, clamp and one-way checkvalve with NRFit male connector	312 cm (123 in)	16 mL	21-7684-24	15/box
CADD yellow-striped high volume administration set, bag spike, flow stop, 0.2-micron air-eliminating filter, clamp and one-way checkvalve with NRFit male connector	340 cm (134 in)	20 mL	21-7665-24	15/box

NRFit is a registered trademark of GEDSA

The product complies with current legislation and has the corresponding CE marking.
For additional information, warnings and / or safety precautions, refer to the manufacturer's user manual.

icumedical
human connections