

Sage Untreated Oral Swab

REF 5602UT-X



Packaging information

Contents:

- (1) Untreated Oral Swab

Indications of use: Aids in the removal of debris and oral secretions. May be used to apply oral care solutions.

Instructions: Moisten swab with water or solution. Clean mouth for 1 minute. Instruct to expectorate. Discard.

Caution:

Do not allow patient to bite down on the oral care tool. Use a bite block if patient has altered levels of consciousness or cannot comprehend commands. Use caution with children or unresponsive individuals. Failure to follow these safety precautions may damage the device and present a choking/aspiration hazard. May not function as intended/potential risk of cross-contamination if device is reused.

For single use.

Not made with natural rubber latex.

Oral Swab: MDD 93/42/EEC

Specifications

Product name	Sage Untreated Oral Swab
Product code	5602UT-X
Description	Applicator swab designed to aid in the removal of debris and oral secretions and apply oral care solutions
Manufacturer accreditation	ISO 13485:2016
Length of Swab	15.9 cm
Case quantity	800

Features and benefits

Swab

Swab	• Helps remove debris and oral secretions between brushings
Soft foam heads	• Gentle to fragile tissue

This document is intended solely for the use of healthcare professionals.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, Perox-A-Mint. All other trademarks are trademarks of their respective owners or holders.

The products depicted are CE marked in accordance with applicable EU Regulations and Directives.

This document is not intended for distribution outside EU and EFTA.



Sage Products LLC
3909 Three Oaks Road
Cary, IL 60013 USA
t: 0808 189 0456

sage-products.co.uk

25670-UK
Copyright © 2018 Stryker