



MiniACE® Button Cecostomy + Appendicostomy

MiniACE® Button

Low Profile Antegrade Continence Enema Button

MiniACE® Button Sizes

New! 10F Stoma Size and Longer-Length 12F Devices!

AMT	10F	12F	14F
LENGTH	ORDER NUMBERS		
0.8 cm	ACE-1008	ACE-1208	ACE-1408
1.0 cm	ACE-1010	ACE-1210	ACE-1410
1.2 cm	ACE-1012	ACE-1212	ACE-1412
1.5 cm	ACE-1015	ACE-1215	ACE-1415
1.7 cm	ACE-1017	ACE-1217	ACE-1417
2.0 cm	ACE-1020	ACE-1220	ACE-1420
2.3 cm	ACE-1023	ACE-1223	ACE-1423
2.5 cm	ACE-1025	ACE-1225	ACE-1425
2.7 cm	ACE-1027	ACE-1227	ACE-1427
3.0 cm	ACE-1030	ACE-1230	ACE-1430
3.5 cm	ACE-1035	ACE-1235	ACE-1435
4.0 cm	ACE-1040	ACE-1240	ACE-1440
4.4 cm	ACE-1044	ACE-1244	ACE-1444
5.0 cm	ACE-1050	ACE-1250	ACE-1450
5.5 cm	ACE-1055	ACE-1255	ACE-1455
6.0 cm	ACE-1060	ACE-1260	ACE-1460
6.5 cm	ACE-1065	ACE-1265	ACE-1465
7.0 cm		ACE-1270	ACE-1470
8.0 cm		ACE-1280	ACE-1480
9.0 cm		ACE-1290	ACE-1490
10.0 cm		ACE-12100	ACE-14100

MiniACE® Kits are packaged with a 5 cc syringe, an introducer, gauze, and a lube packet. A sterile pouched 2-1221 bolus irrigation set is also included in the kit.

Irrigation Set Included in Kit



Complete Portfolio – Irrigation Sets

Adapter Options	Right Angle Connector Options	
	12"	24"
Bolus	2-1221	2-2421

Note: Irrigation sets are packaged as 10/box (non-sterile)

For a complete list of products visit us at www.AppliedMedical.net
800 869 7382 | CS@AppliedMedical.net

LP Initial Placement Cleared

MR Conditional




Brobekkveien 107 * 0582 Oslo
Tlf: 22 72 03 00 Fax: 22 72 03 01
post@kvinto.no www.kvinto.no

© 2019 Applied Medical Technology, Inc. C-4642-B
Patents: <http://www.AppliedMedical.net/Patents/> | Trademarks: <http://www.AppliedMedical.net/Trademarks/>
Chait Trapdoor™ is a trademark of Cook® Medical | 1.Claim based on internal testing of 10F-14F devices.

Balloon Fill Volumes

French Size	Min	Rec.	Max
10F	1 ml	1.5 ml	2 ml
12F	2 ml	2.5 ml	3 ml
14F	3 ml	4 ml	5 ml

Recommended water fill volume is clearly printed on the balloon fill valve of each MiniACE® Button. Consult your healthcare provider before adjusting fill volumes.



River-Kai proudly rocks two AMT devices, the MiniONE® and the MiniACE®!!

[Our surgeons have commented] that they've never seen a stoma look so good. I owe that to AMT for such a great product! I recommend AMT products to any parent who complains about granulation. – River-Kai's Mom

River-Kai's MiniACE®

AppliedMedical.net/River-Kais-MiniONE-and-MiniACE



NEW! 10F



The Only FDA Cleared
Balloon Button Device for
Antegrade Enemas

Cleared for Initial Placement

Superior Leakage Control

Improved Stoma Retention




Innovating. Educating. Changing Lives.™




MiniACE® Button

- The only FDA cleared balloon button device for Antegrade Continence Enema (ACE) procedures
- Cleared for Initial Placement
- **NOW IN 10F!** Minimizes stoma dilation associated with conversion from Chait to MiniACE®
- Cecostomy + Appendicostomy



Cecostomy




Appendicostomy

Patient Comfort First

- Medical grade silicone, gentle to the skin
- Flexible + Soft, easy to clean around the button
- Exclusive AMT balloon, designed to reduce pullouts and leakage
- Balloon bolster shields the internal tip of the MiniACE®, protecting the small bowel


Reduce Leaks, Superior Retention

- AMT's MiniACE® balloon button is 450% more secure within the stoma when compared to the Chait Trapdoor™¹
- Tapered tubing at the external bolster, along with the exclusive AMT Balloon, reduce the risk of leaks.
- Anti-leakage valve helps prevent backflow.



Safe and Secure Irrigation Set Attachment


- The connection between the MiniACE® Button and its irrigation set is 10x stronger than the competition.¹




STABILIZE



ALIGN & INSERT



ROTATE



LOCK



IRRIGATE