


229 Under Buttocks Drape

Product details

Size :	85cm x 120cm
Descriptive feature :	Cuffs, Patched, Pouch
Sterile :	Sterile
Brand :	BARRIER
Equivalences :	9229

Images



Delivered items

229-21

Sales released in: Austria, Bahrain, Belgium, Chile, Colombia, Cyprus, Czechia, Denmark, Estonia, Finland, France, Germany, Greece, Hong Kong, Ireland, Italy, Japan, Kuwait, Latvia, Lithuania, Luxembourg, Malta, Martinique, Mexico, Netherlands, New Zealand, Norway, Panama, Poland, Portugal, Qatar, Saudi Arabia, Singapore, Spain, Sweden, Switzerland, Thailand, United Arab Emirates, United Kingdom of Great Britain and Northern Ireland

Country of origin: Thailand

Shelf life: 5 years

Sterilization method: Irradiation

Production Responsibility: Mölnlycke Health Care (Thailand) Ltd, 160 Moo 17 Bangplee Industrial Estate, Bangna-Trad Road, Km 23 Tambol Bangsaothong, Amphur Bangsaothong, Samutprakarn 10570, Thailand

Packing information: First packaging layer is a peel-open sterile barrier, plastic/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a plastic liner/bag. Third layer is a corrugated board dispenser box. Fourth layer is a corrugated board transport box.

Is suitable for Tray: No

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2022 Mölnlycke Health Care AB. All rights reserved.



Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Piece	1	7332430080851			
Consumer pack	1	7323190294510			
Dispenser box	20	7323190293278			
Transport box	40	7323190293261	397x580x329 mm	75.8 dm3	9.1 / 7.4 kg
Pallet	800	7323190293254	800x1200x1789 mm	1717.4 dm3	207.8 / 148.0 kg

Technical

Dimension

Dimension text	Dimension value
Outer dimension	85 cm x 120 cm
Patch	50 cm x 35 cm
Pouch	82 cm x 79 cm

Disposal instructions

Used Mölnlycke products and sterile barriers should, in the majority of cases, be classified as non-hazardous waste. They contain high amounts of energy and are well suited for incineration. Mölnlycke products do not contain any hazardous substances that can leach out if the products are land filled. Transport- and dispenser boxes are designed to fit existing recovery systems. The Mölnlycke packaging system complies with the EU Directive on packaging and packaging waste.

Storage instructions

Mölnlycke Health Care recommends that Mölnlycke products are stored under normal storage conditions. All layers of packaging should be kept intact until access to the underlying layers is needed. Storage facilities for products only protected by the sterile barrier should be kept under conditions where low level of particulate air contamination prevail, so that it would not constitute a risk to the patient when the sterile barrier is opened and the product is used.

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2022 Mölnlycke Health Care AB. All rights reserved.



Classifications

Regulation type(s)	MDR Class IS	MDR Class IS	Locally Regulated	Locally Regulated
CE Certificate Number :	MDR 722028			
Intended Purpose :	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.
MDR Classification Rule :	1			
Conformity Annexes :	IX			
Measuring Function :	No			
Notified body medical devices/PPE :	BSi 2797	BSi 2797	Local Authority Americas	Local Authority Thailand
Regulatory Released :	Austria, Belgium, Czechia, Denmark, Estonia, Finland, France, Germany, Greece, Ireland, Italy, Latvia, Lithuania, Luxembourg, Netherlands, Norway, Poland, Portugal, Spain, Sweden, Switzerland, United Kingdom of Great Britain and Northern Ireland	Austria, Belgium, Czechia, Denmark, Estonia, Finland, France, Germany, Greece, Ireland, Italy, Latvia, Lithuania, Luxembourg, Netherlands, Norway, Poland, Portugal, Spain, Sweden, Switzerland, United Kingdom of Great Britain and Northern Ireland	Chile, Colombia, Panama	Thailand

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
 The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
 ©2022 Mölnlycke Health Care AB. All rights reserved.



Classifications

Regulation type(s)	Locally Regulated	Locally Regulated	Locally Regulated	Class I
CE Certificate Number :				
Intended Purpose :	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.
MDR Classification Rule :				
Conformity Annexes :				
Measuring Function :				
Notified body medical devices/PPE :	Local Authority Hong Kong	Local Authority Mexico	Local Authority APAC	PMDA (Japan)
Regulatory Released :	Hong Kong	Mexico	New Zealand, Singapore	Japan

Applied standards : The standards presented below is a selection of the most essential standards that are adhered to.

EN 1041, EN 556-1, EN 13795, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 11607-1, EN ISO 11607-2, EN ISO 15223-1, ISO 15223-2, ISO 14001

Removable label

Yes

GMDN Code (Global Medical Device Nomenclature)

47783 Patient surgical drape, single-use

UNSPSC

42132103 Patient barrier drapes

Commodity Code

3005905000 Articles of Wadding, Gauze, Dressings, drapes singlepacked and Sets - textile

CE Responsibility

Molnlycke Health Care AB, Gamlestadsavagen 3c, Box 13080, SE-402 52 Goteborg Sweden

Basic UDIDI

73324300000000177KG

Find out more at www.molnlycke.com

