


84541230 Dislocation Bag

Product details

Size : 120cm x 75cm
Descriptive feature : Adhesive, Wrapped
Sterile : Sterile
Brand : BARRIER
Equivalences : 922709

Images



Delivered items

84541230-21

Sales released in: Algeria, Belgium, Czechia, Finland, France, Germany, Ireland, Lithuania, Luxembourg, Montenegro, Netherlands, Norway, Poland, Qatar, Serbia, Singapore, South Africa, Spain, Sweden, Switzerland, Thailand, United Arab Emirates, United Kingdom of Great Britain and Northern Ireland

Shelf life: 5 years

Sterilization method: Irradiation

Production Responsibility: Mölnlycke Health Care (Thailand) Ltd, 160 Moo 17 Bangplee Industrial Estate, Bangna-Trad Road, Km 23 Tambol Bangsaothong, Amphur Bangsaothong, Samutprakarn 10570, Thailand

Packing information: First packaging layer is a peel-open sterile barrier, plastic/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a corrugated board dispenser box. Third layer is a corrugated board transport box.

Is suitable for Tray: No

Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Piece	1	7332430118714			

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2022 Mölnlycke Health Care AB. All rights reserved.



Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Consumer pack	1	7323190292417			
Dispenser box	16	7323190289530			
Transport box	32	7323190289523	332x428x288 mm	40.9 dm3	5.2 / 4.3 kg
Pallet	1152	7323190289516	800x1200x1872 mm	1797.1 dm3	212.2 / 153.7 kg

84541230-02

Sales released in: Algeria, Belgium, Czechia, Denmark, Finland, France, Germany, Ireland, Lithuania, Luxembourg, Montenegro, Netherlands, Norway, Poland, Qatar, Serbia, Singapore, South Africa, Spain, Sweden, Switzerland, Thailand, United Arab Emirates, United Kingdom of Great Britain and Northern Ireland

Shelf life: 5 years

Sterilization method: Irradiation

Production Responsibility: Mölnlycke Health Care (Thailand) Ltd, 160 Moo 17 Bangplee Industrial Estate, Bangna-Trad Road, Km 23 Tambol Bangsaothong, Amphur Bangsaothong, Samutprakarn 10570, Thailand

Packing information: First packaging layer is a peel-open sterile barrier, plastic/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a corrugated board dispenser box. Third layer is a corrugated board transport box.

Is suitable for Tray: No

Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Piece	1	7332430118714			
Consumer pack	1	7323190039906			
Dispenser box	14	7323190039142			
Transport box	84	7323190039135	366x753x295 mm	81.3 dm3	13.2 / 11.2 kg
Pallet	1512	7323190039128	800x1200x1914 mm		

Technical

Dimension

Dimension text	Dimension value
Outer dimension	120 cm x 75 cm

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2022 Mölnlycke Health Care AB. All rights reserved.



Disposal instructions

Used Mölnlycke products and sterile barriers should, in the majority of cases, be classified as non-hazardous waste. They contain high amounts of energy and are well suited for incineration. Mölnlycke products do not contain any hazardous substances that can leach out if the products are land filled. Transport- and dispenser boxes are designed to fit existing recovery systems. The Mölnlycke packaging system complies with the EU Directive on packaging and packaging waste.

Storage instructions

Mölnlycke Health Care recommends that Mölnlycke products are stored under normal storage conditions. All layers of packaging should be kept intact until access to the underlying layers is needed. Storage facilities for products only protected by the sterile barrier should be kept under conditions where low level of particulate air contamination prevail, so that it would not constitute a risk to the patient when the sterile barrier is opened and the product is used.

Classifications

Regulation type(s)	MDR Class IS	Locally Regulated	Locally Regulated	Locally Regulated
CE Certificate Number :	MDR 722078;MDR 722028;MDR 722028;MDR 722028;MDR 722028;MDR 722028;MDR 722028;MDR 722028;MDR 722028;MDR 722028;MDR 722028;MDR 722028			
Intended Purpose :	The Dislocation Bag is attached to other draped areas and used to maintain a dislocated leg within the sterile field during a surgical procedure.	The Dislocation Bag is attached to other draped areas and used to maintain a dislocated leg within the sterile field during a surgical procedure.	The Dislocation Bag is attached to other draped areas and used to maintain a dislocated leg within the sterile field during a surgical procedure.	The Dislocation Bag is attached to other draped areas and used to maintain a dislocated leg within the sterile field during a surgical procedure.
MDR Classification Rule :	1;1;1;1;1;1;1;1;1;1;1;1;1			
Conformity Annexes :	IX;IX;IX;IX;IX;IX;IX;IX;IX;IX;IX;IX;IX			
Measuring Function :	No			
Notified body medical devices/PPE :	BSi 2797	Local Authority Thailand	Local Authority MEA	Local Authority United Arab Emirates

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2022 Mölnlycke Health Care AB. All rights reserved.



Classifications

Regulation type(s)	Locally Regulated	Locally Regulated
CE Certificate Number :		
Intended Purpose :	The Dislocation Bag is attached to other draped areas and used to maintain a dislocated leg within the sterile field during a surgical procedure.	The Dislocation Bag is attached to other draped areas and used to maintain a dislocated leg within the sterile field during a surgical procedure.
MDR Classification Rule :		
Conformity Annexes :		
Measuring Function :		
Notified body medical devices/PPE : Local Authority APAC		Local Authority EU Non EEC

Applied standards : The standards presented below is a selection of the most essential standards that are adhered to.

EN 1041, EN 556-1, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 11607-1, EN ISO 11607-2, EN ISO 15223-1, ISO 15223-2, ISO 14001

Removable label

Yes

GMDN Code (Global Medical Device Nomenclature)

47783 Patient surgical drape, single-use

UNSPSC

42131701 Surgical drapes

Commodity Code

3926909790 Plastics and articles thereof

CE Responsibility

Mölnlycke Health Care AB, Gamlestadsvägen 3c, Box 13080, SE-402 52 Göteborg Sweden

Basic UDIDI

733243000000000179KL

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2022 Mölnlycke Health Care AB. All rights reserved.

