

# Blood collection devices and accessories

## Safety blood collection sets

### BD Vacutainer® Push Button Blood Collection Set

The BD Vacutainer® Push Button Blood Collection Set with in-vein activation offers split-second protection against needlestick injuries.

- Indication of successful venepuncture:  
When the vein has been successfully penetrated, blood flows immediately into the inspection chamber.<sup>70</sup>
- Versatile:  
For taking blood samples and for short-term infusions of up to 2 hours.
- Single-handed activation possible:  
Activating the safety mechanism with one hand allows you to focus more attention on the patient and the venepuncture site.
- Protection against needlestick injuries:  
When pressing the button, the needle is withdrawn straight from the vein and disappears permanently inside the housing of the blood collection set. This provides an extremely high level of protection<sup>70</sup>.












### Further information

BD recommends for all needles without a pre-attached holder that a BD Vacutainer® One-Use Holder is used.

Clinical and technical information is available on request.

### BD Vacutainer® Push Button Blood Collection Sets

Cat. no.	Size	Needle length	Length of tubing	With Luer adapter	Colour code	Box/Case
367338	21 G (0.8 mm)	19 mm	178 mm	Yes		50/200
367344	21 G (0.8 mm)	19 mm	305 mm	Yes		50/200
367326	21 G (0.8 mm)	19 mm	305 mm	No		50/200
367336	23 G (0.6 mm)	19 mm	178 mm	Yes		50/200
367342	23 G (0.6 mm)	19 mm	305 mm	Yes		50/200
367324	23 G (0.6 mm)	19 mm	305 mm	No		50/200
367335	25 G (0.5 mm)	19 mm	178 mm	Yes		50/200
367341	25 G (0.5 mm)	19 mm	305 mm	Yes		50/200
367323	25 G (0.5 mm)	19 mm	305 mm	No		50/200

# Blood collection devices and accessories

## Safety blood collection sets

### BD Vacutainer® Push Button Blood Collection Set with Pre-Attached Holder

With this Safety Blood Collection Set, the holder is already pre-attached, so it is not necessary to manually assemble the needle and holder. This ready-for-use blood collection set is individually packaged in a sterile blister.





This sterile closed system is ideally suited for taking of samples using the BD BACTEC™ Blood Culture Bottles.

#### Further Information

Technical information is available on request.



### BD Vacutainer® Push Button Blood Collection Sets with Pre-Attached Holder

Cat. no.	Size	Needle length	Length of tubing	Colour code	Box/Case
367355	21 G (0.8 mm)	19 mm	178 mm		20/100
368657	21 G (0.8 mm)	19 mm	305 mm		20/100
367354	23 G (0.6 mm)	19 mm	178 mm		20/100
368658	23 G (0.6 mm)	19 mm	305 mm		20/100

**Ready-  
for-use**

blood  
collection set