

Technical Data Sheet

10/05/2019 - rev. 03

Item	ClinMat 200		
	Three-layered super absorbent mat for operating theaters - antislip on floor		
	Medical Device Class I - Not sterile - according with Directive n. 93/42/CEE CND A0699		
Composition	Top layer	: perforated PE film - light blue	
	Absorbent nonwoven	: 60% Viscose fibers - 40% SAF fibers (200 gr/m2 ± 10%)	
	Bottom layer	: antislip (adhesive) polyolefin film	
Mat weight	300 gr/m2		± 10%
Thickness	2,9 mm		± 10%
Absorbency	5 lt/m2		max.

Available in:

ROLLS

to be cut according to needs

Roll width	: 90 cm	± 5%
Roll length	: 50 or 100 lin.mt	aprox.

other sizes on request

Splices per roll	: max. 4	
Ext. diameter of roll	: max. 48 cm	± 10%

Packaging : PE bag /single roll

Item code	: rotoli cm 90x50 mt	CM200H90L50A
	: rotoli cm 90x100 mt	CM200H90L100A

Instruction for use

At the time of use, the mat has to be unrolled according to the requested size and cut by common scissors. The remaining material – in case of partial use – must be immediately put in its bag and accurately closed. The remaining material – in case of partial use – must be immediately put in its bag and accurately closed.

The mat has to be placed on the floor. The adhesive white side under in contact with the floor; the light blue side, where the holes of the protective film are visible, on top for facilitating the absorption.

CUT IN SIZE

Standard sizes	: 90x45 - 90x76 - 90x78 - 90x116 - 90x150 90x200 cm
----------------	--

other sizes on request

Packaging	primary	: PE bag	/single piece
	secondary	: cardboard box	/quantity
		90x45 cm	55 pcs
		90x76 cm	35 pcs
		90x78 cm	35 pcs

Technical Data Sheet

10/05/2019 - rev. 03

Item

ClinMat 200

Three-layered super absorbent mat for operating theaters - antislip on floor

90x116 cm 23 pcs

90x150 cm 18 pcs

90x200 cm 16 pcs

Item code	: size 90x45 cm	CM2009045A
	: size 90x76 cm	CM2009076A
	: size 90x78 cm	CM2009078A
	: size 90x116 cm	CM20090116A
	: size 90x150 cm	CM20090150A
	: size 90x200 cm	CM20090200A

Instruction for use

The mat has to be placed on the floor. The adhesive white side under in contact with the floor; the light blue side, where the holes of the protective film are visible, on top for facilitating the absorption.

Warnings

This material is extremely sensitive to the humidity. It is recommended to store it well sealed and protected by the contact with air and sun light and by external contaminations until the use.

Pay attention to alterations (reducing of absorbency degree) that can occur in case it is left open in contact with humidity.

The product and its components comply with current rules about air emissions, water protection and content of harmful toxic substances. Dispose of in accordance with local requirements.



This item is produced in Italy