

The S-Cath™ System

- Seldinger based suprapubic catheter insertion system
- 3 stage guidewire ensuring accurate placement along one anaesthetised tract. ^{1, 3, 4, 5, 6}
- Decrease the risk of posterior bladder wall perforations and bowel injuries. ^{1, 3, 4, 5, 6, 9}
- Facilitate use in an outpatient or sub acute setting. ^{1, 3, 4, 5, 7}
- Greater confidence for the user. ^{8, 10}
- 100% silicone Foley catheter with integrated balloon



Data

- **1 - Department of Health Briefing pack**
S-Cath™ System Suprapubic Foley Catheter Introduction Set, August 2012
- **2 - British Association of Urological Surgeons' Suprapubic Catheter Practice Guidelines**
Simon C.W. Harrison*, William T. Lawrence†, Roland Morley‡, Ian Pearce and Joby Taylor§
*Pinderfields General Hospital, Wakefield, †Eastbourne District General Hospital, Eastbourne, ‡Kingston Hospital, Kingston upon Thames, and §Manchester Royal Infirmary, Manchester, UK.
BJU International, Volume 107, Issue 1, pages 77–85, January 2011
- **3 - A dedicated Specialist Nurse led Suprapubic Catheterisation Clinic: The Taunton Experience**
J.Jelski, Xi Cheng, Faith McMeekin, A. MacCormick, R. Macdonagh, Musgrove Park Hospital, Taunton, UK
- **4 - Introducing a New Technique for Suprapubic Catheterisation**
McMeekin F, MacCormick A, Macdonagh R, Musgrove Park Hospital, Taunton, UK
Continence UK, 2010, Vol 4, No 1
- **5 - Optimising Safety of Catheterisation**
The Clinical Services Journal, October 2009
- **6 - How to Why to Guide - The Seldinger Suprapubic Catheter Kit**
NHS National Technology Adoption Centre (NTAC) April 2009
- **7 - Suprapubic catheter insertion is an outpatient procedure: cost savings resultant on closing an audit loop**
Khan A, Abrams P, Bristol Urological Institute, Southmead Hospital, Bristol, UK
BJU Int. 2009 Mar;103(5):640-4. Epub 2008 Oct 24.
- **8 - Suprapubic bladder catheterisation using the Seldinger technique**
Vasdev V, Kachroo N, Mathur S, Pickard R., Department of Urology, Freeman Hospital, Newcastle upon Tyne, UK; BioMed Centre, Bristol Urological Institute, Bristol, UK
The Internet Journal of Urology™ ISSN: 1528-8390, Feb 09
- **9 - A new model for suprapubic catheterization: the MediPlus Seldinger suprapubic catheter**
Mohammed A, Khan A, Shergill IS, Gujral SS. ST1 in Surgery/Urology, King George Hospital, Ilford, London, IG3 8QY, UK
Expert Rev Med Devices. 2008 Nov;5(6):705-7.
- **10 - Suprapubic catheter insertion using the Seldinger technique, with the Mediplus SPC kit**
D.M. Gulur, F. Housami, M.J. Drake, Bristol Urological Institute, Bristol, UK
BJUI, www.bjui.org, Atlas of Surgery, Mar 2008

Product Code	Description	Size	Length	Box Quantity
5758	S-Cath™ System Silicone Suprapubic Foley Catheter Set with Introducer	8Fr.	33cm	5
5760	S-Cath™ System Silicone Suprapubic Foley Catheter Set with Introducer	12Fr.	42cm	5
5754	S-Cath™ System Silicone Suprapubic Foley Catheter Set with Introducer	14Fr.	42cm	5
5757	S-Cath™ System Silicone Suprapubic Foley Catheter (Open Tip) Set with Introducer	14Fr.	42cm	5
5756	S-Cath™ System Silicone Suprapubic Foley Catheter Set with Introducer	16Fr.	42cm	5
5759	S-Cath™ System Silicone Suprapubic Foley Catheter (Open Tip) Set with Introducer	16Fr.	42cm	5

Apr-15/1

mediplus

*Working together
to make the difference*

Mediplus Ltd

Unit 7 The Gateway Centre, Coronation Road, Cressex Business Park, High Wycombe, Bucks, HP12 3SU UK

Tel: +44 (0)1494 551200 | Fax: +44 (0)1494 536333

Email: help@mediplus.co.uk | Web: www.mediplus.co.uk



www.mediplus.co.uk