



Electrosurgical Accessories

EDGE™ COATED ELECTRODES

Standard EDGE™ Coated Electrodes



E1450X

Total length: 6.35 cm (2.5 in.) • Active length: 2.8 cm (1.1 in.)



E1450G

Total length: 7.62 cm (3 in.) • Active length: 2.8 cm (1.1 in.)



E1475X

For use with AccuVac™ smoke evacuation attachment
Total length: 6.99 cm (2.75 in.) • Active length: 2.8 cm (1.1 in.)



E1452

Total length: 7.2 cm (2.84 in.) • Active length: 2.8 cm (1.1 in.)



E1450-4

Total length: 10.16 cm (4 in.) • Active length: 2.54 cm (1 in.)



E1450-6

Total length: 16.51 cm (6.5 in.) • Active length: 2.54 cm (1 in.)



E1452-6

Total length: 16.51 cm (6.5 in.) • Active length: 2.54 cm (1 in.)

Product shown actual size

COVIDIEN, COVIDIEN with logo, Covidien logo and ™ marked brands are internationally registered trademarks of Covidien AG. © 2010 Covidien.

contact www.covidien.com

Insulated EDGE™ Coated Electrodes



E1455

Total Length: 6.99 cm (2.75 in.) • Active Length : 5.1 mm (0.2 in.)



with Safety Sleeve™

E1455B

Total Length: 6.99 cm (2.75 in.) • Active Length : 5.1 mm (0.2 in.)



E1455-4

Total Length: 10.16 cm (4 in.) • Active Length : 5.1 mm (0.2 in.)



with Safety Sleeve™

E1455B-4

Total Length: 10.16 cm (4 in.) • Active Length : 5.1 mm (0.2 in.)



E1455-6

Total Length: 16.51 cm (6.5 in.) • Active Length : 5.1 mm (0.2 in.)



E1465

Total Length: 7.21 cm (2.84 in.) • Active Length : 5.1 mm (0.2 in.)



with Safety Sleeve™

E1465B

Total Length: 7.21 cm (2.84 in.) • Active Length : 5.1 mm (0.2 in.)



E1465-4

Total Length: 10.16 cm (4 in.) • Active Length : 5.1 mm (0.2 in.)



with Safety Sleeve™

E1465B-4

Total Length: 10.16 cm (4 in.) • Active Length : 5.1 mm (0.2 in.)



E1465-6

Total Length: 16.51 cm (6.5 in.) • Active Length : 5.1 mm (0.2 in.)

STAINLESS STEEL ELECTRODES



E1550

Total length: 5.33 cm (2.1 in.) • Active length: 1 cm (0.4 in.) • Diameter: 5 mm (3/16 in.)



E1551G

Total length: 6.2 cm (2.44 in.) • Active length: 2.8 cm (1.1 in.)



E1551X

Total length: 6.2 cm (2.44 in.) • Active length: 2.54 cm (1 in.)



E1551-6

Total length: 16.51 cm (6.5 in.) • Active length: 2.54 cm (1 in.)



E1552

Total length: 7.2 cm (2.84 in.) • Active length: 2.8 cm (1.1 in.)



E1552-6

Total length: 16.51 cm (6.5 in.) • Active length: 2.54 cm (1 in.)



E1502

Straight Electrode Extension - Reusable

Total length: 13 cm (5.1 in.) • Active length: 10.2 cm (4 in.)



E1504

Straight Electrode Extension - Reusable

Total length: 34.3 cm (13.5 in.) • Active length: 31.5 cm (12.4 in.)

Product shown actual size

COVIDIEN, COVIDIEN with logo, Covidien logo and TM marked brands are internationally registered trademarks of Covidien AG. © 2010 Covidien.

contact www.covidien.com


Monopolar, Footswitching Cords



E0510
**Footswitching Monopolar
Laparoscopic Cord**
Requires E0502 Series or E0017 active adapter
Accepts 4 mm male single prong plug
Length: 3 m (10 ft)



E2999
**Footswitching
Laparoscopic Cord - Reusable**
Accepts 4 mm male single prong plug
Length: 3 m (10 ft)



360187
Footswitching Cord
*For use with VALLEYLAB generators
& with Monopolar Forceps 713XXX.*
Length: 3.5 m

Bipolar, Footswitching Cords



E0509
**Footswitching Bipolar
Forceps Cord**
Length: 3.6 m (12 ft)



E0512
**Footswitching Bipolar
Forceps Cord**
For use with Valleylab generators only
Length: 3.6 m (12 ft)



E0018
Handswitching Forceps Cord - Reusable
For use with Valleylab generators only
Length: 3.6 m (12 ft)

Product shown actual size

COVIDIEN, COVIDIEN with logo, Covidien logo and TM marked brands are internationally registered trademarks of Covidien AG. © 2010 Covidien.

contact www.covidien.com

Bipolar, Reusable Cord

for use with E700XXX bipolar Forceps.



Actual size ?

E0022W

Handswitching Forceps Cord

For use with WOLF generators only

Length: 4.5 m (12 ft)



Actual size ?

E0020V

Footswitching Forceps Cord

For use with VALLEYLAB generators only.

Length: 4.5 m (12 ft)



Actual size ?

E0021S

Handswitching Forceps Cord

For use with ERBE/STORZ generators only

Length: 4.5 m (12 ft)

COAGULATION FORCEPS

Bipolar, Footswitching Coagulation Forceps

Use with E0509 or E0512 single use cords. Use with DQ12401 or E0019 reusable cords.



E4051-CT

Scoville-Greenwood Bayonet Forceps

Total length: 19.7 cm (7.75 in.) • Tip: 1.5 mm



E4052-CT

Cushing Forceps

Total length: 15.2 cm (6 in.) • Tip: 0.7 mm



E4053-CT

Cushing Forceps

Total length: 17.8 cm (7 in.) • Tip: 1.5 mm



E4054-CT

Semkin Forceps with Stops

Total length: 19.1 cm (7.5 in.) • Tip: 2.0 mm

Product shown actual size

COVIDIEN, COVIDIEN with logo, Covidien logo and TM marked brands are internationally registered trademarks of Covidien AG. © 2010 Covidien.

contact www.covidien.com

Bipolar, Footswitching Coagulation Forceps

Use with E0509 or E0512 single use cords. Use with DQ12401 or E0019 reusable cords.



E4057-CT

Gerald Bayonet Forceps

Total length: 19.1 cm (7.5 in.) • Tip: 0.7 mm



E4058-CT

Cushing Bayonet Forceps

Total length: 19.1 cm (7.5 in.) • Tip: 0.7 mm



E4059-CT

Hardy Bayonet Forceps with Stops

Total length: 20.9 cm (8.25 in.) • Tip: 0.5 mm

Bipolar, Footswitching Coagulation Forceps

Use with E0509 or E0512 single use cords. Use with DQI2401 or E0019 reusable cords.



E4060-CT

Semkin Forceps with Stops

Total length: 14 cm (5.5 in.) • Tip: 0.5 mm



E4073-CT

**Titanium Bayonet
Forceps Straights**

Total length: 22.2 cm (8.75 in.) • Tip: 1.0 mm



E4055-CT

Jewlers Forceps, Straight Tip

Total length: 10.2 cm (4 in.) • Tip: 0.4 mm



E4062-CT

Curved Iris Forceps

Total length: 8.9 cm (3.5 in.) • Tip: 0.5 mm

Product shown actual size

COVIDIEN, COVIDIEN with logo, Covidien logo and TM marked brands are internationally registered trademarks of Covidien AG. © 2010 Covidien.

contact www.covidien.com

Bipolar, micro Footswitching Coagulation Forceps

Use with E360XXX reusable cord.



E700150

Total length: 16.5 cm (6.5 in.)
straight, pointed



E700151

Total length: 16.5 cm (6.5 in.)
Tip: 1 mm, straight



E700155

Total length: 16.5 cm (6.5 in.)
angled, pointed, delicate tines



E700157

Total length: 16.5 cm (6.5 in.)
Tip: 1 mm, angled



E700158

Total length: 16.5 cm (6.5 in.)
Tip: 0.7 mm, short angled

Bipolar, micro Footswitching Coagulation Forceps

Use with E360XXX reusable cord.



E700167

Total length: 16.5 cm (6.5 in.)
Tip: 1 mm, bayonet



E700170

Total length: 19 cm (7.5 in.)
Tip: 1 mm, pointed



E700171

Total length: 19 cm (7.5 in.)
Tip: 1 mm, straight

Product shown actual size

COVIDIEN, COVIDIEN with logo, Covidien logo and TM marked brands are internationally registered trademarks of Covidien AG. © 2010 Covidien.

contact www.covidien.com

Bipolar, micro Footswitching Coagulation Forceps

Use with E360XXX reusable cord.



E700175

Total length: 19 cm (7.5 in.)
Tip: 1 mm, angled



E700176

Total length: 19 cm (7.5 in.)
Tip: 2 mm, angled



E700181

Total length: 19 cm (7.5 in.)
Tip: 1 mm, bayonet

Bipolar, micro Footswitching Coagulation Forceps

Use with E360XXX reusable cord.



E700184

Total length: 19 cm (7.5 in.)
Tip: 1 mm, bayonet, angled downwards



E700185

Total length: 19 cm (7.5 in.)
Tip: 1 mm, bayonet, angled upwards



E700191

Total length: 23 cm (9.1 in.)
Tip: 0.7 mm, bayonet



E700192

Total length: 23 cm (9.1 in.)
Tip: 1.2 mm, bayonet



E700195

Total length: 23 cm (9.1 in.)
Tip: 1.2 mm, bayonet, angled upwards

Bipolar, micro Footswitching Coagulation Forceps

Use with E360XXX reusable cord.



E700204

Total length: 14.5 cm (5.75 in.)

Tip: 1 mm



E700206

Total length: 16.5 cm (6.5 in.)

Tip: 1 mm



E700209

Total length: 19 cm (7.5 in.)

Tip: 1 mm



E700212

Total length: 22 cm (8.75 in.)

Tip: 1 mm, angled



E700214

Total length: 14.5 cm (5.75 in.)

Tip: 2 mm

Bipolar, micro Footswitching Coagulation Forceps

Use with E360XXX reusable cord.



E700216

Total length: 16.5 cm (6.5 in.)

Tip: 2 mm



E700219

Total length: 19 cm (7.5 in.)

Tip: 2 mm



E700222

Total length: 22 cm (8.75 in.)

Tip: 2 mm



E700240

Total length: 10.5 cm (4.25 in.)

straight, pointed



E700241

Total length: 10.5 cm (4.25 in.)

straight, blunt

Product shown actual size

COVIDIEN, COVIDIEN with logo, Covidien logo and TM marked brands are internationally registered trademarks of Covidien AG. © 2010 Covidien.

contact www.covidien.com

Bipolar, micro Footswitching Coagulation Forceps

Use with E360XXX reusable cord.



E700244

Total length: 10.5 cm (4.25 in.)
short angled, pointed



E700245

Total length: 10.5 cm (4.25 in.)
long angled, pointed



E700246

Total length: 10.5 cm (4.25 in.)
long angled, blunt

Bipolar, micro Footswitching Coagulation Forceps

Use with E360XXX reusable cord.



E700292

Total length: 25 cm (10 in.)
Tip: 1.2 mm, bayonet



E700331

Total length: 11 cm (4.5 in.)
Tip: 1.2 mm, straight



E700346

Total length: 16.5 cm (6.5 in.)
Tip: 2.2 mm, straight



E700349

Total length: 19 cm (7.5 in.)
Tip: 2.2 mm, straight

Product shown actual size

COVIDIEN, COVIDIEN with logo, Covidien logo and TM marked brands are internationally registered trademarks of Covidien AG. © 2010 Covidien.

contact www.covidien.com

Bipolar, micro Footswitching Coagulation Forceps

Use with E360XXX reusable cord.



E700352

Total length: 22 cm (8.75 in.)
Tip: 2.2 mm, straight



E700355

Total length: 25 cm (10 in.)
Tip: 2.2 mm, straight



E700376

Total length: 19 cm (7.5 in.)
Tip: 2.2 mm, angled

Bipolar, micro Footswitching Coagulation Forceps

Use with E360XXX reusable cord.



700760

Total length: 16.5 cm (6.5 in.)
Tip: 1.1 mm, selectal neuro, bayonet



700764

Total length: 20 cm (8 in.)
Tip: 1.1 mm, selectal neuro, bayonet



700768

Total length: 23 cm (9 in.)
Tip: 1.1 mm, selectal neuro, bayonet



700770

Total length: 25 cm (10 in.)
Tip: 1.1 mm, selectal neuro, bayonet

Product shown actual size

COVIDIEN, COVIDIEN with logo, Covidien logo and TM marked brands are internationally registered trademarks of Covidien AG. © 2010 Covidien.

contact www.covidien.com

Monopolar Coagulation Forceps

Use with E36018X reusable cord.



713050

Total length: 14.5 cm
Tip: 1 mm, Adson



713056

Total length: 18 cm
Tip: 2 mm, straight



713058

Total length: 21 cm
Tip: 2 mm, straight



713060

Total length: 24 cm
Tip: 2 mm, straight

Opti²™ Solid Laparoscopic Electrodes



E277036
Straight Spatula Tip
 Total length: 36 cm (14 in.)



E277136
Curved Spatula Tip
 Total length: 36 cm (14 in.)



E277236
Wire J Hook Tip
 Total length: 36 cm (14 in.)



E277336
Wire L Hook Tip
 Total length: 36 cm (14 in.)



E277436
Flat L Hook Tip
 Total length: 36 cm (14 in.)



E277045
Straight Spatula Tip
 Total length: 45 cm (18 in.)



E277345
Wire L Hook Tip
 Total length: 45 cm (18 in.)

Product shown actual size

COVIDIEN, COVIDIEN with logo, Covidien logo and ™ marked brands are internationally registered trademarks of Covidien AG. © 2010 Covidien.

contact www.covidien.com

SUCTION COAGULATORS

Handswitching, Monopolar Suction Coagulators Footswitching



E3310

E3305

Total length:

10 Fr ○
15.24 cm
(6 in.)

12 Fr ○
15.24 cm
(6 in.)



E2608-6

E2506-10FR

Total length:

6 Fr ○
15.24 cm
(6 in.)

10 Fr ○
15.24 cm
(6 in.)



ADAPTERS COAGULATION FORCEPS

Active and Specialty Adapters



E0017

Universal Monopolar Adapter

Accepts most standard active accessory cords up to 6 mm



E0507-B

Multiple Return/S Cord Adapter

Enables use of two patient return electrodes or an S-cord.

For use with all Valleylab REM™-equipped generators



E0504-2

E7506 Return Electrode Adapter

See catalog for a list of compatible generators

* Product no longer available

Product shown actual size

COVIDIEN, COVIDIEN with logo, Covidien logo and ™ marked brands are internationally registered trademarks of Covidien AG. © 2010 Covidien.

contact www.covidien.com

COVIDIEN, COVIDIEN with logo, Covidien logo and "positive results for life" are U.S. and internationally registered trademarks of Covidien AG.
Other brands are trademarks of a Covidien company. © 2011 Covidien. -V-VL-Accessory/FR - 05/2011



COVIDIEN (UK)
COMMERCIAL LTD

154 FAREHAM ROAD
GOSPORT HAMPSHIRE
PO 13 OAS
UNITED KINGDOM (UK)

+44 (0) 2030 271757 [T]
+44 (0) 2030 271758 [F]

WWW.COVIDIEN.COM