


Removal Stylet - Tillbehör till air-Q® larynxmask



air-Q® Removal Stylet & O.E.T.T. Guide

air-Q® storlek	air-Q® patientvikt	Removal Stylet storlek	Max O.E.T.T. storlek
0,5	<4 kg	00	4.0
1,0	4-7 kg	00	4.5
1,5	7-17 kg	00	5.0
2,0	17-30 kg	0	5.5
2,5	30-50 kg	0	6.5
3,5	50-70 kg	1	7.5
4,5	70-100 kg	1	8.5

Orderinformation - Removal Stylet

Artikelnummer	Benämning	Storlek	Passar airQ® strl	Förp
10-1006	Removal stylet enpatientmodell, Small	00	0,5 - 1,5	10 st/pkt
10-1005	Removal stylet enpatientmodell, Medium	0	2,0 & 2,5	10st/pkt
10-1004	Removal stylet enpatientmodell, Large	1	3,5 & 4,5	10st/pkt