



Safety Winged Infusion Set

Product Information Card




Description:

The Surflo® Safety Winged Infusion Set with Protector (= Scalp Vein or SV-set) is intended to introduce or withdraw fluids from the vascular system after vein puncture.

The Surshield® protector is intended to cover the exposed needle immediately after use to prevent inadvertent pricking with the contaminated needle.

Does not contain components made of DEHP or natural rubber latex.

Components:



Description	Raw Materials
1. Surshield protector	Polypropylene
2. Cannula	Stainless steel
3. Wing hub	Polyvinyl chloride
4. Tube	Polyvinyl chloride
5. Adapter	PMMA
6. Adapter cap	Polypropylene
Safety hinge	ABS
Protector	Polyethylene



Product Range:

Product Code	Colour	External Cannula Ø		Cannula Length	Dead Space Volume	Tubing Length
SV*S18NL09	Pink	18G	1.20 mm	19 mm	0.21 ml	90 mm
SV*S19NL09	Yellow	19G	1.10 mm	19 mm	0.21 ml	90 mm
SV*S21NL09	Green	21G	0.80 mm	19 mm	0.20 ml	90 mm
SV*S22NL09	Black	22G	0.70 mm	19 mm	0.20 ml	90 mm
SV*S23NL09	Light Blue	23G	0.60 mm	19 mm	0.20 ml	90 mm
SV*S25NL09	Orange	25G	0.50 mm	19 mm	0.19 ml	90 mm
SV*S18NL30	Pink	18G	1.20 mm	19 mm	0.45 ml	300 mm
SV*S19NL30	Yellow	19G	1.10 mm	19 mm	0.45 ml	300 mm
SV*S21NL30	Green	21G	0.80 mm	19 mm	0.41 ml	300 mm
SV*S22NL30	Black	22G	0.70 mm	19 mm	0.41 ml	300 mm
SV*S23NL30	Light Blue	23G	0.60 mm	19 mm	0.40 ml	300 mm
SV*S25NL30	Orange	25G	0.50 mm	19 mm	0.40 ml	300 mm

Safety Winged Infusion Set



Packaging: *Single unit:* blister lid paper / thermoformable blister film

Product Code	Unit Box	Shipping Carton		
		Quantity	Gross Weight	Dimensions (LxWxH)
SV*S18NL09 / SV*S19NL09 SV*S21NL09 / SV*S22NL09 SV*S23NL09 / SV*S25NL09	50	500	2.9 kg	510 × 257 × 260 mm
SV*S18NL30 / SV*S19NL30	50	500	3.2 kg	510 × 257 × 260 mm
SV*S21NL30 / SV*S22NL30 SV*S23NL30 / SV*S25NL30	50	500	3.1 kg	510 × 257 × 260 mm

Sterilization: With ethylene oxide to a Sterility Assurance Level of at least 10⁻⁶

Storage: Shelf life of 3 years
** The batch number and expiry date are printed on the unit pack*
** Do not store at extreme temperature and humidity*

Precautions for Use: Do not use if package is damaged
 For single use only
 Dispose of safely after single use to avoid risk of infection

Status of Regulation: Class IIa, sterile product following MDD 93/42/EEC

Notified body: TÜV SÜD 0123
 Manufacturer: TERUMO MEDICAL PRODUCTS HANGZHOU (China)
 M4-9-5, Economic & Technological Development Zone
 Hangzhou, 310018
 China
 EC representative: TERUMO EUROPE N.V.
 Interleuvenlaan 40
 B - 3001 Leuven

HO 168GB-10.07TE-1107.16E



Terumo Corporation
 +81 3 3374 8111

Terumo Europe NV
 +32 16 38 12 11

Terumo Medical Products EMEA
 +32 16 38 12 11

EMEA Sales Offices

Terumo Europe NV African Business Division
 +32 16 38 13 08

Terumo Europe NV Benelux Sales Division Belgium
 +32 16 339 25 80
 The Netherlands
 0800 0220396

Terumo Europe NV Emerging Markets Division
 +32 16 38 12 11

Terumo Deutschland GmbH
 +49 6196 80 230

Terumo Deutschland GmbH Switzerland
 +41 56 419 10 10

Terumo Europe España SL
 +34 902 10 12 98

Terumo France S.A.S.
 +33 130 96 13 00

Terumo Italia S.R.L.
 +39 0651 96 14 20

Terumo Russia LLC
 +7 495 988 4740

Terumo Sweden AB
 +46 3174 85 880

Terumo Sweden AB Denmark
 +45 7020 93 80

Terumo Middle East FZE
 +971 4 292 0200

Terumo UK Ltd
 +44 1276 480 440

Terumo BCT Tibbi Cihazlar Dagim ve Hizmetleri A.S.
 +90 216 590 0074

