



Technical Data Sheet

<b>Company Name:</b>	<b>Local company name</b>	<b>Telephone</b>
Cook Medical Europe LTD.	Cook Österreich GMBH	+43 17 956 7121
	Cook Belgium NV/SA	+32 27001633
<b>Address:</b>  O' Halloran Road National Technology Park Limerick, Ireland	Cook Danmark ApS	+45 38 48 76 07
	Cook Medical Distribution Partners	+353 61 239 240
	Cook Medical Finland Oy	+358 9725 19996
	Cook Deutschland GMBH	+49 69 500 72 804
	Cook Medical Hungary KFT.	+36 1 777 9199
	Cook Medical Europe LTD.	+353 61 239 252
	Cook Nederland BV	+31 202013367
	Cook Norge AS	+47 23162968
	Cook Medical SP. Z.O.O.	+48 22 3060159
	Cook Sweden AB	+46 85 876 94 68
	Cook Switzerland AG	+41 44 800 96 09
<b>Web:</b> <a href="http://www.cookmedical.com">www.cookmedical.com</a>	Cook UK LIMITED	+44 207 3654 183

<b>Vigilance Contact:</b>  Sinead Quaid	<b>Tel:</b> +353 61 33 44 40  <b>Email:</b> <a href="mailto:sinead.quaid@cookmedical.com">sinead.quaid@cookmedical.com</a>
---	---

**Product:** Williams Cystoscopic Injection Needle

**Product Code:**

090001  
090001-S26  
090001-S34

## **DIMENSIONS**

### **090001**

Needle

- Gauge: 23 G
- Tip Length: 8 mm

Sheath

- Diameter: 5.0 Fr.
- Length: 35 cm

### **090001-S26**

Needle

- Gauge: 23 G
- Tip Length: 8 mm

Sheath

- Diameter: 5.0 Fr.
- Length: 45 cm


### **090001-S34**

Needle

- Gauge: 23 G
- Tip Length: 8 mm

Sheath

- Diameter: 3.7 Fr.
- Length: 35 cm

**Picture****Components:**

Needle  
Sheath

**Description:**

Radiopaque polyurethane sheath with a stainless steel needle cannula.

**Intended use:**

Used for cystoscopic injection into the urethra, bladder neck, and bladder wall.

**Contraindications:**

None known.

**Precautions:**

None known.

**Recommended storage instructions:**

Store in a clean and dry location, away from temperature extremes.

Upon removal from package, inspect the product to ensure no damage has occurred.

<b>Patient contacting materials</b>	Stainless steel Radiopaque polyurethane
-------------------------------------	--

<b>CE Marked</b>	Yes
<b>Class</b>	IIa
<b>Notified body number</b>	0088
<b>Manufacturer</b>	Cook Incorporated

<b>Single use</b>	Yes
<b>Sterilized</b>	Yes
<b>Address of sterilization</b>	6300 North Matthews Drive, Ellettsville, Indiana 47429 USA
<b>Method of sterilization</b>	EtO
<b>Shelf life</b>	Up to 36 months
<b>Biological origin</b>	No
<b>Latex</b>	<b>No</b>
<b>Biocompatibility</b>	Yes
<b>Quantity supplied</b>	1