



The impregnated gauze dressing for superficial, exuding wounds



CUTICERIN° impregnated gauze dressing for superficial, exuding wounds

The impregnated gauze, CUTICERIN is made of a smooth acetate¹ fabric with a defined mesh size² and water repellent³ ointment coating. This reduces the risk of granulation tissue growing through the dressing4.

Advantages

CUTICERIN is low adherent^{4,5} and may reduce pain on dressing change⁴. It is made from a smooth acetate fabric¹ with a defined mesh size² and water repellent³ ointment coating.

CUTICERIN minimises pain at dressing changes due to the reduced risk of granulation tissue growing through the dressing allowing the dressing to be easily removed4.

CUTICERIN is an impregnated gauze dressing which allows air and wound exudate to pass freely through the dressing6.


Presentation

CUTICERIN is sterile and individually sealed

Indications

For the management of superficial, exuding wounds, such as;

- Burns
- Abrasions
- Split-thickness skin graft donor sites





References

- 1. Laboratory Report Reference AA400168 dated 26.07.96 Claim Smooth Gauze.
- 2. Laboratory Report Reference AA400251 dated 18.09.96 Claim Defined Mesh Size.
- 3. Laboratory Report Reference AA400167 dated 15.07.96 Claim Hydrophobic Ointment Coating.
- 4. Petres J (1988) Moist Wound Dressing in Superficial Skin Defects Zeuschrift fur Haut krankheithen 63 (Supplement 2) 25-29.
- 5. Waymack et al., (1986), An Evaluation of Aquaphor Gauze Dressing in burned children. Burns (1986) 12. 443-448.
- 6. Laboratory Report Reference AA400166 dated 26.07.96 Claim Air Permeability and Secretion Permeability.

www.smith-nephew.com/wound

°Trademark of Smith & Nephew

©May 2016 Smith & Nephew

Ordering Information

Code	Size	Dressings per inner/outer
66045560	7.5cm x 7.5cm	50/1000
66045561	7.5cm x 20.0cm	50/1000
66045562	7.5cm x 7.5cm	10/200
66045563	7.5cm x 20.0cm	10/200
66045502	20.0cm x 40.0cm	25/150
66045503	7.5cm x 20cm (3s)	25/3
66045564	10cm x 10cm	10

Wound Management Smith & Nephew Medical Ltd 101 Hessle Road Hull HU3 2BN

F +44 (0)1482 328326

T +44 (0)1482 225181