



Peel-Away® Sheath Introducer

Used for percutaneous placement of external and internal drainage catheters that are constructed from softer materials.

- The dual knobs allow the sheath to be peeled away and removed.
- The maximum diameter of the instrument or catheter should be determined before its introduction to ensure its passage through the sheath.

Components



| Order Number | Reference Part Number | Introducer Fr | Accepts Wire Guide Diameter inch | Sheath Length cm | Dilator Length cm |
|--------------|-----------------------|---------------|----------------------------------|------------------|-------------------|
| G04897 | PLVW-6.0-35 | 6 | .035 | 15.5 | 20 |
| G00925 | PLVW-7.0-35 | 7 | .035 | 15.5 | 20 |
| G00930 | PLVW-8.0-38 | 8 | .038 | 15.5 | 20 |
| G01697 | PLVW-9.0-38 | 9 | .038 | 15.5 | 20 |
| G01980 | PLVW-10.0-38 | 10 | .038 | 15.5 | 20 |
| G06491 | PLVW-11.0-38 | 11 | .038 | 15.5 | 20 |
| G06492 | PLVW-12.0-38 | 12 | .038 | 15.5 | 20 |
| G06493 | PLVW-13.0-38 | 13 | .038 | 15.5 | 20 |
| G06496 | PLVW-14.0-38 | 14 | .038 | 15.5 | 20 |
| G06497 | PLVW-15.0-38* | 15 | .038 | 15.5 | 20 |
| G02997 | PLVW-16.0-38 | 16 | .038 | 15.5 | 20 |
| G06444 | PLVW-18.0-38 | 18 | .038 | 15.5 | 20 |
| G06445 | PLVW-20.0-38 | 20 | .038 | 13 | 20 |
| G04096 | PLVW-22.0-38 | 22 | .038 | 13 | 20 |
| G06446 | PLVW-24.0-38 | 24 | .038 | 13 | 20 |
| G06447 | PLVW-26.0-38 | 26 | .038 | 13 | 20 |

*Nonstock

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.