

SAPPHIRE
INFUSION PLATFORM
CANADA CATALOGUE 2020




SAPPHIRE
by Eitan Group

CONTENTS

●	<u>INFUSION KITS</u>	4
●	<u>ACCESSORIES</u>	6
●	<u>ADMINISTRATION SETS</u>	10
●	<u>PC APPLICATIONS</u>	13

Sapphire™
Multi-therapy infusion pump kit




Includes:

- Sapphire™ Multi-therapy pump
- Power supply
- Mini cradle
- Mini cradle knob key
- User manual

P/N 17000-028-0075

Sapphire™
Epidural infusion pump kit



Includes:

- Sapphire™ Epidural pump
- Power supply
- Bolus handle
- Mini cradle with splitter
- Mini cradle knob key
- User manual

P/N 17000-031-0036

Sapphire™
Mini Cradle

Sapphire™
Mini Cradle

With Splitter

P/N 17000-026-0002

Without Splitter

P/N 05040-300-0011



Sapphire™
Multi Pump Mounting System



P/N 17000-020-0008

Sapphire™
Splitter



P/N 05040-300-0012

Sapphire™
Integrated Power Supply



P/N 15072-000-0010

Sapphire™
Power Supply



P/N 15072-000-0007

Sapphire™
AC Power Cord - Type B



P/N 15070-112-0005

Sapphire™
External Battery Pack (EBP)



P/N 15079-111-0001

Sapphire™
Power Cord - compatible
with IPS (US and Canada)



P/N 15171-112-0002

Sapphire™
Infusion Pouch 500 mL



P/N 05080-605-0001

Sapphire™
Travel Case



P/N 15020-010-0002

Sapphire™ Homecare Large Backpack



P/N 15088-000-0001

Sapphire™
PCA lockbox 250mL



P/N 15063-000-0001

Sapphire™
PCA lockbox 100mL



P/N 15061-131-0004

Sapphire™ Bolus handle



P/N 15071-000-0005






ADMINISTRATION SETS | STANDARD SETS

NAME	P/N	DESCRIPTION	LENGTH (CM)	PRIME VOLUME	UNITS/ BOX
AP 403	12003-000-0005	Infusion set with vented/non-vented drip chamber and needle y-site	290	~6 mL	30
AP 406	12000-000-0007	Microbore with non vented spike	287	~24 mL	30
AP 407	12000-000-0008	Macro bore with vented/non vented drip chamber	286	~24 mL	20
AP 409	12005-000-0001	Primary infusion set with vented/non-vented drip chamber, back check valve and 2 needless y-site	277	~23 mL	30
AP 413	12003-000-0007	Microbore with needless Y-site, 2 spin male-luer locks	289	~6 mL	30
AP 423	12003-000-0017	Microbore with non-vented spike and a back check valve	292	~7 mL	30

ADMINISTRATION SETS | FILTER SETS

NAME	P/N	DESCRIPTION	LENGTH (CM)	PRIME VOLUME	UNITS/ BOX
AP 201	12000-000-0006	Infusion set with non-vented spike and 0.2µm filter	291	~20 mL	20
AP 204	12003-000-0004	Microbore with non-vented spike, 0.2µm filter and back check valve	295	~8 mL	30
AP 205	12000-000-0004	Infusion set with non-vented spike, 1.2µm filter and needle y-site	297	~23 mL	20
AP 206	12000-000-0005	Macro bore with non vented spike, 1.2µ filter and needless Y-site	295	~23 mL	20
AP 207	12000-000-0009	Macro bore with vented/non vented drip chamber, 1.2µ filter and needle Y-site	291	~26 mL	30
AP 210	12003-000-0010	Microbore with non-vented spike and 1.2µ filter	304	~9 mL	30

ADMINISTRATION SETS | SPECIALTY SETS

NAME	P/N	DESCRIPTION	LENGTH (CM)	PRIME VOLUME	UNITS/BOX
AP 101	12006-000-0002	Blood set with non-vented spike, 200µ filter and needleless Y-site	292	~28 mL	20
					
AP 402	12004-000-0001	Epidural (yellow marked) microbore with non-vented spike	290	~6 mL	30
					
AP 410	12005-000-0002	Secondary extension set with hook	Set: 71Hook: 20	~9 mL	30
					
AP 424	12003-000-0015	PCA Microbore with non-vented spike, side arm with female luer lock and back check valve	294	~7 mL	30
					
AP 504	12000-000-0016	UV-Resistant macrobore with vented/non-vented drip chamber and needleless Y-site	287	~26 mL	20
					

PC APPLICATIONS

The Sapphire modern, software based platform expands beyond the pump to the PC applications suite that complements the pumps;



PUMP LOADER

Enables upgrading or reloading software to the Sapphire pump via a wire connection.



DRUG LIBRARY

Enables to set configuration settings, drug specific limit and more. When complete, the library can be uploaded to multiple Sapphire pumps.



EVENT VIEWER

Enables downloading and viewing the event log of the Sapphire infusion pump on a computer.






ANNUAL CERTIFICATION KIT










(SW + disposable kit [15075-020-0001])
Enables to run performance test (air detection, pressure detection and accuracy). Results are compared to manufacturing database per pump.



ACCESSORIES | PC APPLICATIONS

PC APPLICATION	ITEM	P/N	PICTURE
Annual Certification Test	Annual certification kit	15075-020-0001	
All	Communication cable	05020-110-0213	
All	RS232 to USB 4 port adapter	15077-000-0001	

SPARE PARTS | ACCESSORIES

SPARE PART FOR	ITEM	P/N	PICTURE
Sapphire pump	Li-Ion battery	05020-160-0001	
Mini cradle	connection knob flange Key	05040-380-0019	
Mini cradle	connection knob key	05040-380-0018	
Mini cradle	MC- AVP Connection Bolt Knob Tip	05040-380-0027	
Mini cradle	MC- AVP Press Plane	05040-380-0026	
PCA lockbox 100	Key	05070-190-0015	
PCA lockbox 250	Key	10010-122-0046	
PCA lockbox 250	Release clip (eliminates lockbox locking onto cradle)	15063-121-0029	
PCA Lockbox 500	Key	15062-132-0001	

REF 10050-003-0010 | Rev.01 / 03.2020

Q Core Medical Ltd.
29 Yad Haruzim St., P.O. Box 8639
Netanya 4250529, Israel
Tel: +972.73.2388888
Fax: +972.73.2388800
Email: support@Qcore.com
www.qcore.com

U.S. & Canada | Contact information:
65 Enterprise 3rd Floor
Aliso Viejo, CA 92656
Tel: 877-541-9944 | 949.303.7845
Fax: 949-288-5482
Email: customerservice@eitangroup.com
www.eitangroup.com



SAPPHIRE
by Eitan Group