



Återförsäljare:

MEDI PLAST

Mediplast AB
 Box 1004, 212 10 Malmö
 T 020-78 80 35
 info.se@mediplast.com
 www.mediplast.com



Technical sheet

Surgical clipper with single use Ergosafe blades

Reference(s)

The surgical clipper is used for hair removal from the operative site prior to surgery.
 DYND70800EU – Clipper handle and charger, EU plug
 DYND70800UK – Clipper handle and charger, UK plug
 DYND70801EU – Clipper handle
 DYND70802EU – Charger, EU plug
 DYND70802UK – Charger, UK plug

Description

Blade attachments references	CUB1 – Universal Ergosafe blade	CNB2 – Narrow Ergosafe blade	CCB3 – Coarse Ergosafe blade
Recommended use	General purpose for body hair	Urology – Obstetrics Gynecology	Neurosurgery, Thoracic surgery
Areas	Body	Perinea and sensitive areas	Scalp and other thick coarse hair
Clipper charging time	12 hours (Electromagnetic, NiMH battery)		
Clipper run time on a fully charged battery	90 – 120 minutes uninterrupted running		
Submergibility	Submersible in water no deeper than 1 meter for up to 30 min.		
Cleaning/disinfections	Soap and water - Disinfecting swab / antimicrobial soap Rinse with hot or cold water then dry		

Storage

Temperature : -10°C to +60°C
 Humidity: 10% to 100% RH
 Pressure: 500hPa – 1060hPa

Main materials

Clipper body: Polyoxymethylene
Clipper charger: ABS
Blades: Polyoxymethylene and Stainless Steel
Products do not contain phthalates. Latex free.

Pack factor

Blades: 1ea / blister - 50 ea/ box
Clipper : 1ea / box

Manufacturer

Medline Industries LP, Three Lakes Drive, Northfield, IL 60093 - EN ISO 13485

Regulation

Compliant to Medical Device Regulation (EU) 2017/745 and IEC 60601-1, IPX7
Class : I
Please refer to complete labeling and instruction for use before using the product.

