

ETHICON

a Johnson & Johnson company



MONOCRYL*
(Poliglecaprone 25) Suture

MONOCRYL*

(Poliglecaprone 25) Monofilament

SUTURE

Description: **Monofilament** For smooth passage through tissue.
Synthetic For minimal tissue reaction.
Absorbable By hydrolysis. Essentially complete in 90 -120 days.

Tensile strength: **Post implantation** **Approximate % original strength remaining**

7 days	60% (violet) 50% (undyed)
14 days	30% (violet) 20% (undyed)
21 days	0% (undyed)
28 days	0% (violet)

Color: Violet or undyed.

Range: 6-0 to 1 (USP). Supplied as needled sutures and ligatures.

Indications: MONOCRYL* sutures are intended for use in general soft tissue approximation and/or ligation where an absorbable material is indicated. Typical areas of use include:

- Subcuticular sutures
- Small intestine anastomoses
- Urological anastomoses

Full details in the Instruction for Use included in every package.

