

Disposables – the Original by B. Braun

Flexibility and Accuracy for All Pumps



Automated Infusion Systems

Infusomat® Space Lines

A System Dedicated to Flexibility for All Therapies and All Wards

Easy Start

Easy and secure container connection due to special spike design



Design Assisted Loading

- Long term accuracy up to 96 h
- Anti-freeflow clamp for automated freeflow protection when removing the line from the pump

Find out more online:

1. Test Result IV-Study
2. HemoIyse Test Result
3. Test Result Infection Control

Either scan the 2D barcode with your smartphone or follow the link to find out more
www.space.bbraun.com



Clinically proven
Bio-compatibility¹



**Infusomat® Space Line
Type Transfusion?**
Same hemolytic effect as gravity application



**Eurofix®
Injection-Port**
Ergonomic plate
to avoid needle-
stick injuries

- Safeflow®³**
- Needle-free Injection Port
 - Swabable – state of the art in infection control
 - Minimal priming volume (0.08 ml)



100 % light protection



DEHP-free

Free of leachables



PVC-free

User Benefits

Golden Standard:

- ① High long term accuracy
- ① Active set based freeflow protection
- ① Clinically proven materials and components


Therapy diversity:

- ① Different tubing materials
- ① High specialized components offering enhanced safety
- ① Lines with therapy specific configuration

Original Perfusor® Syringes and Lines

A System Dedicated to Accuracy

100 % light protection



Original Perfusor® and Omnifix® Syringes

- Excellent start-up properties
- Optimal gliding characteristics
- Minimal residual volume
- Available in all sizes from 3 ml to 50/60 ml




Find out more online:

1. Test Result IV-Study

Either scan the 2D barcode with your smart-phone or follow the link to find out more

www.space.bbraun.com





PVC-free

DEHP-free

DEHP-free

DEHP-free

Original Perfusor® Line PCA

- With BCV (Back-Check Valve) to prevent reflux of drugs
- Slide clamps for convenient flow stop

Original Perfusor® Lines

- Available in different lengths and materials
- Minimal residual volume
- Clinically proven Bio-compatibility¹

Original Perfusor® Line with Safsite® valve

- Automates freeflow protection
- Clinically proven Bio-compatibility¹

User Benefits

- ➊ Original Perfusor® Syringes & Lines set the benchmark for accuracy and long term consistency in infusion therapy
- ➋ A Landmark for High Performance

Infusomat® Space Line Integrates ...

Ready to Use Solutions



Multiconnector compatible with most ready to use containers



PVC-free


Y-Port for additional drug applications

Dedicated Enteral Patient connector

- Male Luer Lock avoids mix up with IV access ports
- Scaled cone adapter for gastro tubes

... Enteral Nutrition in Multifunctional System

Patient Tailored Nutrition Solutions



Infusomat® Space Line for Enteral Nutrition
With integrated nutrition bag for easy admixture

Y-Port with integrated back-check valve for additional drug applications

User Benefits


- ➊ Maximal flexibility with universal connector and integrated nutrition bag
- ➋ Enhanced safety with dedicated enteral nutrition patient connector
- ➌ Integration of EN in a multifunctional system

Dedicated Solutions for Smallest Patients



Dosifix[®] Burette

- For safe and convenient drug admixture
- Needle-free Injection Port





Find out more online:

3. Test Result Infection Control
Either scan the 2D barcode with your smart-
phone or follow the link to find out more
www.space.bbraun.com



Safeflow³

- Needle-free Injection Port
- Swabable – state of the art in infection control
- Minimal priming volume (0.08 ml)



Original Perfusor® Line PE

- Minimal residual volume
- Different sizes available from 50 - 300 cm
- Leachable-free

Original Perfusor® Line with Sterifix

- Particle and air filter (0.2 µm pore size)
- Minimal residual volume

Omnifix®

- Small syringes for short alarm reaction time
- Smooth start-up
- Long term consistency



User Benefits

- 1 Always DEHP-free
- 2 Components and tubes have minimal priming volume
- 3 Needle-free connectors
- 4 Original Perfusor® Syringes in small sizes allow short alarm reaction time

CytoSet Infusomat® Space

Oncology – We Make Things Safe So You Can Have It Easy

Cyto-Set® Line

- For preparation on ward
- PrimeStop for automated priming

Acoustic feedback indicates safe connection

Safe Preparation

Acoustic feedback indicates safe connection

Cyto-Set® Line

- With integrated back-check valve and 0.2 µm filter

Acoustic feedback indicates safe connection

Find out more online:





Either scan the 2D barcode with your smart-phone or follow the link to find out more
www.space.bbraun.com







Product Specifications

Original Infusomat® Space Lines Type IV-Standard



Infusomat® Space Lines Type IV-Standard	Add on	PVC (DEHP-free)	PVC-free	Length (cm)	Sales Units / Pcs.	Code No. (REF)
		■	☐	250	100	8700036 SP
	Ward pack	■	☐	250	100	8700435 SP
	Extra long	■	☐	300	100	8270350 SP
	Eurofix Y-Port for needle access	■	☐	270	100	8700087 SP
		☐	■	250	100	8250731 SP
	Safeflow injection port	☐	■	300	100	8700110 SP

Therapy Specific Variants



Infusomat® Space Lines Type piggy-back Neutrapur®	Add on	PVC (DEHP-free)	PVC-free	Length (cm)	Sales Units / Pcs.	Code No. (REF)
	Primary line with BCV and Safeflow injection port	-	■	300	25	8250710 SP
	Secondary line patient connector with integrated BCV	-	■	75	25	4062877
Infusomat® Space Lines Type Transfusion	Add on	PVC (DEHP-free)	PVC-free	Length (cm)	Sales Units / Pcs.	Code No. (REF)
	Drip chamber with integrated blood filter (200 µm pore size)	■	-	250	100	8270066 SP
Infusomat® Space Lines UV-Protected	Add on	PVC (DEHP-free)	PVC-free	Length (cm)	Sales Units / Pcs.	Code No. (REF)
	Opaque tube material	■	-	250	100	8700125 SP
	Opaque tube material with needle based injection port	■	-	250	100	8250430 SP

Product Specifications


Original Infusomat® Space Lines Type Enteral Nutrition


Infusomat® Space Lines Type Enteral Nutrition Neutrapur®	Add on	PVC (DEHP-free)	PVC-free	Length (cm)	Sales Units / Pcs.	Code No. (REF)
	With multiconnector, needle-free injection port, female Luer Lock patient connector and scaled cone adapter	-	■	320	40	8250857 SP
Infusomat® Space Lines Type Enteral Nutrition	Add on	PVC (DEHP-free)	PVC-free	Length (cm)	Sales Units / Pcs.	Code No. (REF)
	With integrated nutrition bag, needle-free injection port, female Luer Lock patient connector and scaled cone adapter	■	-	230	25	8250839 SP

Original Infusomat® Space Lines Type Neonatology

Infusomat® Space Lines Type IV-Standard	Add on	PVC (DEHP-free)	PVC-free	Length (cm)	Sales Units / Pcs.	Code No. (REF)
	With Dosifix® Burette 150 ml and needle free safeflow Y-Port	-	■	270	40	8250245 SP
	With Dosifix® Burette 150 ml and needle based Y-Port	■	-	270	25	8250294 SP

Oncology Application

Infusomat Space Lines Type Oncology	Add on	PVC (DEHP-free)	PVC-free	Length (cm)	Sales Units / Pcs.	Code No. (REF)
	With Sterifix®, particle and air filter 0,2 µm	—	■	250	20	8700095 SP

Cyto-Set® Infusomat Space	Add on	PVC (DEHP-free)	PVC-free	Sales Units / Pcs.	Code No. (REF)
	with 3 Safsite®	—	■	20	8250910 SP
	with 5 Safsite®	—	■	20	8250812 SP
	with 5 Safsite® and 0.2 µm filter	—	■	20	8250413 SP
	with 3 Safsite® and BCV	—	■	20	8250913 SP
	with 3 Safsite® and BCV, Opaque	■	—	20	8250915 SP
	with 5 Safsite® and BCV	—	■	20	8250813 SP
	with 5 Safsite® and BCV, Opaque	■	—	20	8250818 SP

Preparation - Cyto-Set® Mix







Cyto-Set® Mix	Add on	PVC (DEHP-free)	PVC-free	Sales Units / Pcs.	Code No. (REF)
	with 1 Safsite®	—	■	20	A 2901 N
	with 1 Safsite® and BCV	—	■	20	A 2902 N
	with 1 Safsite® and BCV, Opaque	■	—	20	A 2904 N
	with 1 Safsite®, with filter and BCV	—	■	20	A 29023 N

Preparation - Cyto-Set® Line

Cyto-Set® Line	Add on	PVC (DEHP-free)	PVC-free	Sales Units / Pcs.	Code No. (REF)
	without Safsite®	—	■	20	A 2581 NF
	without Safsite®, with BCV	—	■	20	A 25811 NF
	with 1 Safsite®, with filter and BCV	—	■	20	A 25812 NF




Product Specifications

Omnifix® Syringes

Omnifix® Syringes	Volume	Cone	Graduation	Style of syringe cone	Sales Units / Pcs.	Code No. (REF)
	3 ml	centric	0.1 ml	Luer-Lock	2500 (25 x 100)	4617022 V
	5 ml	centric	0.2 ml	Luer-Lock	2000 (20 x 100)	4617053 V
	10 / 12 ml	centric	0.5 ml	Luer-Lock	1200 (12 x 100)	4617100 V
	20 ml	centric	1.0 ml	Luer-Lock	800 (8 x 100)	4617207 V
	30 ml	centric	1.0 ml	Luer-Lock	600 (6 x 100)	4617304 F
	50 / 60 ml	centric	1.0 ml	Luer-Lock	100 (1 x 100)	4617509 F



Original Perfusor® Syringes

Original-Perfusor® Syringe 20 ml	Needle outer ø	Needle length	Volume	Cone	Max. Volume	Gradua- tion	Style of syringe cone	Sales Units / Pcs.	Code No. (REF)
			20 ml	centric	20 ml	1 ml	Luer-Lock	100	8728615
	2 mm	30 mm	20 ml	centric	20 ml	1 ml	Luer-Lock	100	8728623
Original-Perfusor® Syringe 50 ml	Needle outer ø	Needle length	Volume	Cone	Max. Volume	Gradua- tion	Style of syringe cone	Sales Units / Pcs.	Code No. (REF)
			50 ml	centric	50 ml	1 ml	Luer-Lock	100	8728844 F
	2 mm	30 mm	50 ml	centric	50 ml	1 ml	Luer-Lock	100	8728810 F
	with filter								
	2 mm	28 mm	50 ml	centric	50 ml	1 ml	Luer-Lock	100	8728852 F
Original-Perfusor® Syringe Light protection	Needle outer ø	Needle length	Volume	Cone	Max. Volume	Gradua- tion	Style of syringe cone	Sales Units / Pcs.	Code No. (REF)
	Opaque with filter								
	2 mm	28 mm	50 ml	centric	50 ml	1 ml	Luer-Lock	100	8728828 F






Product Specifications

Original Perfusor® Lines

Original Perfusor® Lines	PVC (DEHP-free)	Inner diameter (mm)	Length (cm)	Sales Units / Pcs.	Code No. (REF)
	▪	1.5	300	100	8255253
	▪	1.5	250	100	8255490
	▪	1.5	200	100	8722862
	▪	1.5	150	100	8722960
	▪	1.5	50	100	8255172
Original Perfusor® Line	PVC (DEHP-free)	Inner diameter (mm)	Length (cm)	Sales Units / Pcs.	Code No. (REF)
	▪	0.9	75	100	8722870
	▪	0.9	150	100	8255504
Original Perfusor® Line with Sterifix® 0.2 µm particle & air filter	PVC (DEHP-free)	Inner diameter (mm)	Length (cm)	Sales Units / Pcs.	Code No. (REF)
	▪	1.5	200	100	8723001

Original Perfusor® Lines

Original Perfusor® Line Type PCA with BCV and slide clamp	PVC (DEHP-free)	Inner diameter (mm)	Length (cm)	Sales Units / Pcs.	Code No. (REF)
	▪	1.5	150	100	8726019
Original Perfusor® Lines	PVC-Free (PE)	Inner diameter (mm)	Length (cm)	Sales Units / Pcs.	Code No. (REF)
	▪	0.9	50	100	8255059
	▪	0.9	100	100	8255067
	▪	0.9	300	100	8250146
	▪	0.9	250	100	8272565
	▪	0.9	200	100	8723060
	▪	0.9	150	100	8722935
Original Perfusor® Line	PVC-Free (PE)	Inner diameter (mm)	Length (cm)	Sales Units / Pcs.	Code No. (REF)
	▪	0.9	150	100	8723010

