



D-fend™ Water Traps

D-fend water traps help protect your GE Healthcare respiratory monitoring equipment from humidity, secretions, viral and bacterial contamination. The D-fend water trap is indicated for use with respiratory modules in both the anesthesia and critical care environments.

Features

- Suitable for GE Healthcare respiratory modules for accurate measurement
- Tested by an independent third-party lab for viral filtration efficiency (VFE) against a surrogate infectious material known as Phi X 174.

PRODUCT COMPATIBILITY	Anesthesia		Critical Care		Anesthesia or Critical Care
	D-fend 876446-HEL	D-fend Pro M1182629	D-fend+ 881319-HEL	D-fend Pro+ M1200227	Mini D-fend 8001274
E-CAiO	X				
E-CAiOV	X				
E-CAiOVX	X				
E-sCAiO		X			
E-sCAiOV		X			
E-sCAiOVX		X			
E-sCAiOE		X			
E-sCAiOVE		X			
N-CAiO		X			
E-CO			X		
E-COV			X		
E-COVX			X		
E-sCO				X	
E-sCOV				X	
E-sCOVX				X	
E-miniC					X

Our equipment.
Our accessories.
Simply compatible.

You can trust GE Healthcare accessories and consumables to keep your clinical devices working efficiently.

- Simplified ordering process
- Fast delivery
- One contact point for expert advice
- Quality support at every step

PRODUCT SPECIFICATIONS

Product use

D-fend water traps are disposable.

D-fend+ (881319-HEL),
D-fend Pro+ (M1200227): Critical Care, Maximum 24 hours, single-patient-use

D-fend (876446-HEL),
D-fend Pro (M1182629): Anesthesia, Maximum 2 months.

Mini D-fend (8002174): Anesthesia, Maximum 2 months. Critical Care, Maximum 24 hours, single-patient-use

Materials

	D-fend Pro, D-fend Pro+	D-fend, D-fend+, Mini D-fend
Body	Acrylonitrile Butadiene Styrene (ABS) and Polycarbonate (PC)	Polyamide (PA)
Top Part	Polycarbonate (PC) Thermoplastic Elastomer (TPE)	N/A
Gasket	N/A	Silicone rubber
Container	Polycarbonate (PC) Thermoplastic Elastomer (TPE)	Polycarbonate (PC)
O-Ring	N/A	Nitrile rubber
Latch	Polyoxymethylene (POM)	N/A
Filter	Polytetrafluoroethylene (PTFE)	GORE microfiltration media

Materials composition

Not made with natural rubber latex, PVC, DEHP or other phthalates.

Product may not be available in all countries and regions. Contact a GE Healthcare Representative for more information. Please visit www.gehealthcare.com.

© 2020 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram and D-fend are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

DOC2390764 10/20

Maintenance

Do not wash, dry, sterilize or reuse exhausted/occluded water traps. When cleaning the outer surface of the module, keep the water trap and gas sampling line connected to prevent cleaning agent from entering the water trap or gas module.

Disposal considerations

Dispose as clinical waste, in accordance with hospital policy, local guidelines and regulations.

Packaging materials

Primary: Polyethylene (PE)

Compliance

CE Mark

Bacterial and Viral Filtration Efficiency

Part	BFE	VFE
D-fend D-fend+ Mini D-fend	>99.99998%	99.99997 - 99.99999%
D-fend Pro D-fend Pro+	>99.999970%	>99.9988%

PART LIST

Part Number	Description	QTY
876446-HEL	D-fend Water Trap, Black	10
M1182629	D-fend Pro Water Trap, Dark Steel Blue	10
881319-HEL	D-fend+ Water Trap, Green	10
M1200227	D-fend Pro+ Water Trap, Green	10
8002174	Mini D-fend Water Trap	10

Not available in all countries