

# MONOCRYL\*

(Poliglecaprone 25) Monofilament

## SUTURE

**Description:** **Monofilament** For smooth passage through tissue.  
**Synthetic** For minimal tissue reaction.  
**Absorbable** By hydrolysis. Essentially complete in 90 -120 days.

**Tensile strength:** **Post implantation** **Approximate % original**  
**7 days** **length remaining**  
60% (violet)  
50% (undyed)

14 days 30% (violet)  
20% (undyed)

21 days 0% (undyed)

28 days 0% (violet)

**Color:** Violet or undyed.













**Range:** 6-0 to 1 (USP). Supplied as needed sutures and ligatures.

**Indications:** MONOCRYL\* sutures are intended for use in general soft tissue approximation and/or ligation where an absorbable material is indicated. Typical areas of use include:





- Subcuticular sutures
- Small intestine anastomoses
- Urological anastomoses




Full details in the Instruction for Use included in every package.






|  <b>1/2 Circle Taper Point Needles</b>     |       |              |             |                               |                |         |  |
|---|-------|--------------|-------------|-------------------------------|----------------|---------|--|
| Needle  | Gauge | Suture Color | Length (cm) | Other information             | Number per box | Code    |  |
|  TF<br>13 mm                               | 6-0   | ●            | 70          |                               | 36             | Y432H   |  |
|   | 5-0   | ●            | 70          |                               | 36             | Y925H   |  |
|  RB-1<br>17 mm                             | 6-0   | ●            | 70          |                               | 36             | Y212H   |  |
|   | 5-0   | ●            | 70          |                               | 36             | Y213H   |  |
|   | 4-0   | ●            | 70          |                               | 36             | Y214H   |  |
|   | 3-0   | ○            | 70          |                               | 36             | T4215AH |  |
|   | 3-0   | ●            | 70          |                               | 36             | Y215H   |  |
|  JRB-1<br><b>VISI-BLACK</b><br>17 mm       | 4-0   | ●            | 70          | VISI-BLACK*                   | 36             | Y3040H  |  |
|  JRB-1 JRB-1<br><b>VISI-BLACK</b><br>17 mm | 5-0   | ●            | 70          | VISI-BLACK*                   | 36             | Y3213H  |  |
|   | 4-0   | ●            | 70          | VISI-BLACK*                   | 36             | Y3214H  |  |
|  SH-1<br>22 mm                             | 4-0   | ●            | 70          |                               | 36             | Y218H   |  |
|   | 3-0   | ●            | 70          |                               | 36             | Y219H   |  |
|  JB-1<br><b>VISI-BLACK</b><br>22 mm      | 4-0   | ●            | 70          | VISI-BLACK*                   | 36             | Y3100H  |  |
|   | 3-0   | ●            | 70          | VISI-BLACK*                   | 36             | Y3110H  |  |
|   | 2-0   | ●            | 70          | VISI-BLACK*                   | 36             | Y3120H  |  |
|  JB-1 JB-1<br><b>VISI-BLACK</b><br>22 mm | 4-0   | ●            | 70          | VISI-BLACK*                   | 36             | Y3218H  |  |
|   | 4-0   | ●            | 90          | VISI-BLACK*                   | 36             | Y5218H  |  |
|   | 3-0   | ●            | 70          | VISI-BLACK*                   | 36             | Y3219H  |  |
|  SH<br>26 mm                             | 4-0   | ●            | 70          |                               | 36             | Y415H   |  |
|   | 3-0   | ○            | 70          |                               | 12             | NP0001  |  |
|   | 3-0   | ●            | 70          |                               | 36             | Y416H   |  |
|   | 2-0   | ●            | 70          |                               | 36             | Y417H   |  |
|  JB<br><b>VISI-BLACK</b><br>26 mm        | 4-0   | ●            | 70          | VISI-BLACK*                   | 36             | Y3150H  |  |
|   | 3-0   | ●            | 70          | VISI-BLACK*                   | 36             | Y3160H  |  |
|   | 3-0   | ●            | 8x45        | CONTROL RELEASE*, VISI-BLACK* | 12             | Y3864G  |  |
|   | 2-0   | ●            | 70          | VISI-BLACK*                   | 36             | Y3170H  |  |
|  JB JB<br><b>VISI-BLACK</b><br>26 mm     | 3-0   | ●            | 70          | VISI-BLACK*                   | 36             | Y3416H  |  |
|   | 3-0   | ●            | 90          | VISI-BLACK*                   | 36             | Y5270H  |  |
|   | 2-0   | ●            | 70          | VISI-BLACK*                   | 36             | Y3417H  |  |
|  MH<br>36 mm                             | 2-0   | ●            | 70          |                               | 36             | Y227H   |  |

**MONOCRYL\***  
(Poliglecaprone 25) Monofilament







| 1/2 Circle Taper Point Needles (Heavier Wire Diameter)  |       |              |             |                   |                |       |
|---|-------|--------------|-------------|-------------------|----------------|-------|
| Needle  | Gauge | Suture Color | Length (cm) | Other information | Number per box | Code  |
|  CT-2<br>26 mm | 2-0   | ●            | 90          |                   | 36             | Y326H |
|   | 0     | ●            | 70          |                   | 36             | Y270H |
|  CT-1<br>36 mm | 2-0   | ●            | 70          |                   | 36             | Y250H |
|   | 0     | ●            | 70          |                   | 36             | Y260H |
|   | 0     | ●            | 90          |                   | 36             | Y946H |
|   | 1     | ●            | 70          |                   | 36             | Y261H |
|  CT<br>40 mm   | 1     | ●            | 90          |                   | 36             | Y950H |
|  CTX<br>48 mm  | 2-0   | ●            | 90          |                   | 12             | Y979G |
|   | 1     | ●            | 90          |                   | 12             | Y977G |

| 5/8 Circle Taper Point Needles  |       |              |             |                   |                |       |
|---|-------|--------------|-------------|-------------------|----------------|-------|
| Needle  | Gauge | Suture Color | Length (cm) | Other information | Number per box | Code  |
|  UR-6<br>26 mm | 3-0   | ●            | 70          |                   | 12             | W3625 |
|   | 2-0   | ●            | 70          |                   | 36             | Y605H |
|  UR-6<br>26 mm | 3-0   | ●            | 70          |                   | 36             | Y602H |
|   | 2-0   | ●            | 70          |                   | 36             | A1H   |
|  UR-5<br>36 mm | 2-0   | ●            | 70          |                   | 36             | Y375H |


| 1/2 Circle TAPERCUT* Needles  |       |              |             |                   |                |       |
|---|-------|--------------|-------------|-------------------|----------------|-------|
| Needle  | Gauge | Suture Color | Length (cm) | Other information | Number per box | Code  |
|  V-7<br>26 mm  | 2-0   | ●            | 70          |                   | 36             | Y999H |
|  V-34<br>36 mm | 2-0   | ●            | 70          |                   | 12             | W3488 |
|   | 0     | ●            | 70          |                   | 12             | W3489 |
|   | 1     | ●            | 90          |                   | 36             | Y919H |

| Compound Curve TAPERCUT* Needles  |       |              |             |                   |                |        |
|---|-------|--------------|-------------|-------------------|----------------|--------|
| Needle  | Gauge | Suture Color | Length (cm) | Other information | Number per box | Code   |
|  DA-1<br>13 mm | 5-0   | ○            | 45          |                   | 12             | Y4950G |




**3/8 Circle Reverse Cutting PRIME\* Needles**



| Needle   | Gauge | Suture Color | Length (cm) | Other information | Number per box | Code     |
|--|-------|--------------|-------------|-------------------|----------------|----------|
|  P-1    | 11 mm | 6-0          | 45          | MULTIPASS*        | 36             | MPY489H  |
|  |       | 6-0          | 45          |                   | 36             | Y489H    |
|  |       | 5-0          | 45          | MULTIPASS*        | 36             | MPY490H  |
|  P-3    | 13 mm | 6-0          | 45          | MULTIPASS*        | 36             | MPY492H  |
|  |       | 6-0          | 45          |                   | 36             | Y492H    |
|  |       | 5-0          | 45          | MULTIPASS*        | 36             | MPY493H  |
|  |       | 5-0          | 45          |                   | 36             | Y491H    |
|  |       | 5-0          | 45          |                   | 36             | Y493H    |
|  |       | 4-0          | 45          | MULTIPASS*        | 36             | MPY494H  |
|  |       | 4-0          | 45          |                   | 36             | Y494H    |
|  PS-3   | 16 mm | 5-0          | 45          | MULTIPASS*        | 36             | MPY500H  |
|  |       | 5-0          | 45          |                   | 12             | W3204    |
|  |       | 4-0          | 45          | MULTIPASS*        | 36             | MPY501H  |
|  |       | 4-0          | 45          |                   | 12             | W3205    |
|  PS-2 | 19 mm | 5-0          | 45          | MULTIPASS*        | 36             | MPY495H  |
|  |       | 4-0          | 45          | MULTIPASS*        | 36             | MPY496H  |
|  |       | 4-0          | 45          |                   | 12             | W3206    |
|  |       | 4-0          | 70          | MULTIPASS*        | 36             | MPY426H  |
|  |       | 3-0          | 45          | MULTIPASS*        | 36             | MPY497H  |
|  |       | 3-0          | 45          |                   | 12             | W3207    |
|  |       | 3-0          | 70          |                   | 36             | C427     |
|  |       | 3-0          | 70          | MULTIPASS*        | 36             | MPY427H  |
|  PS-1 | 24 mm | 4-0          | 45          |                   | 36             | T4682H   |
|  |       | 4-0          | 70          | MULTIPASS*        | 36             | MPY935H  |
|  |       | 4-0          | 70          |                   | 36             | Y495H    |
|  |       | 3-0          | 70          |                   | 12             | W3208    |
|  |       | 3-0          | 70          | MULTIPASS*        | 36             | MPY683H  |
|  |       | 2-0          | 70          | MULTIPASS*        | 36             | MPY684H  |
|  |       | 2-0          | 70          |                   | 36             | Y684H    |
|  PS   | 26 mm | 4-0          | 70          | MULTIPASS*        | 36             | MPY3212H |
|  |       | 4-0          | 70          |                   | 12             | W3212    |
|  |       | 3-0          | 70          | MULTIPASS*        | 36             | MPY3213H |
|  |       | 3-0          | 70          |                   | 12             | W3213    |


**Compound Curve Reverse Cutting PRIME\* Needles**

| Needle  | Gauge | Suture Color | Length (cm) | Other information | Number per box | Code   |
|---|-------|--------------|-------------|-------------------|----------------|--------|
|  PS-4C | 16 mm | 5-0          | 45          |                   | 12             | NP0127 |

**MONOCRYL\***  
(Poliglecaprone 25) Monofilament

| ▼ 3/8 Circle Reverse Cutting Needles  |       |              |             |                   |                |       |
|---|-------|--------------|-------------|-------------------|----------------|-------|
| Needle  | Gauge | Suture Color | Length (cm) | Other information | Number per box | Code  |
|  FS-2<br>19 mm | 4-0   | ○            | 45          |                   | 36             | Y292H |
|   | 4-0   | ○            | 70          |                   | 36             | Y422H |
|   | 3-0   | ○            | 45          |                   | 36             | Y293H |
|   | 3-0   | ○            | 70          |                   | 36             | Y423H |
|  FS-1<br>24 mm | 4-0   | ○            | 70          |                   | 36             | Y291H |
|   | 3-0   | ○            | 70          |                   | 36             | Y442H |
|   | 2-0   | ○            | 70          |                   | 36             | Y443H |
|  FS<br>28 mm   | 3-0   | ○            | 70          |                   | 12             | W3326 |

| ▲ Straight Cutting Needles  |       |              |             |                   |                |       |
|---|-------|--------------|-------------|-------------------|----------------|-------|
| Needle  | Gauge | Suture Color | Length (cm) | Other information | Number per box | Code  |
|  TS<br>48 mm   | 4-0   | ○            | 70          |                   | 36             | C1041 |
|  KS<br>60 mm | 4-0   | ○            | 70          |                   | 12             | W3649 |
|   | 3-0   | ○            | 70          |                   | 12             | W3650 |

| ▼ 1/4 Circle Spatula Needles   |       |              |             |                   |                |       |
|--|-------|--------------|-------------|-------------------|----------------|-------|
| Needle   | Gauge | Suture Color | Length (cm) | Other information | Number per box | Code  |
|  S-14<br>8 mm | 6-0   | ○            | 45          |                   | 12             | W3552 |

| Ligatures |       |              |             |                   |                |       |
|-----------|-------|--------------|-------------|-------------------|----------------|-------|
| Needle    | Gauge | Suture Color | Length (cm) | Other information | Number per box | Code  |
| LIGAPAK*  | 3-0   | ●            | 250         | Reel              | 10             | W3020 |
|           | 2-0   | ●            | 250         | Reel              | 10             | W3021 |
| SUTUPAK*  | 2-0   | ●            | 150         |                   | 12             | W3113 |
|           | 1     | ●            | 150         |                   | 12             | W3115 |
|           | 0     | ●            | 150         |                   | 12             | W3114 |
|           | 2-0   | ●            | 2x70        |                   | 36             | Y625H |
|           | 3-0   | ○            | 3x45        |                   | 36             | Y634H |