


## 906544 Aperture Drape

### Product details

Size :	50cm x 60cm, Aperture 6x8cm
Descriptive feature :	Centered ellipse aperture
Sterile :	Sterile
Brand :	BARRIER
Equivalences :	904736

### Images



### Delivered items

906544-60

**Country of origin:** Thailand

**Shelf life:** 5 years

**Sterilization method:** Irradiation

**Production Responsibility:** Mölnlycke Health Care (Thailand) Ltd, Amata Nakorn (Bang Pakong), Industrial Estate, 700/461 Moo 7 Bangna-Trad Road, Km 57 Tambol Donhuaroh, Amphur Muang Chonburi 20000, Thailand

**Packing information:** First packaging layer is a peel-open sterile barrier, paper/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a corrugated board dispenser box. Third layer is a corrugated board transport box.

**Is suitable for Tray:** No

Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Piece	1	7323190258369			
Consumer pack	1	7323190274963			

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.  
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2022 Mölnlycke Health Care AB. All rights reserved.



Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Dispenser box	70	7323190262236			
Transport box	140	7323190262229	332x428x288 mm	40.9 dm <sup>3</sup>	4.4 / 3.5 kg
Pallet	5040	7323190262212	800x1200x1872 mm	1797.1 dm <sup>3</sup>	184.5 / 108.0 kg

## Technical

### Dimension

Dimension text	Dimension value
Outer dimension	50 cm x 60 cm
Aperture	6 cm x 8 cm

### Disposal instructions

Used Mölnlycke products and sterile barriers should, in the majority of cases, be classified as non-hazardous waste. They contain high amounts of energy and are well suited for incineration. Mölnlycke products do not contain any hazardous substances that can leach out if the products are land filled. Transport- and dispenser boxes are designed to fit existing recovery systems. The Mölnlycke packaging system complies with the EU Directive on packaging and packaging waste.

### Storage instructions

Mölnlycke Health Care recommends that Mölnlycke products are stored under normal storage conditions. All layers of packaging should be kept intact until access to the underlying layers is needed. Storage facilities for products only protected by the sterile barrier should be kept under conditions where low level of particulate air contamination prevail, so that it would not constitute a risk to the patient when the sterile barrier is opened and the product is used.

### Classifications

Regulation type(s)	MDR Class IS	Locally Regulated
CE Certificate Number :	MDR 722028	
Intended Purpose :	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.	The surgical patient drapes are intended to be used as protective patient covering to create a sterile field for a surgical procedure.
MDR Classification Rule :	1	
Conformity Annexes :	IX	
Measuring Function :	Yes	
Notified body medical devices/PPE :	BSi 2797	Local Authority Mexico
510(k) clearance number :		203300CT070363

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.  
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2022 Mölnlycke Health Care AB. All rights reserved.



**Applied standards :** The standards presented below is a selection of the most essential standards that are adhered to.

EN 1041, EN 556-1, EN 13795, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 11607-1, EN ISO 11607-2, EN ISO 15223-1, ISO 15223-2, ISO 14001

### **Removable label**

Yes

### **GMDN Code (Global Medical Device Nomenclature)**

47783 Patient surgical drape, single-use

### **UNSPSC**

42131701 Surgical drapes

### **Commodity Code**

3005100000 Wadding, Gauze, Dressings, drapes singlepacked - adhesive articles, Sets mainly consisting thereof

### **CE Responsibility**

Mölnlycke Health Care AB, Gamlestadsvägen 3c, Box 13080, SE-402 52 Goteborg Sweden

### **Basic UDIDI**

733243000000000160JX

**Find out more at [www.molnlycke.com](http://www.molnlycke.com)**

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.  
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2022 Mölnlycke Health Care AB. All rights reserved.

