


You want to extend or connect?

More mobility for your patients

Prolonged bed rest and immobilization are the standards of the past. Nowadays, we know that it is usually better for patients to get back on their feet as soon as possible.

Early mobilization effectively promotes healing. Today's mobile patients need flexible solutions to keep their infusions running.



Spiral-Tubings	Length (cm)	Priming volume (ml)	Inner ø (mm)	Material	Sales unit/ Pcs.	Code No. (REF)
	150	1.2	1.0	PUR	50	4092945
	150	4.7	2.0	PUR	50	4090365
	300	2.3	1.0	PUR	50	4092937
	300	9.4	2.0	PUR	50	4090373

Perfusor® lines	Length (cm)	Priming volume (ml)	Inner ø (mm)	Material	Sales unit/ Pcs.	Code No. (REF)
Extension and connection lines	30	0.2	0.9	PVC (DEHP-free)	100	9500049
	30	0.3	1.2	PVC (DEHP-free)	100	4256000 ¹
	30	0.3	1.2	PVC (DEHP-free)	100	4256034
Perfusor® lines	300	5.3	1.5	PVC (DEHP-free)	100	8255253
	250	4.4	1.5	PVC (DEHP-free)	100	8255490
	200	3.5	1.5	PVC (DEHP-free)	100	8722862
	150	2.6	1.5	PVC (DEHP-free)	100	8722960
	50	0.9	1.5	PVC (DEHP-free)	100	8255172
	200	3.3	1.5	PVC (DEHP-free)	100	8723001 ²
	150	2.6	1.5	PVC (DEHP-free)	100	8726019 ³
	150	2.7	1.5	PVC (DEHP-free)	100	8722820
	150	2.7	1.5	PVC (DEHP-free)	100	8722919 ⁴
	75	0.5	1.0	PVC (DEHP-free)	100	8722870
	150	0.9	1.0	PVC (DEHP-free)	100	8255504
	50	0.4	1.0	PE	100	8255059
	100	0.8	1.0	PE	100	8255067
	300	2.4	1.0	PE	100	8250146
	250	2.0	1.0	PE	100	8272565
200	1.6	1.0	PE	100	8723060	
150	1.2	1.0	PE	100	8722935	
150	1.2	1.0	PE	100	8723010 ⁴	
150	1.2	1.0	PE	100	8723017 ⁵	