


832000 Fluid Collection Pouch W/Incise

Product details

Product group name :	Fluid collection pouch w. incise
Size :	80cm x 100cm, Incise 38x32cm
Descriptive feature :	Incise film
Colour :	Transparent
Sterile :	Sterile
Brand :	BARRIER®
Equivalences :	983200

Images



Delivered items

832000-78

Country of origin: Thailand

Shelf life: 5 years

Sterilization method: Irradiation

Packing information: First packaging layer is a peel-open sterile barrier, plastic/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a plastic liner/bag. Third layer is a corrugated board dispenser box. Fourth layer is a corrugated board transport box.

Is suitable for Tray: No

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.

The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2024 Mölnlycke Health Care AB. All rights reserved.

Packing level	Quantity	GTIN Code	UDI-DI	Width x Length x Height	Volume	Weight gross / net
Piece	1	7332430072306				352.0 g / -
Consumer pack	1	07333350111618	07333350111618			
Dispenser box	10	07333350767945	07333350767945	246x366x569 mm		
Transport box	20	07333350755706	07333350755706	373x588x505 mm	110.8 dm3	7.0 / 5 kg
Pallet	240	7313661374926		800x1200x1659 mm	1592.6 dm3	109.5 kg / -

832000-01

Country of origin: Thailand**Shelf life:** 5 years**Sterilization method:** Irradiation

Packing information: First packaging layer is a peel-open sterile barrier, plastic/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a plastic liner/bag. Third layer is a corrugated board dispenser box. Fourth layer is a corrugated board transport box.

Is suitable for Tray: No

Packing level	Quantity	GTIN Code	UDI-DI	Width x Length x Height	Volume	Weight gross / net
Piece	1	7332430072306			5034.5 cm3	308.0 g / -
Consumer pack	1	07333350579593	07333350579593			
Dispenser box	11	07333350852245	07333350852245			
Transport box	22	07333350383398	07333350383398	373x588x505 mm	110.8 dm3	6.8 / 4.2 kg
Pallet	264	7313661157789		800x1200x1659 mm	1592.6 dm3	96.3 kg / -

832000-00

Country of origin: Thailand**Shelf life:** 5 yearsFind out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.

The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2024 Mölnlycke Health Care AB. All rights reserved.



Sterilization method: Irradiation

Production Responsibility: Mölnlycke Health Care (Thailand) Ltd, 160 Moo 17 Bangplee Industrial Estate, Bangna-Trad Road, Km 23 Tambol Bangsaothong, Amphur Bangsaothong, Samutprakarn 10570, Thailand

Packing information: First packaging layer is a peel-open sterile barrier, plastic/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a plastic liner/bag. Third layer is a corrugated board dispenser box. Fourth layer is a corrugated board transport box.

Is suitable for Tray: No

Packing level	Quantity	GTIN Code	UDI-DI	Width x Length x Height	Volume	Weight gross / net
Consumer pack	1	7310791268903	7310791268903			
Dispenser box	11	7310791268897	7310791268897			
Transport box	22	7310791268880	7310791268880	373x588x505 mm	110.8 dm3	6.8 / 4.2 kg
Pallet	264	7332430444912		800x1200x1659 mm		

Technical

Dimension

Dimension text	Dimension value
Outer dimension	100 cm x 80 cm
Incise	38 cm x 32 cm

Disposal instructions

Used Mölnlycke products and sterile barriers should, in the majority of cases, be classified as non-hazardous waste. They contain high amounts of energy and are well suited for incineration. Mölnlycke products do not contain any hazardous substances that can leach out if the products are land filled. Transport- and dispenser boxes are designed to fit existing recovery systems. The Mölnlycke packaging system complies with the EU Directive on packaging and packaging waste.

Storage instructions

Mölnlycke Health Care recommends that Mölnlycke products are stored under normal storage conditions. All layers of packaging should be kept intact until access to the underlying layers is needed. Storage facilities for products only protected by the sterile barrier should be kept under conditions where low level of particulate air contamination prevail, so that it would not constitute a risk to the patient when the sterile barrier is opened and the product is used.

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.

The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2024 Mölnlycke Health Care AB. All rights reserved.



Classifications

Regulation type(s)	MDR Class IS	Locally Regulated	Locally Regulated	Locally Regulated
CE Certificate Number :	UKCA 760670			
Intended Purpose :	The pouches are attached to other draped areas and are used for fluid management during surgical procedures.	The pouches are attached to other draped areas and are used for fluid management during surgical procedures.	The pouches are attached to other draped areas and are used for fluid management during surgical procedures.	The pouches are attached to other draped areas and are used for fluid management during surgical procedures.
MDR Classification Rule :	1			
Conformity Annexes :	V			
Measuring Function :	No			
Notified body medical devices/PPE :	BSi 2797	Local Authority Americas	Local Authority Thailand	Local Authority APAC

Classifications

Regulation type(s)	Locally Regulated	Locally Regulated
CE Certificate Number :		
Intended Purpose :	The pouches are attached to other draped areas and are used for fluid management during surgical procedures.	The pouches are attached to other draped areas and are used for fluid management during surgical procedures.
MDR Classification Rule :		
Conformity Annexes :		
Measuring Function :		
Notified body medical devices/PPE :	Local Authority MEA	Local Authority United Arab Emirates

Applied standards : The standards presented below is a selection of the most essential standards that are adhered to.

EN 1041, EN 556-1, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 11607-1, EN ISO 11607-2, EN ISO 15223-1, ISO 15223-2, ISO 14001

Removable label

Yes

GMDN Code (Global Medical Device Nomenclature)

56731 Surgical tissue/fluid collection bag, sterile

Find out more at www.molnlycke.com



UNSPSC

42293507 Surgical suction reservoirs

Commodity Code

3926909790 Plastics and articles thereof

CE Responsibility / Legal Manufacturer

Mölnlycke Health Care AB, Gamlestadvägen 3c, Box 13080, SE-402 52 Goteborg Sweden

Basic UDIDI

733243000000000034JR

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.

The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2024 Mölnlycke Health Care AB. All rights reserved.

