








Clear, concise labeling on the package and pads guides experienced and infrequent caregivers alike through pad selection, placement, and connection.



Technical specifications - ALS

	989803166021	M3713A	M3716A	M3718A	M3501A	M3717A	M3719A	M3504A	
Description/name	Adult/child preconnect	Adult/child Plus	Adult/child radiolucent	Adult/child radiotransparent	Adult/child with barrel connector	Infant Plus	Infant radiotransparent	Infant with barrel connector	
Pads per case	10	10	10	10	10	5	5	5	
Indications	Monitoring, external pacing, defibrillation, and cardioversion								
Patient						Infant	Infant	Infant	
Infant < 10 kg (22 lbs)									
Child 10 kg – 25 kg	Child	Child	Child	Child	Child				
Adult > 25 kg (55 lbs)	Adult	Adult	Adult	Adult	Adult				
Conductor gel surface area (cm²)	75	75	75	97	75	32	50	32	
Approximate lead wire length (cm)	120	120	120	120	60	60	60	60	
Number of hours on skin	24	24	24	24	24	24	24	24	
Number of defibrillation shocks	50	50	50	50	50	50	50	50	
Duration of pacing (hours) at 100mA, 80 bpm	8	8	8	3	8	8	1.5	8	
Temperature range for storage	15° – 35°C	15° – 35°C	15° – 35°C	15° – 35°C	15° – 35°C	15° – 35°C	15° – 35°C	15° – 35°C	
Shelf life (in months) from date of manufacture (may be less upon receipt)	24	24	24	24	24	24	18	24	
Philips adapter compatibility	05-10000 05-10100 05-10200	05-10000 05-10100 05-10200	05-10000 05-10100 05-10200	05-10000 05-10100 05-10200		N/A	05-10000 05-10100	05-10000 05-10100	N/A

Defibrillator compatibility

Intrepid* 867172; M3508A plug connector cable	●	●	●	●		●	●	
Intrepid* 867172; M3507A barrel connector cable					●			●
DFM100* 866199; M3508A plug connector cable		●	●	●		●	●	
DFM100* 866199; M3507A barrel connector cable					●			●
XL+** 861290; M3508A plug connector cable	●	●	●	●		●	●	
XL+** 861290; M3507A barrel connector cable					●			●
MRx** M3535A, M3536A; M3508A plug connector cable	○ ²	○ ²	○ ²	○ ²		○ ¹	○ ¹	
MRx** M3535A, M3536A; M3507A barrel connector cable					○ ²			○ ¹
FR3 861388, 861389	○ ^{3,4}	○ ^{3,4}	○ ^{3,4}					
FRx 861304, 861305	○ ⁵	○ ⁵	○ ⁵					
FR2/FR2+ M3840A, M3841A, M3860A, M3861A, 861458, 861459	○ ^{3,4}	○ ^{3,4}	○ ^{3,4}	○ ^{3,4}				


* Not available for sale in the USA.

** MRx and XL+ have been discontinued and no longer for sale.


Not all products are available in all geographies. Please contact your Philips representative.

● Fully compatible
○ Partially compatible


Technical specifications - AEDs



989803158211 (1 set)
989803158221 (5 sets)



989803139261



989803149981 (1 set)
989803149991 (5 sets)



M3870A

Description/name	Defibrillator pads	SMART Pads II	SMART Pads III	FR2 Infant/child
Pads per case	1 or 5	1	1 or 5	1
Indications	Monitoring, external pacing, defibrillation, and cardioversion			Defibrillation
Patient				
Infant < 10 kg (22 lbs)	Infant	Infant	Infant	Infant
Child 10 kg – 25 kg	Child	Child	Child	Child
Adult > 25 kg (55 lbs)	Adult	Adult	Adult	
Conductor gel surface area (cm²)	85	80	80	85
Approximate lead wire length (cm)	120	120	120	120
Number of hours on skin	24	24	24	24
Number of defibrillation shocks	50	50	50	50
Duration of pacing (hours) at 100mA, 80 bpm	8	8	8	N/A
Temperature range for storage	0° – 50°C	0° – 50°C	0° – 50°C	0° – 32°C
Shelf life (in months) from date of manufacture (may be less upon receipt)	30	30	30	30
Philips adapter compatibility	05-10000 05-10100 05-10200	05-10000 05-10100 05-10200	05-10000 05-10100 05-10200	N/A

Defibrillator compatibility

Intrepid* 867172; M3508A plug connector cable	●	●	●	
Intrepid* 867172; M3507A barrel connector cable				
DFM100* 866199; M3508A plug connector cable	●	●	●	
DFM100* 866199; M3507A barrel connector cable				
XL** 861290; M3508A plug connector cable	●	●	●	
XL** 861290; M3507A barrel connector cable				
MRx** M3535A, M3536A; M3508A plug connector cable	○ ²	○ ²	○ ²	
MRx** M3535A, M3536A; M3507A barrel connector cable				
FR3 861388, 861389	○ ⁴	●	●	
FRx 861304, 861305		●		
FR2/FR2+ M3840A, M3841A, M3860A, M3861A, 861458, 861459	○ ^{3,4}	○ ^{3,4}	○ ^{3,4}	○ ⁴

* Not available for sale in the USA.

** MRx and XL+ have been discontinued and no longer for sale.

Not all products are available in all geographies. Please contact your Philips representative.

1 Not in AED mode, as the energy automatically delivered is 150 Joules

2 Not for use with infants/children (10–25 kg / 22–55 lbs) in AED Mode

3 Infant/Child Key is required for use with infants/children (10–25 kg / 22–55 lbs) in Manual and AED Mode

4 Not preconnectable and cannot be left attached to device after use, when in Standby Mode

5 Infant/Child Key is required for use with infants/children (10–25 kg / 22–55 lbs) in AED Mode

© 2021 Koninklijke Philips N.V. All rights are reserved.

● Fully compatible
○ Partially compatible

4522 991 60851 * JAN 2021