PHILIPS

Hospital respiratory care

Confidence in your patient's comfort

Philips Respironics patient interfaces and circuits for noninvasive ventilation

Delivering exceptional NIV

Philips Respironics masks are designed to support enhanced NIV. From comfort for your patients to ease of use for you, our interfaces are specially designed to meet the rigors of acute care environments.

Table of contents

| 9 | Total face masks | 4 |
|---|--------------------------|----|
| | Oro-nasal masks | 8 |
| 4 | Nasal masks | 16 |
| | Packs and accessories | 20 |
| Ð | Circuits and connections | 28 |
| | Trilogy accessories | 30 |

Total face masks

Total face masks allow quick fitting and eliminate nasal bridge challenges by sealing around the perimeter of the face, where patients have less pressure sensitivity and smoother facial contours.

Mask selection criteria

- Rapid application
- Mouth breather
- Facial abnormalities
- Lack of teeth



Philips Respironics Total



Key benefits

- Clear, unobstructed view
- Perimeter facial contact increases comfort and simplifies fitting
- Digital Auto-Trak compatibility enhances NIV performance

For headgear, see page 25

Total

Philips Respironics Total multi-use

| Description | Quantity | Product number | Part number |
|----------------------------|----------|----------------|--------------|
| Total face mask | 1 | 302433 | 989805609231 |
| Total face mask head strap | 1 | 1003052 | 453561517111 |
| Total face mask head strap | 9 | 1010609 | 989805610031 |

Philips Respironics PerforMax single-use



Key benefits

- Available for patients 7 years or older and > 20 kg (40 lb)
- Perimeter facial contact increases comfort and simplifies fitting
- Clear, unobstructed view

Compatible with a leak 2 entrainment elbow (EE), standard elbow (SE), and bronchoscopy elbow (BE) For accessories, see page 24 For headgear, see page 25

PerforMax EE leak 2

Philips Respironics PerforMax single-use

| Description | Quantity | Product number | Part number |
|--|----------|----------------|--------------|
| PerforMax with EE leak 2 elbow, small*† | 1 | 1061619 | 989805627831 |
| PerforMax with EE leak 2 elbow, large † | 1 | 1061620 | 989805627841 |
| PerforMax with EE leak 2 elbow, extra large [†] | 1 | 1080426 | 989805634831 |
| PerforMax with EE leak 2 elbow, small*† | 10 | 1052535 | 989805611281 |
| PerforMax with EE leak 2 elbow, large † | 10 | 1052536 | 989805611291 |
| PerforMax with EE leak 2 elbow, extra large [†] | 5 | 1080424 | 989805634811 |
| PerforMax with SE standard elbow, small*† | 1 | 1061621 | 989805627851 |
| PerforMax with SE standard elbow, large ⁺ | 1 | 1061622 | 989805627861 |
| PerforMax with SE standard elbow, extra large ⁺ | 1 | 1080517 | 989805634841 |
| PerforMax with SE standard elbow, small*† | 10 | 1052547 | 989805611301 |
| PerforMax with SE standard elbow, large † | 10 | 1052548 | 989805611311 |
| PerforMax with SE standard elbow, extra large [†] | 5 | 1080425 | 989805634821 |

* Intended for use by patients 7 years or older and > 20 kg (40 lb)

⁺ Not available in all markets, please contact your Philips sales consultant.

Philips Respironics PerforMax multi-use



Key benefits

- Available for patients 7 years or older and > 20 kg (40 lb)
- Perimeter facial contact supports comfort and simplifies fitting
- Clear, unobstructed view

Compatible with a leak 2 entrainment elbow (EE), standard elbow (SE), and bronchoscopy elbow (BE) For accessories, see page 24 For headgear, see page 25

PerforMax SE

Philips Respironics PerforMax multi-use

| Description | Quantity | Product number | Part number |
|--|-------------------|----------------|--------------|
| PerforMax with EE leak 2 and SE elbow, small* | 1 mask/4 headgear | 1052553 | 989805611361 |
| PerforMax with EE leak 2 and SE elbow, large | 1 mask/4 headgear | 1052554 | 989805611371 |
| PerforMax with EE leak 2 and SE elbow, extra large | 1 mask/4 headgear | 1080420 | 989805634771 |
| PerforMax with EE leak 2, small* | 1 mask/4 headgear | 1052555 | 989805611381 |
| PerforMax with EE leak 2, large | 1 mask/4 headgear | 1052556 | 989805611391 |
| PerforMax with EE leak 2, extra large | 1 mask/4 headgear | 1080421 | 989805634781 |
| PerforMax with SE standard elbow, small* | 1 mask/4 headgear | 1052568 | 989805611401 |
| PerforMax with SE standard elbow, large | 1 mask/4 headgear | 1052569 | 989805611411 |
| PerforMax with SE standard elbow, extra large | 1 mask/4 headgear | 1080422 | 989805634791 |
| PerforMax with EE leak 2 and SE elbow, small* | 1 mask/1 headgear | 1044385 | 989805611091 |
| PerforMax with EE leak 2 and SE elbow, large | 1 mask/1 headgear | 1044386 | 989805611101 |
| PerforMax with EE leak 2 and SE elbow, extra large | 1 mask/1 headgear | 1080423 | 989805634801 |

Oro-nasal masks

Oro-nasal masks are the standard of care for hospital NIV, and the Philips Respironics family includes a comprehensive selection of systems with features that address patient comfort, mask fitting, and ventilator compatibility.



Mask selection criteria

- Mouth breather
- Long-term NIV
- Two cushion styles for mask rotation

Patient comfort made easy

CapStrap headgear provides excellent stability and simplifies initial setup by keeping the mask in place while you adjust the patient's straps. When you need to interrupt treatment for oral care and medications, the CapStrap headgear allows quick removal and simple reapplication. CapStrap headgear is available on several Philips Respironics masks:

- AF541
- AF531
- AF811
- PerformaTrak



AF541 with four-point headgear, EE leak 1

Key benefits

- The under-the-nose cushion can be used to off load mask pressure at bridge of nose
- Push-button forehead pad adjustment, easy-to-use CapStrap headgear, and snap-in-place elbows and talon clips contribute to efficient workflow
- Discharge kit configures mask to work with home NIV systems
- CleanClip system supports your infection control efforts
- Discharge kit for patient transition to lower acuity setting

For accessories, see page 24 For headgear, see page 25

Philips Respironics AF541

| Description | Quantity | Product number | Part number |
|--|----------|----------------|--------------|
| AF541 with EE leak 1 elbow, CapStrap headgear, small* | 10 | 1120925 | 989805653931 |
| AF541 with EE leak 1 elbow, CapStrap headgear, medium* | 10 | 1120926 | 989805653941 |
| AF541 with EE leak 1 elbow, CapStrap headgear, large | 10 | 1120927 | 989805653951 |
| AF541 with EE leak 1 elbow, CapStrap headgear, extra large | 10 | 1120928 | 989805653961 |
| AF541 with EE leak 2 elbow, CapStrap headgear, small* | 10 | 1120935 | 989805654011 |
| AF541 with EE leak 2 elbow, CapStrap headgear, medium* | 10 | 1120936 | 989805654021 |
| AF541 with EE leak 2 elbow, CapStrap headgear, large | 10 | 1120937 | 989805654031 |
| AF541 with EE leak 2 elbow, CapStrap headgear, extra large | 10 | 1120938 | 989805654041 |
| AF541 with SE standard elbow, CapStrap headgear, small* | 10 | 1120945 | 989805654091 |
| AF541 with SE standard elbow, CapStrap headgear, medium* | 10 | 1120946 | 989805654101 |
| AF541 with SE standard elbow, CapStrap headgear, large | 10 | 1120947 | 989805654111 |
| AF541 with SE standard elbow, CapStrap headgear, extra large | 10 | 1120948 | 989805654121 |
| AF541 with EE leak 1 elbow, four-point headgear, small* | 10 | 1120920 | 989805653891 |
| AF541 with EE leak 1 elbow, four-point headgear, medium* | 10 | 1120921 | 989805653901 |
| AF541 with EE leak 1 elbow, four-point headgear, large | 10 | 1120922 | 989805653911 |
| AF541 with EE leak 1 elbow, four-point headgear, extra large | 10 | 1120923 | 989805653921 |
| AF541 with EE leak 2 elbow, four-point headgear, small* | 10 | 1120930 | 989805653971 |
| AF541 with EE leak 2 elbow, four-point headgear, medium* | 10 | 1120931 | 989805653981 |
| AF541 with EE leak 2 elbow, four-point headgear, large | 10 | 1120932 | 989805653991 |
| AF541 with EE leak 2 elbow, four-point headgear, extra large | 10 | 1120933 | 989805654001 |
| AF541 with SE standard elbow, four-point headgear, small* | 10 | 1120940 | 989805654051 |
| AF541 with SE standard elbow, four-point headgear, medium* | 10 | 1120941 | 989805654061 |
| AF541 with SE standard elbow, four-point headgear, large | 10 | 1120942 | 989805654071 |
| AF541 with SE standard elbow, four-point headgear, extra large | 10 | 1120943 | 989805654081 |
| AF541 size A under-the-nose cushion | 10 | 1120980 | 989805654271 |
| AF541 size B under-the-nose cushion | 10 | 1120981 | 989805654281 |
| AF541 size C under-the nose cushion | 10 | 1120982 | 989805654291 |
| Bronchoscopy elbow | 5 | 1120984 | 989805654301 |
| AF541 discharge kit, CapStrap headgear, small/medium | 5 | 1120986 | 989805654311 |
| AF541 discharge kit, CapStrap headgear, large/extra large | 5 | 1120988 | 989805654331 |
| AF541 discharge kit, four-point headgear, small/medium | 5 | 1120987 | 989805654321 |
| AF541 discharge kit, four-point headgear, large/extra large | 5 | 1120989 | 989805654341 |





AF531 with CapStrap

AF531 EE leak 1 elbow

Philips Respironics AF531 with CapStrap single-use

Key benefits

- CapStrap headgear provides an excellent fit and simple reapplication
- Supports your infection control efforts with the CleanClip system
- You can adapt one mask to multiple ventilators with interchangeable elbows

Compatible with a leak 1 or leak 2 entrainment elbow (EE), EE NIVO elbow (EE/NE), standard elbow (SE), SE NIVO elbow (SE/NE), and bronchoscopy elbow (BE) For accessories, see page 24 For headgear, see page 25

| Description | Quantity | Product number | Part number |
|--|----------|----------------|--------------|
| AF531 with EE leak 1 elbow, CapStrap headgear, small* | 1 | 1061727 | 989805628611 |
| AF531 with EE leak 1 elbow, CapStrap headgear, medium | 1 | 1061728 | 989805628621 |
| AF531 with EE leak 1 elbow, CapStrap headgear, large | 1 | 1061729 | 989805628631 |
| AF531 with EE leak 1 elbow, CapStrap headgear, small* | 10 | 1072627 | 989805628731 |
| AF531 with EE leak 1 elbow, CapStrap headgear, medium | 10 | 1072628 | 989805628741 |
| AF531 with EE leak 1 elbow, CapStrap headgear, large | 10 | 1072629 | 989805628751 |
| AF531 with SE standard elbow, CapStrap headgear, small* | 1 | 1061730 | 989805628641 |
| AF531 with SE standard elbow, CapStrap headgear, medium | 1 | 1061731 | 989805628651 |
| AF531 with SE standard elbow, CapStrap headgear, large | 1 | 1061732 | 989805628661 |
| AF531 with SE standard elbow, CapStrap headgear, small* | 10 | 1072630 | 989805628761 |
| AF531 with SE standard elbow, CapStrap headgear, medium | 10 | 1072631 | 989805628771 |
| AF531 with SE standard elbow, CapStrap headgear, large | 10 | 1072632 | 989805628781 |
| AF531 with EE leak 1/NE nebulizer elbow, CapStrap headgear, small* | 10 | 1072636 | 989805635001 |
| AF531 with EE leak 1/NE nebulizer elbow, CapStrap headgear, medium | 10 | 1072637 | 989805635011 |
| AF531 with EE leak 1/NE nebulizer elbow, CapStrap headgear, large | 10 | 1072638 | 989805635021 |

* Intended for use by patients 7 years or older and > 20 kg (40 lb)







Key benefits

- Supports your infection control efforts with the CleanClip system
- You can adapt one mask to multiple ventilators with interchangeable elbows

For accessories, see page 24 For headgear, see page 25

AF531 with four-point headgear, EE leak 1

AF531 Bronchoscopy elbow*

AF531 SE standard elbow

Philips Respironics AF531 with four-point headgear single-use

| Quantity | Product number | Part number |
|----------|---|---|
| 1 | 1061718 | 989805628521 |
| 1 | 1061719 | 989805628531 |
| 1 | 1061720 | 989805628541 |
| 10 | 1072621 | 989805628671 |
| 10 | 1072622 | 989805628681 |
| 10 | 1072623 | 989805628691 |
| 1 | 1061721 | 989805628551 |
| 1 | 1061722 | 989805628561 |
| 1 | 1061723 | 989805628571 |
| 10 | 1072624 | 989805628701 |
| 10 | 1072625 | 989805628711 |
| 10 | 1072626 | 989805628721 |
| 10 | 1061633 | 989805634971 |
| 10 | 1061634 | 989805634981 |
| 10 | 1061635 | 989805634991 |
| | 1 1 1 10 10 10 10 1 1 1 1 10 10 | 1 1061718 1 1061719 1 1061720 10 1072621 10 1072622 10 1072623 10 1072623 1 1061721 1 1061722 1 1061723 10 1072624 10 1072625 10 1072626 10 1072626 10 1072626 10 1061633 10 1061633 10 1061634 |

* Bronchoscopy elbow sold separately, see page 24

** Intended for use by patients 7 years or older and > 20 kg (40 lb)

Philips Respironics AF531 with leak 2 elbow



Key benefits

- Leak built into the mask
- Can be used without an additional exhalation port
- \cdot Available with four point headgear
- and CapStrap headgear
- Convenient ten-packs

For accessories, see page 24 For headgear, see page 25

AF531 EE leak 2

Philips Respironics AF531 with leak 2 elbow single-use

| Description | Quantity | Product number | Part number |
|---|----------|----------------|--------------|
| AF531 with EE leak 2 elbow, four-point headgear, small* | 1 | 1061724 | 989805628581 |
| AF531 with EE leak 2 elbow, four-point headgear, medium | 1 | 1061725 | 989805628591 |
| AF531 with EE leak 2 elbow, four-point headgear, large | 1 | 1061726 | 989805628601 |
| AF531 with EE leak 2 elbow, four-point headgear, small* | 10 | 1100393 | 989805645651 |
| AF531 with EE leak 2 elbow, four-point headgear, medium | 10 | 1100394 | 989805645661 |
| AF531 with EE leak 2 elbow, four-point headgear, large | 10 | 1100395 | 989805645671 |
| AF531 with EE leak 2 elbow, CapStrap headgear, small* | 10 | 1100390 | 989805645621 |
| AF531 with EE leak 2 elbow, CapStrap headgear, medium | 10 | 1100391 | 989805645631 |
| AF531 with EE leak 2 elbow, CapStrap headgear, large | 10 | 1100392 | 989805645641 |
| | | | |

* Intended for use by patients 7 years or older and > 20 kg (40 lb)



Key benefits

- Cushion uses a gel pad with a silicone flap for increased comfort
- CapStrap headgear system simplifies mask application and routine patient care
- Available in three sizes to fit a wide patient population

Entrainment elbow (EE) Also available with standard elbow (SE) For headgear, see page 25

AF811 EE

Philips Respironics AF811 single-use

| Description | Quantity | Product number | Part number |
|--|----------|----------------|--------------|
| AF811 oro-nasal mask EE, CapStrap headgear, small | 1 | 1063047 | 989805626321 |
| AF811 oro-nasal mask EE, CapStrap headgear, medium | 1 | 1063016 | 989805626311 |
| AF811 oro-nasal mask EE, CapStrap headgear, large | 1 | 1063015 | 989805626301 |
| AF811 oro-nasal mask EE, CapStrap headgear, small | 5 | 1063053 | 989805626381 |
| AF811 oro-nasal mask EE, CapStrap headgear, medium | 5 | 1063052 | 989805626371 |
| AF811 oro-nasal mask EE, CapStrap headgear, large | 5 | 1063051 | 989805626361 |
| AF811 oro-nasal mask SE, CapStrap headgear, small | 1 | 1063050 | 989805626351 |
| AF811 oro-nasal mask SE, CapStrap headgear, medium | 1 | 1063049 | 989805626341 |
| AF811 oro-nasal mask SE, CapStrap headgear, large | 1 | 1063048 | 989805626331 |
| AF811 oro-nasal mask SE, CapStrap headgear, small | 5 | 1063056 | 989805626411 |
| AF811 oro-nasal mask SE, CapStrap headgear, medium | 5 | 1063055 | 989805626401 |
| AF811 oro-nasal mask SE, CapStrap headgear, large | 5 | 1063054 | 989805626391 |
| Aron oro-nasar mask SE, Capstrap neaugear, targe | 5 | 1003034 | 989803020391 |

Philips Respironics PerformaTrak



PerformaTrak

Key benefits

- Dual-density foam forehead pads provide patient comfort and mask stability
- With only three sizes, clinicians can quickly fit a wide range of patients
- Swivel clips simplify reapplication by reducing the need to refit after removal

Entrainment elbow (EE) Also available with standard elbow (SE) For headgear, see page 25

Philips Respironics PerformaTrak single-use

| Description | Quantity | Product number | Part number |
|--|----------|----------------|--------------|
| PerformaTrak oro-nasal mask EE leak 1, four-point headgear, small | 1 | 1012573 | 989805610111 |
| PerformaTrak oro-nasal mask EE leak 1, four-point headgear, medium | 1 | 1012574 | 989805610121 |
| PerformaTrak oro-nasal mask EE leak 1, four-point headgear, large | 1 | 1012635 | 989805610151 |
| PerformaTrak oro-nasal mask EE leak 1, four-point headgear, small | 10 | 1018589 | 989805610481 |
| PerformaTrak oro-nasal mask EE leak 1, four-point headgear, medium | 10 | 1018590 | 989805610491 |
| PerformaTrak oro-nasal mask EE leak 1, four-point headgear, large | 10 | 1018591 | 989805610501 |
| PerformaTrak oro-nasal mask SE, four-point headgear, small | 1 | 1012639 | 989805610191 |
| PerformaTrak oro-nasal mask SE, four-point headgear, medium | 1 | 1012640 | 989805610201 |
| PerformaTrak oro-nasal mask SE, four-point headgear, large | 1 | 1012641 | 989805610211 |
| PerformaTrak oro-nasal mask SE, four-point headgear, small | 10 | 1018605 | 989805610581 |
| PerformaTrak oro-nasal mask SE, four-point headgear, medium | 10 | 1018606 | 989805610591 |
| PerformaTrak oro-nasal mask SE, four-point headgear, large | 10 | 1018595 | 989805610541 |
| PerformaTrak oro-nasal mask with CapStrap EE, small | 10 | 1019492 | 989805610691 |
| PerformaTrak oro-nasal mask with CapStrap EE, medium | 10 | 1019498 | 989805610701 |
| PerformaTrak oro-nasal mask with CapStrap EE, large | 10 | 1019507 | 989805610711 |
| PerformaTrak oro-nasal mask with CapStrap SE, small | 10 | 1019540 | 989805610811 |
| PerformaTrak oro-nasal mask with CapStrap SE, medium | 10 | 1019541 | 989805610821 |
| PerformaTrak oro-nasal mask with CapStrap SE, large | 10 | 1019542 | 989805610831 |

Philips Respironics PerformaTrak multi-use

| Description | Quantity | Product number | Part number |
|--|----------|----------------|--------------|
| PerformaTrak oro-nasal mask EE leak 1, four-point headgear, small | 1 | 1022980 | 989805610871 |
| PerformaTrak oro-nasal mask EE leak 1, four-point headgear, medium | 1 | 1023022 | 989805610881 |
| PerformaTrak oro-nasal mask EE leak 1, four-point headgear, large | 1 | 1023023 | 989805610891 |
| PerformaTrak SE oro-nasal mask, four-point-headgear, small | 1 | 1023025 | 989805610901 |
| PerformaTrak SE oro-nasal mask, four-point-headgear, medium | 1 | 1023026 | 989805610911 |
| PerformaTrak SE oro-nasal mask, four-point-headgear, large | 1 | 1023027 | 989805610921 |
| PerformaTrak four-point headgear assembly with forehead pad | 5 | 1022453 | 989805610861 |





AF421 with four-point headgear, SE

AF421 EE leak 1



AF421 EE leak 2

Key benefits

- Supports your infection control efforts with CleanClip
- \cdot Oral access feature for oral care
- $\boldsymbol{\cdot}$ Leak symbols for ease of ventilator setup
- Available with EE leak 1, SE standard elbow, and EE leak 2 elbow for wide ventilator application
- $\boldsymbol{\cdot}$ Available in four sizes to fit a wide patient population

Philips Respironics AF421

| Description | Quantity | Product number | Part number |
|---|----------|----------------|--------------|
| AF421 with EE leak 1 elbow, small* | 10 pack | 1106170 | 989805648331 |
| AF421 with EE leak 1 elbow, medium | 10 pack | 1106171 | 989805648341 |
| AF421 with EE leak 1 elbow, large | 10 pack | 1106172 | 989805648351 |
| AF421 with EE leak 1 elbow, extra large | 10 pack | 1106173 | 989805648361 |
| AF421 with SE standard elbow, small* | 10 pack | 1106174 | 989805648371 |
| AF421 with SE standard elbow, medium | 10 pack | 1106175 | 989805648381 |
| AF421 with SE standard elbow, large | 10 pack | 1106176 | 989805648391 |
| AF421 with SE standard elbow, extra large | 10 pack | 1106177 | 989805648401 |
| AF421 with EE leak 2 elbow, small* | 10 pack | 1106178 | 989805648411 |
| AF421 with EE leak 2 elbow, medium | 10 pack | 1106179 | 989805648421 |
| AF421 with EE leak 2 elbow, large | 10 pack | 1106180 | 989805648431 |
| AF421 with EE leak 2 elbow, extra large | 10 pack | 1106181 | 989805648441 |
| Accessory, AF421 forehead pad (3/4 in, 19.1 mm) | 10 pack | 1118318 | 989805652341 |
| RP, S/M AF421 headgear with clips | 5 pack | 1109915 | 989805648851 |
| RP, L/XL AF421 headgear with clips | 5 pack | 1109916 | 989805648861 |
| RP, AF421 sizing gauge | 25 pack | 1109917 | 989805648871 |
| | | | |

* Intended for use by patients 7 years or older and > 20 kg (40 lb)

Nasal masks

When chronically ill patients enter the hospital in respiratory distress, they typically want an NIV interface as comfortable as what they use at home. They want an unobtrusive system to maintain quality of life during their hospital stay.

Mask selection criteria

- Oral access
- Long-term NIV



Philips Respironics PerformaTrak



Key benefits

- CapStrap headgear provides an excellent fit and simple reapplication
- With only three sizes, clinicians can quickly fit a wide range of patients
- Hook-n-loop grab tabs allow quick, easy removal and adjustment of headgear

For headgear, see page 25

PerformaTrak with CapStrap

Philips Respironics PerformaTrak

| Description | Quantity | Product number | Part number |
|---|----------|----------------|--------------|
| PerformaTrak single-use nasal mask, CapStrap headgear, petite | 10 | 1048709 | 989805611181 |
| PerformaTrak single-use nasal mask, CapStrap headgear, small | 10 | 1048708 | 989805611171 |
| PerformaTrak single-use nasal mask, CapStrap headgear, medium/large | 10 | 1048710 | 989805611191 |

Philips Respironics Contour Deluxe



Key benefits

- Dual-density foam forehead pads provide patient comfort and mask stability
- With only three sizes, clinicians can quickly fit a wide range of patients
- Hook-n-loop grab tabs allow quick, easy removal and adjustment of headgear

For headgear, see page 25

Contour Deluxe

Philips Respironics Contour Deluxe single-use

| Description | Quantity | Product number | Part number |
|--|----------|----------------|--------------|
| Contour Deluxe nasal mask, four-point headgear, petite | 5 | 1016980 | 989805610311 |
| Contour Deluxe nasal mask, four-point headgear, small | 5 | 1016697 | 989805610291 |
| Contour Deluxe nasal mask, four-point headgear, medium/large | 5 | 1016695 | 989805610271 |
| Contour Deluxe nasal mask, four-point headgear, small | 10 | 1016698 | 989805610301 |
| Contour Deluxe nasal mask, four-point headgear, medium/large | 10 | 1016696 | 989805610281 |



Key benefits

- Easy sizing and fitting
- Complete kit for quick setup
- Clear field of vision
- Cradle style nasal cushion

For headgear, see page 25

AP111

Philips Respironics AP111 single-use

| Description | Quantity | Product number | Part number |
|---|----------|----------------|--------------|
| AP111 interface with cradles* leak 3 | 1 | 1062128 | 989805626481 |
| AP111 interface with cradles* leak 3, circuit, and main flow filter | 1 | 1062130 | 989805626491 |
| AP111 interface with cradles* leak 3, circuit, and main flow filter | 5 | 1062132 | 989805626501 |



* Each Philips Respironics AP111 interface includes all four cushion sizes

Procedure packs and accessories

Get everything you need in one bundle with our procedure packs, designed to simplify the ordering process and help you work more efficiently.

These accessories help you customize your care to suit each individual patient, large or small.

Procedure packs

Description

Quantity Product number Part number

1

AF531 mask with EE leak 1 nebulizer elbow (no nebulizer), CapStrap headgear, and FEP circuit with main flow filter (#1)

| Small* | 5 | 1072642 | 989805634651 |
|--------|---|---------|--------------|
| Medium | 5 | 1072643 | 989805634661 |
| Large | 5 | 1072644 | 989805634671 |

AF531 mask with EE leak 1 entrainment elbow, four-point

headgear, and FEP circuit with main flow filter*

| Small* | 5 | 1072639 | 989805626511 |
|--------|---|---------|--------------|
| Medium | 5 | 1072640 | 989805628811 |
| Large | 5 | 1072641 | 989805628821 |

AF811 gel mask with EE leak 1 entrainment elbow, CapStrap

headgear, and FEP circuit with main flow filter and exhalation filter $^{\scriptscriptstyle \dagger}$

| Small | 5 | 1074839 | 989805633171 |
|--------|---|---------|--------------|
| Medium | 5 | 1074838 | 989805633161 |
| Large | 5 | 1074837 | 989805633151 |

AF811 gel mask with EE leak 1 entrainment elbow, CapStrap

headgear, and FEP circuit with main flow filter⁺

| Small | 5 | 1074061 | 989805633201 |
|--------|---|---------|--------------|
| Medium | 5 | 1074060 | 989805633191 |
| Large | 5 | 1074059 | 989805633181 |

PerformaTrak mask with EE leak 1 entrainment elbow,

four-point headgear, and FEP circuit with main flow filter⁺


| Small | 5 | 1074055 | 989805633111 |
|--------|---|---------|--------------|
| Medium | 5 | 1074054 | 989805633101 |
| Large | 5 | 1074053 | 989805633091 |

PerformaTrak mask with EE leak 1 entrainment elbow,

CapStrap headgear, and FEP circuit with main flow filter⁺

| Small | 5 | 1074058 | 989805633141 |
|--------|---|---------|--------------|
| Medium | 5 | 1074057 | 989805633131 |
| Large | 5 | 1074056 | 989805633121 |

* Intended for use by patients 7 years or older and > 20 kg (40 lb) [†] Not illustrated Note: FEP – Filter Exhalation Port



Procedure packs (continued)

| Description | Quantity | Product number | Part number |
|--|----------|----------------|--------------|
| PerformaTrak mask, NIV circuit (9624C) with main flow filter, O ₂ enrichment port ⁺ | | | |
| Small | 10 | 1004228 | 989805618031 |
| Medium | 10 | 1004229 | 989805618041 |
| Large | 10 | 1004253 | 989805618051 |
| PerformaTrak mask, heated single limb circuit with main flow filter, auto-feed humidification chamber ⁺ | | | |
| Small | 10 | 1004198 | 989805618001 |
| Medium | 10 | 1004199 | 989805618011 |

10

1004227

989805618021

Starter packs

Large

| Description | Quantity | Product number | Part number |
|--|----------|----------------|--------------|
| Philips EE starter kit, oro-nasal and full face, multi-use $^{\scriptscriptstyle +}$ | | 1074044 | 989805626551 |
| PerformaTrak EE, small/medium/large mask | 3 | | |
| PerformaTrak with CapStrap headgear, EE, small/medium/ large mask | 3 | | |
| AF811 with CapStrap headgear, EE, small/medium/large mask | 3 | | |
| PerforMax multi-use, EE/SE, four-point headgear, small/large | 2 | | |
| NIV circuit with DEP filter and main flow filter | 10 | | |
| Philips basic start kit, oro-nasal (1 mask) ⁺ | | 1074050 | 989805626611 |
| PerformaTrak SE, small/medium/large mask | 3 | | |


⁺ Not illustrated

Aerogen Pro-X controllers and systems

NIVO Pro-X controller kits*

| Description | Quantity | Product number | Part number |
|--|----------|----------------|--------------|
| NIVO Pro-X controller kit – Eng EU (English, EU plug) | 1 | 1079421 | 989805634331 |
| NIVO Pro-X controller kit – Northern EU (German, French, Dutch, Italian) | 1 | 1082693 | 989805634471 |
| NIVO Pro-X controller kit – Southern EU (Spanish, Italian, Portuguese, Greek) | 1 | 1082694 | 989805634481 |
| NIVO Pro-X controller kit – Scand EU (Norwegian, Swedish, Danish) | 1 | 1082695 | 989805634491 |
| NIVO Pro-X controller kit – UK/IE (English) |) 1 | 1079422 | 989805634341 |
| NIVO Pro-X controller kit – AU/NZ (English) | 1 | 1079423 | 989805634351 |
| NIVO Pro-X controller kit – English/INTL (English) | . 1 | 1079424 | 989805634361 |
| NIVO Pro-X controller kit – INTL (Spanish, Italian, Portuguese, Greek) | 1 | 1082696 | 989805634501 |

NIVO controller and system components



- AC/DC adapterPower source or charge controllerPro-X control moduleCan use AC/DC adapter or internal battery
- (up to 45 minutes)3 Equipment mount Use for pole or humidifier slot mounting configuration
- 4 AF531 mask^{**} Mask to be used
- 5 Control module cable Connect controller to NIVO nebulizer
- 6 Nebulizer elbows (NE) Entrainment elbow (EE) or standard elbow (SE) with nebulizer port
- 7 NIVO nebulizer Aerosol generator

1

2

* Pro-X nebulizer kit includes #1, 2, 3, and 5 shown above

** NE elbow can be added to an existing AF531 or mask can be purchased with NE elbow installed.

Mask accessories

Philips Respironics AF531 accessories and replacement parts Philips Respironics AF541 accessories and replacement parts

| Description | Quantity | Product number | Part number |
|---|----------|----------------|--------------|
| EE leak 1 entrainment elbow* | 1 | 1072645 | 989805628841 |
| EE leak 1 entrainment elbow* | 5 | 1072646 | 989805628831 |
| EE leak 1 nebulizing elbow (no nebulizer) | * 1 | 1072649 | 989805634531 |
| EE leak 1 nebulizing elbow (no nebulizer) | * 5 | 1072650 | 989805634541 |
| SE nebulizing elbow (no nebulizer)* | 1 | 1076755 | 989805634561 |
| SE nebulizing elbow (no nebulizer)* | 5 | 1076754 | 989805634551 |
| NIVO nebulizer | 5 | 1076302 | 989805634301 |
| CapStrap headgear with clips | 1 | 1072654 | 989805628871 |
| AF531 sizing gauge | 25 | 1078287 | 989805628901 |
| Swivel | 1 | 7041 | 453561517051 |
| RP-PMax and AF531 headgear clips | 2 | 1060913 | 989805649921 |
| Nasogastric sealing pad | 10 | 452040 | 453561517091 |
| CleanClip bracket assembly | 1 | 1076430 | 989805627181 |
| CleanClip bracket thumb screw | 1 | 1076197 | 989805629911 |
| Bronchoscopy elbow* | 1 | 1072651 | 989805635091 |
| Bronchoscopy elbow* | 5 | 1072652 | 989805635101 |

Philips Respironics PerforMax accessories and replacement parts

| Description | Quantity | Product number | Part number |
|------------------------------|----------|----------------|--------------|
| EE leak 2 entrainment elbow* | 1 | 1045313 | 989805611111 |
| EE leak 2 entrainment elbow* | 5 | 1072647 | 989805628851 |
| Bronchoscopy elbow* | 1 | 1072651 | 989805635091 |
| Bronchoscopy elbow* | 5 | 1072652 | 989805635101 |
| SE standard elbow* | 1 | 1043210 | 989805611081 |
| SE standard elbow* | 5 | 1072648 | 989805628861 |


| Description | Quantity | Product number | Part number |
|--|----------|----------------|--------------|
| AF541 over-the-nose cushion, small | 10 | 1120960 | 989805654131 |
| AF541 over-the-nose cushion, medium | 10 | 1120961 | 989805654141 |
| AF541 over-the-nose cushion, large | 10 | 1120962 | 989805654151 |
| AF541 over-the-nose cushion, extra large | 10 | 1120963 | 989805654161 |
| EE leak 1 elbow with NIVO port** | 5 | 1120967 | 989805654171 |
| EE leak 2 elbow with NIVO port** | 5 | 1120968 | 989805654181 |
| SE standard elbow with NIVO port** | 5 | 1120969 | 989805654191 |
| AF541 EE leak 1, without cushion, four-point headgear | 5 | 1120975 | 989805654201 |
| AF541 EE leak 1, without cushion, CapStrap headgear | 5 | 1120976 | 989805654211 |
| AF541 four-point headgear with forehead pad, small/medium | 5 | 1120970 | 989805654221 |
| AF541 four-point headgear with forehead pad, large/extra large | 5 | 1120971 | 989805654231 |
| AF541 CapStrap headgear, small/medium | า 5 | 1126500 | 989805654251 |
| AF541 CapStrap headgear, large/extra large | e 5 | 1120972 | 989805654241 |
| AF541 sizing gauge | 25 | 1120977 | 989805654261 |
| AF541 size A under-the-nose cushion | 10 | 1120980 | 989805654271 |
| AF541 size B under-the-nose cushion | 10 | 1120981 | 989805654281 |
| AF541 size C under-the nose cushion | 10 | 1120982 | 989805654291 |
| Bronchoscopy elbow | 5 | 1120984 | 989805654301 |
| AF541 discharge kit, CapStrap headgear, small/medium | 5 | 1120986 | 989805654311 |
| AF541 discharge kit, four-point headgear, small/medium | 5 | 1120987 | 989805654321 |
| AF541 discharge kit, CapStrap headgear, large/extra large | 5 | 1120988 | 989805654331 |
| AF541 discharge kit, four-point headgear, large/extra large | 5 | 1120989 | 989805654341 |

* Elbow – Conical type: Push and twist to install. Pull and twist to remove. Refer to instructions with compatible mask or elbow for proper usage. Conical type elbows used on AF531, PerforMAX and PerforMAX pediatric mask.

** Elbow – Snap/Squeeze type: Push for audible click to install and squeeze tabs to remove. Refer to instructions with compatible mask or elbow for proper usage. Snap/Squeeze type used on AF541 mask.


Headgear


| Image | Description | Size | Compatibility | Quantity | Product number | Part number |
|-------|--|-------------------|------------------------------------|----------|----------------|--------------|
| 1 | CapStrap headgear with ball and socket clips | Adult (one size) | AF811, PerformaTrak | 5 | 1019543 | 989805610841 |
| 1 | CapStrap headgear with clips | Adult (one size) | AF531 | 1 | 1072654 | 989805628871 |
| 2 | Four-point headgear with ball and socket clips | Adult (one size) | AF811, PerformaTrak | 1 | 1015788 | 989805610221 |
| 3 | Four-point headgear, no clips | Petite | AF811, PerformaTrak Contour Deluxe | 5 | 1018472 | 989805610351 |
| 3 | Four-point headgear, no clips | Adult (one size) | AF811, PerformaTrak Contour Deluxe | 5 | 1007365 | 989805609961 |
| 4 | Four-point headgear | Adult (one size) | Total | 1 | 1003052 | 453561517111 |
| 4 | Four-point headgear | Adult (one size) | Total | 9 | 1010609 | 989805610031 |
| 5 | Four-point headgear | Adult (one size) | Vinyl nasal | 10 | 312106 | 989805609291 |
| 6 | Four-point headgear with crown strap, clips | Small | PerforMax | 5 | 1051857 | 989805611231 |
| 6 | Four-point headgear with crown strap, clips | Large | PerforMax | 5 | 1051816 | 989805611221 |
| 6 | Four-point headgear with crown strap, clips | Extra large | PerforMax | 5 | 1080518 | 989805634851 |
| 6 | Four-point headgear with crown strap, clips | Small | PerforMax | 9 | 1045325 | 989805611131 |
| 6 | Four-point headgear with crown strap, clips | Large | PerforMax | 9 | 1045314 | 989805611121 |
| 6 | Four-point headgear with crown strap, clips | Extra large | PerforMax | 9 | 1080519 | 989805634861 |
| 6 | Four-point headgear with crown strap, clips | Small | AF531 | 1 | 1072655 | 989805628891 |
| 6 | Four-point headgear with crown strap, clips | Large | AF531 | 1 | 1072656 | 989805628881 |
| 6 | Four-point headgear with crown strap, clips | Extra large | AF531 | 5 | 1083188 | 989805634751 |
| 7 | Two-point headgear with chin strap | Adult (one size) | AP111 | 5 | 1062134 | 453561516551 |
| 8 | AF421 forehead pad kit with two pads | 10 mm | AF421 | 10 | 1119366 | 989805652521 |
| * | Four-point headgear assembly with forehead pad | Adult (one size) | PerformaTrak multi-use only | 5 | 1022453 | 989805610861 |
| 9 | AF541 four-point headgear with forehead pad | Small/Medium | AF541 | 5 | 1120970 | 989805654221 |
| 9 | AF541 four-point headgear with forehead pad | Large/Extra large | AF541 | 5 | 1120971 | 989805654231 |
| 10 | AF541 CapStrap headgear | Small/Medium | AF541 | 5 | 1126500 | 989805654251 |
| 10 | AF541 CapStrap headgear | Large/Extra large | AF541 | 5 | 1120972 | 989805654241 |




Circuit and mask components and accessories

| Image | Description | Quantity | Product number | Part number |
|-------|---|----------|----------------|--------------|
| 1 | Plateau exhalation valve (reusable) | 1 | 302312 | 989805617941 |
| 2 | Plateau exhalation valve diaphragm | 5 | 302310 | 989805609221 |
| 3 | Whisper Swivel II | 1 | 332113 | 989805617951 |
| 4 | Disposable exhalation port with filter exhalation port with swivel (no filter) | 10 | 1065775 | 453561517211 |
| 4a | Replacement port caps (silicone) | 5 | 1061825 | 453561517201 |
| 4b | Disposable exhalation port with filter exhalation port and proximal pressure line (7 ft, 2.1 m) | 10 | 1116891 | 989805652291 |
| 4b | Disposable exhalation port with filter exhalation port and proximal pressure line (7 ft, 2.1 m) | 1 | 1116892 | 989805652301 |
| 5 | Disposable exhalation port | 10 | 312149 | 989805609391 |
| 6 | Oxygen enrichment attachment | 1 | 312710 | 989805609421 |
| 6 | Oxygen enrichment attachment | 10 | 312010 | 989805609241 |
| 7 | Bacteria/virus filter (single use) | 1 | 1014047 | 989805618161 |
| 7 | Bacteria/virus filter (single use) | 10 | 342077 | 989805609521 |
| 8 | Proximal airway filter | 1 | 1002362 | 453561517101 |
| 9 | Tube management clip | 2 | 1039967 | 453561517351 |


























Circuit and mask components and accessories

| Image | Description | Quantity | Product number | Part number |
|-------|---|----------|----------------|--------------|
| 10 | Swivel replacement kit (22 mm male-male) | 1 | 7041 | 453561517051 |
| 11 | CleanClip bracket assembly | 1 | 1076430 | 989805627181 |
| 11 | CleanClip bracket thumb screw | 1 | 1076197 | 989805629911 |
| 12 | Nasogastric tube sealing pad | 10 | 452040 | 453561517091 |
| 13 | PerforMax sizing gauge, small, large, and extra large | 1 | 1095097 | 989805648901 |
| 14 | Replacement ball and socket | 2 | 1009723 | 989805652511 |




















Circuits and connections

Philips offers a variety of single-use circuits, including heated-wire and non-heated-wire options, to accompany our noninvasive ventilation systems.

Circuits

| Image | Description | Product number | Part number | Quantity | V60 | BiPAP Vision | Trilogy 202 | BiPAP Focus | BiPAP S/T-D |
|-------|--|-------------------|--------------|----------|-----|-----------------|----------------|----------------|----------------|
| 1 | Invasive circuit with exhalation port, water trap and temperature probe ports, proximal airway filter and humidifier coupling tube, hanger, and hose clips | 652002 | 989805609681 | 10 | • | • | • | | |
| 1 | Invasive circuit with exhalation port, water trap and temperature probe ports, proximal airway filter and humidifier coupling tube, hanger, and hose clips | 652001 | 989805609671 | 20 | • | • | • | | |
| 2 | Noninvasive circuit with main flow filter, disposable swivel exhalation port with filter port and filter, proximal pressure line, hanger, and hose clips | 1065830 | 989805621311 | 10 | • | • | • | • | |
| 3 | Noninvasive circuit with main flow filter, disposable swivel exhalation port with filter port, proximal pressure line, hanger, and hose clips (exhalation filter not included) | 1065832 | 989805621321 | 10 | • | • | • | • | |
| 4 | Noninvasive circuit with disposable swivel exhalation port with filter port, proximal pressure line, hanger, and hose clips (main flow and exhalation filters not included) | 1069210 | 989805634871 | 10 | • | • | • | • | |
| 5 | Noninvasive circuit with disposable exhalation port, proximal pressure line, hanger, and hose clips | 582073 | 989805609611 | 10 | • | • | • | • | |
| 6 | BiPAP S/T-D circuit with exhalation port, proximal pressure line, hanger, and hose clips | 312107 | 989805609301 | 10 | | | | | • |
| 7 | Tubing 6 ft (22 mm) (no dep, proximal pressure line, or hanger) | 532316 | 989805647181 | 1 | | | | | |



Humidification circuits and accessories

| Description | Quantity | Product number | Part number |
|--|----------|----------------|--------------|
| Auto-feed humidification chamber | 24 | 1035609 | 989805619481 |
| Heated single-limb circuit with main flow filter | 10 | 8624C | 989805619521 |
| Heated single-limb circuit with main flow filter, auto-feed humidification chamber | 10 | 1045770 | 989805619511 |

Dry circuits and accessories

| Description | Quantity | Product number | Part number |
|---|----------|----------------|--------------|
| NIV circuit with main flow filter | 10 | 9624 | 989805619411 |
| NIV circuit with main flow filter, O ₂ , enrichment port | 10 | 9624C | 989805619531 |

Trilogy accessories

The innovative technology of Trilogy sets it apart from the rest. The support, service, and upkeep behind it keep it there.

Easy to use and maintain

- Trilogy circuits can be used for a variety of pediatric and adult patient applications
- Passive circuits for NIV applications
- Active exhalation valve for IV application
- \cdot Circuits are available in packs of 10 and 20





Circuit adapters

| Image | Description | Compatibility | Product number | Part number |
|-------|---|--|--|------------------------------|
| 1 | Active exhalation valve with PAP assembly | | 1065659 – 1/pk | 989805628221 |
| 2 | Proximal flow sensor (Trilogy202) | | 1050408 – 1/pk 1050534 – 10/pk | 989805621211 989805645511 |
| 3 | Disposable swivel passive exhalation port | | 1084428 – box of 10 1084429 – box of 20 | 989805643161 989805646961 |
| 4 | Flexible trach adapter with 15 mm cuff | Connects to the flow sensor or the active exhalation valve | 1073863 – box of 10 1073903 – box of 20 | 989805633921 989805633941 |
| 5 | Flexible trach adapter with 22 mm cuff | Connects to Whisper Swivel II and swivel and fixed disposable exhalation ports | 1073902 – box of 10 1073904 – box of 20 | 989805633931 989805633951 |

The ability to use either an active valve or passive exhalation port gives you the flexibility to choose the best available circuit and interface for your patient. Trilogy's unique passive exhalation port reduces extra tubing connections.

5







3





Miscellaneous


| Imag | ge Product | Product number | Part number |
|------|-------------------------------|----------------|--------------|
| 6 | Oxygen inlet adapter – 2/pk | 1040390 | 989805647691 |
| 7 | Tube port adapter, active PAP | 1061398 | 989805621291 |

Trilogy support arm system

8


4

| Image | Product | Product number | Part number |
|-------|------------------------------|-----------------------|--------------|
| 8 | Tubing hanger (22 mm) – 1/pk | 311107 | 989805648891 |
| * | Support arm | 332497 | 989805617961 |
| * | Support arm mount assembly | 1052863 | 989805605731 |












* Image not shown


Trilogy disposable circuits


| Image | Description | Components | Quantity | Product number | Part number |
|-------|--|---|------------------------|--------------------|------------------------------|
| 1 | Adult passive circuit with fixed exhalation port | Passive exhalation port. Patient tube: 6 ft length, 22 mm diameter with 22 mm cuff. | Box of 10 Box of 20 | 1069467 1069468 | 989805621551 989805621561 |
| 2 | Adult passive circuit with swivel design exhalation port | Passive exhalation port. Patient tube: 6 ft length, 22 mm diameter with 22 mm cuff. | Box of 10 Box of 20 | 1073228 1073236 | 989805633861 989805633911 |
| 3* | Adult passive circuit with water trap | Standard 90° degree elbow adapter. Passive fixed exhalation port. Temperature probe port. Patient tube, 2 ft length, 22 mm diameter with 22 mm cuff. Water trap. Humidifier connection tube, 2 ft length, 22 mm diameter with 22 mm cuff. | Box of 10 Box of 20 | 1069443 1069444 | 989805621531 989805621541 |
| 4* | Adult passive circuit, heated wire | Passive swivel exhalation port. Temperature probe port. Heated wire patient tube, 5 ft length, 22 mm diameter with 22 mm cuff. Heated wire adapter. Humidifier connection tube, 2 ft length, 22 mm in diameter with 22 mm cuff. | Box of 10 Box of 20 | 1073227 1073235 | 989805646441 989805646481 |
| 5 | Adult active circuit, no water trap | Standard 90° elbow trach adapter. Active exhalation valve. Exhalation tubing. Proximal pressure tubing. Patient tube, 6 ft length, 22 mm diameter with 22 mm cuff. Tube management clips. | Box of 10 Box of 20 | 1073223 1073232 | 989805633841 989805633881 |
| 6* | Adult active circuit, heated wire | Standard 90° elbow trach adapter. Active exhalation valve. Exhalation tubing. Proximal pressure tubing. Heated wire patient tube, 5 ft length, 22 mm diameter with 22 mm cuff, with heated wire connector. Heated wire adapter. Temperature probe port. Humidifier connection tube, 2 ft length, 15 mm diameter with 22 mm cuff. Tube management clips. | Box of 10 Box of 20 | 1073222 1073231 | 989805646421 989805646461 |

















* Chamber not included

Trilogy disposable circuits


| Image | Description | Components | Quantity | Product number | Part number |
|-------|--|---|------------------------|-----------------------|------------------------------|
| 7 | Pediatric passive circuit with water trap | Passive fixed exhalation port. Patient tube, 6 ft length, 15 mm diameter with 22 mm cuff and water trap. | Box of 10 Box of 20 | 1069445 1069446 | 453561517231 453561517241 |
| 8* | Pediatric passive circuit with water trap | Passive swivel exhalation port. Patient tube, 2.5 ft length, 15 mm diameter with 22 mm cuff. Temperature probe port. Water trap. Humidifier connection tube, 2 ft length, 15 mm diameter with 22 mm cuff. | Box of 10 Box of 20 | 1073224 1073233 | 989805633851 989805633901 |
| 9* | Pediatric passive circuit, heated wire | Passive swivel exhalation port. Heated wire patient tube, 5 ft length, 15 mm diameter with 22 mm cuff, with heated wire connector. Heated wire adapter. Temperature probe port. Humidifier connection tube, 2 ft length, 15 mm diameter with 22 mm cuff. | Box of 10 Box of 20 | 1073225 1073234 | 989805646431 989805646471 |
| 10* | Pediatric active circuit with water trap | Standard 90° elbow trach adapter. Active exhalation valve. Exhalation tubing. Proximal pressure tubing. Patient tube, 2.5 ft length, 15 mm diameter with 22 mm cuff. Temperature probe port. Water trap. Humidifier connection tube, 2 ft length, 15 mm diameter with 22 mm cuff. Tube management clips (one shown). | Box of 10 Box of 20 | 1073220 1073229 | 989805633831 989805633871 |
| 11* | Pediatric active circuit, heated wire | Standard 90° elbow trach adapter. Active exhalation valve. Exhalation tubing. Proximal pressure tubing. Heated wire patient tube, 5 ft length, 15 mm diameter with 22 mm cuff, with heated wire connector. Heated wire adapter. Temperature probe port. Tube management clips. Humidifier connection tube, 2 ft length, 15 mm diameter with 22 mm cuff. | Box of 10 Box of 20 | 1073221 1073230 | 989805646411 989805646451 |
| 12 | MPV (mouthpiece ventilation) circuit, single use | Bacteria filter, 22 mm x 15 mm reducer cuff, 6 ft length, 15 mm tubing, coupler, miniature and angled mouthpiece. | Box of 10 | 1102866 | 989805646981 |








11



12



* Chamber not included

Notes

Notes

©2016 Koninklijke Philips N.V. All rights are reserved. Philips reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.



www.philips.com/NIVmasks

Printed in the Netherlands. 4522 991 16521 * AUG 2016